

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT	A1
(Equipment locations, heat loads, component weights, environmental specs)	
STRUCTURAL LAYOUT	S1
(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	
STRUCTURAL DETAILS	S2
(Floor and Ceiling loading information)	
ELECTRICAL LAYOUT	E1
(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	
ELECTRICAL SPECIFICATIONS	E2
(Maximum wiring run lengths, interconnect diagram, system power specifications)	
ELECTRICAL DETAILS	E3
EQUIPMENT DETAILS	D1

These equipment IS drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the IS and operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

*** REQUIRED REFERENCE ***

Precision RXi
Pre Installation Manual

5123449-100

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



R/F Site Planning



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 18					
<i>Before using this document ensure you have the latest Rev from MyWorkshop on DCCC429762</i>					
GEHC Global Order # :		Customer:			
GEHC PMI :		FE / Installer:			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.					
Item #	Inspection Date:	Storage Is Item Ready?	PM Is Item Ready?	FE Is Item Ready?	Comments If 'N', enter comments or action plan
1					Pre-Project Delivery Requirements: Ensure oxygen venting system is designed and installed with objective evidence that it is compliant with the GEHC Pre-Installation Manual (PIM) requirements, exhaust system is installed and operational, 480V power, and chilled water supply is available that meets system cooling requirements. Electrical connectivity is available for magnet monitoring and phone service is available during delivery.
2					Pre-Install Screen Room Requirements: If Screen Room is tested with objective evidence that it is compliant with GEHC specifications. Do not Bolt installed using 2 post anchors. For rfp system, slotted beam cannot be installed by SF vendor using 2 post anchors.
3					Site Regulatory Requirements: Site Drawing Requirements: final version of equipment installation drawings (including red lined version) verified to match actual room and has been provided to installer. X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO & WA.
4					Site Drawing Requirements: Final version of equipment installation drawings (including red lined version) verified to match actual room and has been provided to installer.
5					Surface Penetration Requirements: Customer/contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls. DR surface penetration permit available and posted in the room when GEHC will perform the work.
6					Delivery Route Requirements: The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protectors, fork lift, rollback truck, etc.)
7					Finished Room Requirements: Rooms that will contain equipment, including storage areas, and clean rooms, are dust free. Precautions taken to maintain rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility. For Storage: Room must meet PIM requirements for storage.
8					Electrical Requirements: Main Disconnect Panel (MDP) is installed and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and lead-side wires can be installed at time of system installation.
9					HVAC Requirements: The HVAC/Chilled Water systems designed to maintain the environment are running and appear to provide the desired environmental conditions (temperature and humidity) for system operation.
10					Flooring Requirements: Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications.
11					Ceiling Requirements: Unit for equivalent location, levelness and spacing is measured per vendor confirmed and consistent with the requirements of the installation drawings. Ceiling grids installed. Permanent lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PIM direction.

GE Healthcare

IS Services Design Center

Minneapolis, Wisconsin
Copyright © 2009, General Electric Company. Proprietary to GE.

SHEET TITLE: SITE READINESS
MODALITY TYPE: PRECISION RXi

THIS PLAN IS SUBMITTED TO CURRENT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST REVISED DRAWINGS. THIS PLAN IS NOT TO BE USED FOR ANY CONSTRUCTION PURPOSES. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL FINAL DRAWINGS
3-50F

PROJECT	REVISION
3-50F	09
DATE:	20.Apr.11
DRAWN BY:	JDR
CHECKED BY:	REK

REVISION HISTORY:

SHEET
C1

RQ - 117726 PM R12

GE EQUIPMENT LISTING

EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER : NEITHER A QUOTE OR CON WAS ISSUED AT THE DATE OF THESE DRAWINGS

NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.

ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN
1	1		X-RAY TUBE ON OVERHEAD TUBE SUSPENSION SYSTEM	496 lbs	511 btu	BD115F C7509	---	DTS S
2	1		DTS STATIONARY RAILS	198 lbs		BD115F	BD115G	C
3	1		TV MONITOR ON MOBILE MONITOR CART	44 lbs	34 btu	C7508	---	TV S
4	1		TABLE CABINET	374 lbs	1245 btu	BD115E	---	PC S
5	1		PRECISION RX: ELEVATING TABLE	2200 lbs	211 btu	BD115A BD115C	BD115C	RXI S
6	1		INTEGRATED CONSOLE	365 lbs	296 btu	BD115M	---	IC S
7	1		PRECISION RX: GENERATOR	235 lbs	75 btu	BD115D	---	GEN S
8	1		PRECISION RX: NON-TILTING WALLSTAND	374 lbs	126 btu	BD115L	---	WB S


EQUIPMENT CROSS REFERENCE CHART
 P = PREAPPROVAL
 C = CALCULATIONS/PENDING APPROVAL
 S = SPECIFICATIONS ONLY

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

1	1		STEP DOWN TRANSFORMER	198 lbs		4502KP1	---	SDT S
2	1		MAIN DISCONNECT CONTROL PANEL GEMS CAT. NO. E4502KP	169 lbs		E4502KP	---	MDC S

EQUIPMENT LAYOUT RECOMMENDED CEILING HEIGHT = 9'-10"

SCALE: 1/4" = 1'-0"
 This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



Drawn by: JOYDEL ROELKE Octel no.: 5603733
 GE Installation
 Project Manager: VINSON MARTIN
 Telephone no.: (281) 852-8641

ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	COUNTER TOP WITH SINK, BASE AND WALL CABINETS
61	COUNTER TOP FOR EQUIPMENT - MINIMUM DEPTH 24 IN. AND ADDITIONAL SHELVING MAY BE REQUIRED BELOW COUNTER TOP FOR PC TOWER. PROVIDE GROMMETTED OPENINGS AS REQUIRED TO ROUTE CABLES.
62	CONTROL WALL, 7 FT. HIGH WITH LEAD GLASS VIEWING WINDOW.
63	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 48 IN. W X 80 IN. H (1219mm X 2032mm). CONTINGENT ON A 78 IN. (1981mm) CORRIDOR WIDTH.
64	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL: 800-200-9760 GE CAT. NO. WX1ABW-WF-XIU
65	DOOR LIMIT SWITCH (NEEDED ONLY IF REQUIRED BY STATE/LOCAL CODES)

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL. REFERENCE JUNCTION POINT "XRL6" ON SHEET "E1" FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.
----	--

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.

SITE ENVIRONMENT SPECIFICATIONS

- AMBIENT OPERATING TEMPERATURE: 59 TO 75 DEGREES (F), MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 15 DEGREES (F)/HOUR.
- HUMIDITY: REFER TO PREINSTALLATION MANUAL FOR THE EQUIPMENT ILLUSTRATED ON THIS DRAWING.
- ALTITUDE: NOT TO EXCEED 8,000 FT. ABOVE SEA LEVEL.
- THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
- DO NOT RESTRICT THE AIR INTAKE AT THE LOWER FRONT OR AIR EXHAUST AT THE TOP OF THE ELECTRONICS CABINETS.

MAGNETIC INTERFERENCE SPECIFICATIONS

IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.
 X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.
 SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.
 OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

GE Healthcare
 IS Services Design Center
 Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT LAYOUT
 MODALITY TYPE: PRECISION RX:
 THIS PLAN IS SUBMITTED TO GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN. EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST GEHC INSTRUCTIONS. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL FINAL DRAWINGS
3-50F

PROJECT	REVISION
3-50F	09


DATE: 20.Apr.11
 DRAWN BY: JDR
 CHECKED BY: REK

REVISION HISTORY:

SHEET
A1


TYPICAL WALL SUPPORT ELEVATIONS

S115




SUPPORT FOR TABLE CABINET (NOT TO SCALE)

S110



SUPPORT FOR STEP DOWN TRANSFORMER (NOT TO SCALE)

S113




SUPPORT FOR MAIN DISCONNECT CONTROL (NOT TO SCALE)

SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-10"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING STATIONARY RAILS. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 27.5" AND REQUIRE 450 LBS (197 LBS IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
2	SUPPORT BACKING. LOCATE AS SHOWN. REFER TO ELEVATION DETAIL S115, FOR TABLE CABINET. WALL BACKING MUST SUPPORT 110lb (50kg) PER ANCHOR LOAD.
3	FLOOR CONTACT AREA FOR PRECISION RX.
4	SUPPORT BACKING. LOCATE AS SHOWN. REFER TO ELEVATION DETAIL S110, FOR STEP DOWN TRANSFORMER.
5	SUPPORT BACKING. LOCATE AS SHOWN. REFER TO ELEVATION DETAIL S113, FOR MAIN DISCONNECT PANEL.
6	FLOOR CONTACT AREA FOR CHEST UNIT

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1.58mm (1/16") DEFLECTION. (10) 12.7mm (1/2") DIA. x 38.1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12.7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6.35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3.17mm (1/8") IN 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT
MODALITY TYPE: PRECISION RX

PROJECT TITLE:
TYPICAL FINAL DRAWINGS
3-50F

PROJECT	REVISION
3-50F	09

DATE: 20.Apr.11
DRAWN BY: JDR
CHECKED BY: REK


REVISION HISTORY:

SHEET
S1

Drawn by: JOYDEL ROELKE Octel no.: 5603733
GE Installation
Project Manager: VINSON MARTIN
Telephone no.: (281) 852-8641


PIM R12
RQ - 117726

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST REVISED AND INSTALLED TO DATE SPECIFICATIONS. THE COMPANY CANNOT ACCEPT LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.


GE Healthcare
 IS Services Design Center
 Milwaukee, Wisconsin

SUPPORT DETAIL
PRECISION RXi RADIOGRAPHIC SUSPENSION MOUNTING


B0115G
REV. DATE: 12/07/04



DETAIL NOT TO SCALE

SUPPORT DETAIL
PRECISION RXi RADIOGRAPHIC SUSPENSION MOUNTING

B0115H
REV. DATE: 2/10/05




UNISTRUT SYSTEM IS NOT SUPPLIED OR INSTALLED BY GENERAL ELECTRIC MEDICAL SYSTEMS.

DETAIL NOT TO SCALE

SUPPORT DETAIL
PRECISION RXi RADIOGRAPHIC SUSPENSION MOUNTING


B0115H1
REV. DATE: 2/10/05



DETAIL NOT TO SCALE

FLOOR MOUNTING DETAIL: PRECISION RXi ELEVATING TABLE INSTALLATION


B0115C
REV. DATE: 02/02/06



DETAIL NOT TO SCALE

FLOOR MOUNTING DETAIL: PRECISION RXi CHEST UNIT INSTALLATION METHODS

B0115J
REV. DATE: 11/23/09



DETAIL NOT TO SCALE

GE Healthcare
IS Services Design Center
Minneapolis, Wisconsin

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: PRECISION RXi

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND REGULATIONS. THE COMPANY CANNOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL FINAL DRAWINGS 3-50F

PROJECT	REVISION
3-50F	09

DATE: 20.Apr.11
DRAWN BY: JDR
CHECKED BY: REK

REVISION HISTORY:

SHEET
S2

SCALE: 1/4" = 1'-0"

ELECTRICAL PLAN

RECOMMENDED CEILING HEIGHT = 9'-10"

ELECTRICAL OUTLET LEGEND
CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS. HEIGHT ABOVE FLOOR DETERMINED BY LOCAL CODES UNLESS OTHERWISE SPECIFIED.

- △ DEDICATED TELEPHONE LINE(S) (SEE ELECTRICAL DETAIL ELEC-1 OR ELEC-87)
- ▲ NETWORK OUTLET (SEE ELECTRICAL DETAILS ELEC-83 AND ELEC-84 OR ELEC-87)

DUCT HATCHING LEGEND





- ▨ ABOVE CEILING DUCT
- ▩ UNDER FLOOR DUCT
- ▧ TRENCH DUCT (FLUSH FLOOR)
- ▦ SURFACE FLOOR DUCT
- ▥ CABLE TRAY
- ABOVE CEILING CONDUIT
- BELOW FLOOR CONDUIT

- JUNCTION POINT NOTES**
- o ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER'S ELECTRICAL CONTRACTOR.
 - o CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS.
 - o CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
 - o CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
 - o ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
 1. DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
 2. DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL POWER PURPOSES.
 3. DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
 4. PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
 - o ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMER'S CONTRACTOR.
 - o GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
 - o 10 FOOT PIGTAILS AT ALL JUNCTION POINTS.
 - o ALL WIRING MUST BE THIN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT INSULATION. **ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.**
 - o GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.


JUNCTION POINT DESCRIPTIONS

THE FOLLOWING MATERIALS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER'S ELECTRICAL CONTRACTOR

DESCRIPTION	QTY.	HARDWARE	DETAIL NO., SHT. E3
DLK1 DOOR SWITCH IN FRAME - NORMALLY OPEN (24V) REQUIRED BY STATE/LOCAL CODES?	1	ROOM DOOR INTERLOCK LIMIT SWITCH	ELEC-5
GEN GENERATOR	1	32 IN. OF GROMMET MATERIAL FOR AN 8 X 8 IN. OPENING IN DUCT COVER	ELEC-6
IC INDEPENDENT CONSOLE	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
MDC 100-AMP PANEL INCLUDED IN ORDER	1	100-AMP CIRCUIT BREAKER PANEL GENS CAT. NO. E4502RP EMERGENCY PUSHBUTTON STATION IS INCLUDED.	ELEC-6
DTS X-RAY TUBE HANGER	1	32 IN. OF GROMMET MATERIAL FOR AN 8 X 8 IN. OPENING IN DUCT COVER	ELEC-5
PC TABLE CONTROL CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8 X 8 IN. OPENING IN DUCT COVER	ELEC-6
RML1 ROOM LIGHTS AVAILABLE FROM GE. CALL: 800-558-5102	1	COVERPLATE SINGLE GANG BOX #E4502SS 24V X-RAY ROOM WARNING LIGHT AND ROOM LIGHT CONTROLLER OR EQUIVALENT.	ELEC-17
RXI X-RAY TABLE	3	3 1/2 IN. NIPPLES, 1 1/2 IN. LONG	ELEC-25
SDT STEP DOWN TRANSFORMER	1	FITTINGS AS REQUIRED	
SEO EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8 IN. DEEP, FLUSH MTD. WALL BOX	ELEC-16
TV TV MONITOR	1	BLANK COVERPLATE 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
WB CHEST UNIT	1	BLANK COVERPLATE 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-79
XRL1 WARNING LIGHT	1	SINGLE GANG BOX X-RAY ON-INCANDESCENT LIGHT FIXTURE. 24V, 8 AMP OR LESS LOW VOLTAGE SOURCE. DO NOT USE FLUORESCENT FIXTURES.	ELEC-17
XRLC WARNING LIGHT CONTROLLER AVAILABLE FROM GE. CALL: 800-558-5102 OR LOCAL GE. INSTALLATION PROJECT MGR.	1	E4502SS WARNING LIGHT & ROOM LIGHT CONTROL	ELEC-17



PLEASE SEE BELOW FOR ADDITIONAL REQUIRED CONDUIT RUNS AND SIZES.



FEEDER TABLE INDICO 100 80 kW REV. DATE: 05/Apr/11

o CALCULATIONS BASED UPON NOMINAL VOLTAGE, WIRE SIZE IN AWG.
o RECOMMENDED FEEDER SIZES FROM DISTRIBUTION TRANSFORMER TO THE POWER CABINET.
o NEUTRAL MUST BE TERMINATED INSIDE THE MAIN DISCONNECT PANEL AND NOT AT ANY GE CABINET.
o THE GROUNDING CONDUCTOR () WILL BE OF SAME SIZE AS THE FEEDER WIRES WITH A 1/0 MINIMUM. THIS GROUND WILL RUN FROM THE EQUIPMENT BACK TO THE FACILITY POWER SOURCE/MAIN GROUNDING POINT AND ALWAYS TRAVEL IN THE SAME CONDUIT WITH THE FEEDERS AND NEUTRAL.
o IF THE GENERAL ELECTRIC EQUIPMENT IS BEING FED BY A DELTA SECONDARY, IT IS RECOMMENDED THAT THE B PHASE ON THE SECONDARY BE CONNECTED TO GROUND TO PREVENT DAMAGE TO THE SYSTEM.
o * MINIMUM SIZE FOR CIRCUIT BREAKER, NEC ARTICLE S17-73.
o FOR A FULL SYSTEM UPS, REFER TO ELECTRICAL DETAILS FOR UPS FEEDER WIRES.

RUN LENGTH IN FEET	POWER SUPPLY VOLTAGE	
	360-440 400	432-528 480
50	* 2 (1/0)	* 2 (1/0)
100	* 2 (1/0)	* 2 (1/0)
150	2 (1/0)	* 2 (1/0)
200	1/0 (1/0)	2 (1/0)
250	2/0 (2/0)	1 (1/0)
300	3/0 (3/0)	1/0 (1/0)
350	4/0 (4/0)	2/0 (2/0)
400	4/0 (4/0)	2/0 (2/0)

ADDITIONAL CONDUIT RUNS FOR PRECISION RXI (BY CONTRACTOR)

CONDUITS REQUIRED FOR BASE SYSTEM (CONDUITS ARE LOCATED ABOVE CEILING)

REV DATE: 10/18/10

XRLC TO RML1	ONE 1/2" CND.
XRLC TO XRL1	ONE 1/2" CND.
XRLC TO GEN	ONE 1/2" CND.
XRLC TO 120-V 1Ø POWER	CND. AS REQ'D
DLK1 TO GEN	ONE 1/2" CND.
MDC TO GEN	ONE CND. AS REQ'D.
MDC TO GEN/PC	ONE 2" CND.
MDC TO OTS	ONE 2" CND.
MDC TO FEEDER	ONE CND. AS REQ'D
MDC TO SEO	ONE 1/2" CND.
MDC TO SDT	ONE CND. AS REQ'D

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

Drawn by: JOYDEL ROELKE Octel no.: 5603733
GE Installation
Project Manager: VINSON MARTIN
Telephone no.: (281) 852-8641

CONTRACTOR SUPPLIED AND INSTALLED WIRING
ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS.

WIRE RUN, FROM - TO	QUANTITY, WIRE SIZE/COLOR
XRLC > 1 PHASE	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
XRLC > RML1	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
XRL1 > XRLC	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
DLK1 > GEN	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
XRLC > GEN	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
MDC > GEN	3-BLACK, 1 GREEN - REFER TO FEEDER TABLE
480-V > MDC	3-BLACK, 1-WHITE, 1-GREEN - SIZE AS REQUIRED
MDC > TRANS	2-ND. 12 BLACK, 1-ND. 12 GREEN
TRANS > MDC	2-ND. 8 BLACK, 1-ND. 8 WHITE, 1-ND. 8 GREEN

NOTE: IF SUPPLIED CABLE (16.4 FT.) IS NOT LONG ENOUGH, SUPPLY THE FOLLOWING:
3-BLACK, 1 GREEN - REFER TO FEEDER TABLE

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL LAYOUT
MODALITY TYPE: PRECISION RXI
THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND REQUIREMENTS OF THE PROJECT. IT IS THE RESPONSIBILITY OF THE USER TO VERIFY ALL CONSTRUCTION DETAILS AND REQUIREMENTS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL FINAL DRAWINGS
3-50F


PROJECT	REVISION
3-50F	09

DATE: 20.Apr.11
DRAWN BY: JDR
CHECKED BY: REK

REVISION HISTORY:

SHEET
E1

INTERCONNECT DIAGRAM



REV DATE 06 Apr 11

POWER SPECIFICATIONS

INDICO 100 SERIES 80 kW GENERATOR SYSTEM

VOLTAGE PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS.
 RANGE OF LINE VOLTAGES
 NOMINAL LINE VOLTAGE OF 400 & 480, 3 PHASE, 50 OR 60 Hz
 RECOMMENDED POWER SUPPLY: DELTA OR WYE-CONNECTED.
 MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

TABLE A ALLOWABLE INPUT VOLTAGES/CURRENT DEMAND

NOMINAL VOLTAGE	NORMAL RANGE ±10 PERCENT	CURRENT (AMPS)		MINIMUM STANDARD OVERCURRENT PROTECTION
		MAX. MOMENTARY	CONTINUOUS	
400	360-440	152	3	100-A
480	432-528	126	3	100-A

NOTE LOW LINE CONDITIONS MAY INHIBIT SOME HIGH KVp TECHNIQUES. THE GENERATOR AUTOMATICALLY ESTABLISHES THESE INHIBITS BASED ON ACTUAL LINE CONDITIONS AND SYSTEM REGULATION.

PHASE-BALANCE. PHASE-TO-PHASE VOLTAGES MUST BE WITHIN +2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED LINE VOLTAGE AT A MAXIMUM DURATION OF 5 CYCLES AND FREQUENCY OF 10 TIMES PER HOUR.

POWER DEMAND CONTINUOUS POWER DEMAND = 2.4 KVA

TABLE B MAXIMUM MOMENTARY POWER DEMAND.

DEMAND	INDICO
kVA * POWER FACTOR AT	105 n/a
mA	1000
kVp	120

* DEMAND INCLUDES POWER FOR ENTIRE GENERATOR SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 10 PERCENT.

DISTRIBUTION TRANSFORMER FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 150 KVA.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

DIAGRAM KEY

- CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN ADEQUATE CONDUIT OR RACEWAY.
- GE FURNISHED CABLE RUNS. ROUTE IN EMPTY CONDUIT OR RACEWAY.
- 59' [18M] MAXIMUM RUN LENGTH BETWEEN JUNCTION POINTS. Feet [Meters]

SHEET TITLE: ELECTRICAL SPECIFICATIONS
 MODALITY TYPE: PRECISION RXi

THIS PLAN IS SUBMITTED TO SURVEY LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST REVISIONS OF THE NATIONAL ELECTRICAL CODE. THE COMPANY CANNOT ACCEPT LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
 TYPICAL FINAL DRAWINGS
 3-50F

PROJECT	REVISION
3-50F	09

DATE: 20.Apr.11
 DRAWN BY: JDR
 CHECKED BY: REK

REVISION HISTORY:

SHEET
 E2

GE Healthcare

IS Services Design Center
 Milwaukee, Wisconsin

RQ - 117726 PIM R12

ELECTRICAL DETAIL HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5
REV. DATE: 03/19/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL VERTICAL WALL DUCT (TYPICAL)

ELEC-6
REV. DATE: 03/19/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL FLUSH FLOOR DUCT (TYPICAL)

ELEC-25
REV. DATE: 4/01/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL INSITE CONNECTION (TYPICAL)

ELEC-1
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.

A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.

OR

B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND ONE A DEDICATED DATA LINE.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL BOX WITH COVERPLATE AND NETWORK JACK

ELEC-83
REV. DATE: 10/06/98

DETAIL NOT TO SCALE

ELECTRICAL DETAIL NETWORK CONNECTION (TYPICAL)

ELEC-84
REV. DATE: 03/06/04

FOR NUCLEAR SYSTEMS A DIRECT NETWORK CONNECTION IS TO BE MADE BETWEEN THE SYSTEM AND THE REVIEW WORKSTATION.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL J.B. / WALL DUCT DETAIL (TYPICAL)

ELEC-2
REV. DATE: 09/30/94

DETAIL NOT TO SCALE

ELECTRICAL DETAIL X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL

ELEC-17
REV. DATE: 05/14/09

E4502SS
X-RAY ROOM WARNING LIGHT / ROOM LIGHTING CONTROL PANEL

TO GE X-RAY ON SIGNAL: MAXIMUM 24-VAC

TO GE ROOM LIGHT SIGNAL: MAXIMUM 24-VAC

120-VAC 10A MAXIMUM

X-RAY WARNING LIGHT: 0-VAC

ROOM LIGHTS: 0-VAC

10A MAX

10A MAX

THE R4500AL IS RECOMMENDED IF "X-RAY ON" WARNING LIGHT AND ROOM LIGHT CONTROL ARE UTILIZED

E4502RL
X-RAY ROOM WARNING LIGHT CONTROL PANEL

TO GE X-RAY ON SIGNAL: MAXIMUM 24-VAC

120-VAC 10A MAXIMUM

X-RAY WARNING LIGHT: 0-VAC

CONTROL PANEL CAN BE LOCATED ABOVE THE CEILING NEAR THE WARNING LIGHT

UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER, ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

ELECTRICAL DETAIL EMERGENCY OFF BUTTON

ELEC-16
REV. DATE: 05/14/09

DETAIL NOT TO SCALE

ELECTRICAL DETAIL BOX WITH DIVIDER AND SPLIT COVERPLATE (TYPICAL)

ELEC-79
REV. DATE: 04/06/04

DETAIL NOT TO SCALE


ELECTRICAL DETAIL BOX WITH COVERPLATE (TYPICAL)

ELEC-8
REV. DATE: 09/30/94

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi X-RAY TUBE SUSPENSION


B0115F
REV. DATE: 12/07/04



PLAN VIEW
DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi X-RAY TUBE SUSPENSION


C75-09
REV. DATE: 11/23/09



DETAIL NOT TO SCALE

EQUIPMENT DETAIL
LCD MONITOR MOBILE CART


C75-08
REV. DATE: 11/23/09



DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi TABLE CABINET

B0115E
REV. DATE: 02/14/05



FRONT VIEW


PLAN VIEW

SIDE VIEW

**CUSTOMER SUPPLIED HARDWARE:
QTY. (2) 12mm DIA. WALL ANCHORS
50mm LONG AND QTY. (2) 6mm DIA.
CABINET BOLTS WITH NUTS 50mm LONG
DRILL WALL ANCHOR HOLE 12mm DIA.
60mm DEEP.

EQUIPMENT DETAIL
PRECISION RXi TABLE-ELEVATING


B0115
REV. DATE: 11/16/09



FRONT VIEW
DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi TABLE-ELEVATING


B0115A
REV. DATE: 11/16/09



SIDE VIEW
DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi - CONSOLE COMPONENTS

B0115M
REV. DATE: 11/20/09



LCD MONITOR


GENERATOR TOUCH SCREEN

POSITIONER CONSOLE

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi GENERATOR CABINET

B0115D
REV. DATE: 11/20/09



FRONT VIEW

RIGHT SIDE VIEW


REAR VIEW

GENERATOR CLEARANCES

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
STEP-DOWN TRANSFORMER


E4502KP1
REV. DATE: 11/16/09



DETAIL NOT TO SCALE

EQUIPMENT DETAIL
SYSTEM MAIN DISCONNECT


E4502KP
REV. DATE: 08/20/04



DETAIL NOT TO SCALE

EQUIPMENT DETAIL
PRECISION RXi NON-TILTING WALLSTAND

B0115L
REV. DATE: 11/23/09



DETAIL NOT TO SCALE

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: PRECISION RXi

PROJECT TITLE:
TYPICAL FINAL DRAWINGS
3-50F

PROJECT	REVISION
3-50F	09

DATE: 20.Apr.11
DRAWN BY: JDR
CHECKED BY: REK

REVISION HISTORY:

SHEET
D1

PIM R12
RQ - 117726

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin