

Discovery CT590 RT  
Optima CT580  
Discovery RT

## Pre-Installation Manual

OPERATING DOCUMENTATION



5366636-1EN  
Rev 18  
© 2009-2017, GE Company  
All Right Reserved.

## Effectivity

The information in this manual applies to the following RT CT Systems:

- Discovery CT590 RT
- Optima CT580
- Discovery RT

# IMPORTANT PRECAUTIONS

## LANGUAGE

|                               |  |
|-------------------------------|--|
| <b>ПРЕДУПРЕЖДЕНИЕ</b><br>(BG) | <p>Това упътване за работа е налично само на английски език.</p> <ul style="list-style-type: none"><li>• Ако доставчикът на услугата на клиента изиска друг език, задължение на клиента е да осигури превод.</li><li>• Не използвайте оборудването, преди да сте се консултирали и разбрали упътването за работа.</li><li>• Неспазването на това предупреждение може да доведе до нараняване на доставчика на услугата, оператора или пациента в резултат на токов удар, механична или друга опасност.</li></ul>   |
| <b>警告</b><br>(ZH-CN)          | <p>本维修手册仅提供英文版本。</p> <ul style="list-style-type: none"><li>• 如果客户的维修服务人员需要非英文版本，则客户需自行提供翻译服务。</li><li>• 未详细阅读和完全理解本维修手册之前，不得进行维修。</li><li>• 忽略本警告可能对维修服务人员、操作人员或患者造成电击、机械伤害或其他形式的伤害。</li></ul>   |
| <b>警告</b><br>(ZH-HK)          | <p>本服務手冊僅提供英文版本。</p> <ul style="list-style-type: none"><li>• 倘若客戶的服務供應商需要英文以外之服務手冊，客戶有責任提供翻譯服務。</li><li>• 除非已參閱本服務手冊及明白其內容，否則切勿嘗試維修設備。</li><li>• 不遵從本警告或會令服務供應商、網絡供應商或病人受到觸電、機械性或其他危險。</li></ul>   |
| <b>警告</b><br>(ZH-TW)          | <p>本維修手冊僅有英文版。</p> <ul style="list-style-type: none"><li>• 若客戶的維修廠商需要英文版以外的語言，應由客戶自行提供翻譯服務。</li><li>• 請勿試圖維修本設備，除非 您已查閱並瞭解本維修手冊。</li><li>• 若未留意本警告，可能導致維修廠商、操作員或病患因觸電、機械或其他危險而受傷。</li></ul>  |
| <b>UPOZORENJE</b><br>(HR)     | <p>Ovaj servisni priručnik dostupan je na engleskom jeziku.</p> <ul style="list-style-type: none"><li>• Ako davatelj usluge klijenta treba neki drugi jezik, klijent je dužan osigurati prijevod.</li><li>• Ne pokušavajte servisirati opremu ako niste u potpunosti pročitali i razumjeli ovaj servisni priručnik.</li><li>• Zanimajte li ovo upozorenje, može doći do ozljede davatelja usluge, operatera ili pacijenta uslijed strujnog udara, mehaničkih ili drugih rizika.</li></ul>  |
| <b>VÝSTRAHA</b><br>(CS)       | <p>Tento provozní návod existuje pouze v anglickém jazyce.</p> <ul style="list-style-type: none"><li>• V případě, že externí služba zákazníkům potřebuje návod v jiném jazyce, je zajištění překladu do odpovídajícího jazyka úkolem zákazníka.</li><li>• Nesnažte se o údržbu tohoto zařízení, aniž byste si přečetli tento provozní návod a pochopili jeho obsah.</li><li>• V případě nedodržování této výstrahy může dojít k poranění pracovníka prodejního servisu, obslužného personálu nebo pacientů vlivem elektrického proudu, respektive vlivem mechanických či jiných rizik.</li></ul> |

|                             |  |
|-----------------------------|--|
| <b>ADVARSEL</b><br>(DA)     | Denne servicemanual findes kun på engelsk. <ul style="list-style-type: none"><li>• Hvis en kundes tekniker har brug for et andet sprog end engelsk, er det kundens ansvar at sørge for oversættelse.</li><li>• Forsøg ikke at servicere udstyret uden at læse og forstå denne servicemanual.</li><li>• Manglende overholdelse af denne advarsel kan medføre skade på grund af elektrisk stød, mekanisk eller anden fare for teknikeren, operatøren eller patienten.</li></ul>  |
| <b>WAARSCHUWING</b><br>(NL) | Deze onderhoudshandleiding is enkel in het Engels verkrijgbaar. <ul style="list-style-type: none"><li>• Als het onderhoudspersoneel een andere taal vereist, dan is de klant verantwoordelijk voor de vertaling ervan.</li><li>• Probeer de apparatuur niet te onderhouden alvorens deze onderhoudshandleiding werd geraadpleegd en begrepen is.</li><li>• Indien deze waarschuwing niet wordt opgevolgd, zou het onderhoudspersoneel, de operator of een patiënt gewond kunnen raken als gevolg van een elektrische schok, mechanische of andere gevaren.</li></ul>                               |
| <b>WARNING</b><br>(EN)      | This service manual is available in English only. <ul style="list-style-type: none"><li>• If a customer's service provider requires a language other than english, it is the customer's responsibility to provide translation services.</li><li>• Do not attempt to service the equipment unless this service manual has been consulted and is understood.</li><li>• Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.</li></ul>  |
| <b>HOIATUS</b><br>(ET)      | See teenindusjuhend on saadaval ainult inglise keeles <ul style="list-style-type: none"><li>• Kui klienditeeninduse osutaja nõuab juhendit inglise keelest erinevas keeles, vastutab klient tõlketeenuse osutamise eest.</li><li>• Ärge üritage seadmeid teenindada enne eelnevalt käesoleva teenindusjuhendiga tutvumist ja sellest aru saamist.</li><li>• Käesoleva hoiatuse eiramine võib põhjustada teenuseosutaja, operaatori või patsiendi vigastamist elektrilöögi, mehaanilise või muu ohu tagajärjel.</li></ul>   |
| <b>VAROITUS</b><br>(FI)     | Tämä huolto-ohje on saatavilla vain englanniksi. <ul style="list-style-type: none"><li>• Jos asiakkaan huoltohenkilöstö vaatii muuta kuin englanninkielistä materiaalia, tarvittavan käännöksen hankkiminen on asiakkaan vastuulla.</li><li>• Älä yritä korjata laitteistoa ennen kuin olet varmasti lukenut ja ymmärtänyt tämän huolto-ohjeen.</li><li>• Mikäli tätä varoitusta ei noudateta, seurauksena voi olla huoltohenkilöstön, laitteiston käyttäjän tai potilaan vahingoittuminen sähköiskun, mekaanisen vian tai muun vaaratilanteen vuoksi.</li></ul>                                   |
| <b>ATTENTION</b><br>(FR)    | Ce manuel d'installation et de maintenance est disponible uniquement en anglais. <ul style="list-style-type: none"><li>• Si le technicien d'un client a besoin de ce manuel dans une langue autre que l'anglais, il incombe au client de le faire traduire.</li><li>• Ne pas tenter d'intervenir sur les équipements tant que ce manuel d'installation et de maintenance n'a pas été consulté et compris.</li><li>• Le non-respect de cet avertissement peut entraîner chez le technicien, l'opérateur ou le patient des blessures dues à des dangers électriques, mécaniques ou autres.</li></ul> |

|                                       |   |
|---------------------------------------|---|
| <p><b>WARNUNG</b><br/>(DE)</p>        | <p>Diese Serviceanleitung existiert nur in englischer Sprache.</p> <ul style="list-style-type: none"> <li>Falls ein fremder Kundendienst eine andere Sprache benötigt, ist es Aufgabe des Kunden für eine entsprechende Übersetzung zu sorgen.</li> <li>Versuchen Sie nicht diese Anlage zu warten, ohne diese Serviceanleitung gelesen und verstanden zu haben.</li> <li>Wird diese Warnung nicht beachtet, so kann es zu Verletzungen des Kundendiensttechnikers, des Bedieners oder des Patienten durch Stromschläge, mechanische oder sonstige Gefahren kommen.</li> </ul>  |
| <p><b>ΠΡΟΕΙΔΟΠΟΙΗΣΗ</b><br/>(EL)</p>  | <p>Το παρόν εγχειρίδιο σέρβις διατίθεται μόνο στα αγγλικά.</p> <ul style="list-style-type: none"> <li>Εάν ο τεχνικός σέρβις ενός πελάτη απαιτεί το παρόν εγχειρίδιο σε γλώσσα εκτός των αγγλικών, αποτελεί ευθύνη του πελάτη να παρέχει τις υπηρεσίες μετάφρασης.</li> <li>Μην επιχειρήσετε την εκτέλεση εργασιών σέρβις στον εξοπλισμό αν δεν έχετε συμβουλευτεί και κατανοήσει το παρόν εγχειρίδιο σέρβις.</li> <li>Αν δεν προσέξετε την προειδοποίηση αυτή, ενδέχεται να προκληθεί τραυματισμός στον τεχνικό σέρβις, στο χειριστή ή στον ασθενή από ηλεκτροπληξία, μηχανικούς ή άλλους κινδύνους.</li> </ul>                     |
| <p><b>FIGYELMEZTETÉS</b><br/>(HU)</p> | <p>Ezen karbantartási kézikönyv kizárólag angol nyelven érhető el.</p> <ul style="list-style-type: none"> <li>Ha a vevő szolgáltatója angoltól eltérő nyelvre tart igényt, akkor a vevő felelőssége a fordítás elkészíttetése.</li> <li>Ne próbálja elkezdni használni a berendezést, amíg a karbantartási kézikönyvben leírtakat nem értelmezték.</li> <li>Ezen figyelmeztetés figyelmen kívül hagyása a szolgáltató, működtető vagy a beteg áramütés, mechanikai vagy egyéb veszélyhelyzet miatti sérülését eredményezheti.</li> </ul>  |
| <p><b>AÐVÖRUN</b><br/>(IS)</p>        | <p>Þessi þjónustuhandbók er aðeins fáanleg á ensku.</p> <ul style="list-style-type: none"> <li>Ef að þjónustuveitandi viðskiptamanns þarfnast annas tungumáls en ensku, er það skylda viðskiptamanns að skaffa tungumálþjónustu.</li> <li>Reynið ekki að afgreiða tækið nema að þessi þjónustuhandbók hefur verið skoðuð og skilin.</li> <li>Brot á sinna þessari aðvörun getur leitt til meiðsla á þjónustuveitanda, stjórnanda eða sjúklings frá raflosti, vélrænu eða öðrum áhættum.</li> </ul>  |
| <p><b>AVVERTENZA</b><br/>(IT)</p>     | <p>Il presente manuale di manutenzione è disponibile soltanto in lingua inglese.</p> <ul style="list-style-type: none"> <li>Se un addetto alla manutenzione richiede il manuale in una lingua diversa, il cliente è tenuto a provvedere direttamente alla traduzione.</li> <li>Procedere alla manutenzione dell'apparecchiatura solo dopo aver consultato il presente manuale ed averne compreso il contenuto.</li> <li>Il mancato rispetto della presente avvertenza potrebbe causare lesioni all'addetto alla manutenzione, all'operatore o ai pazienti provocate da scosse elettriche, urti meccanici o altri rischi.</li> </ul> |
| <p><b>警告</b><br/>(JA)</p>             | <p>このサービスマニュアルには英語版しかありません。</p> <ul style="list-style-type: none"> <li>サービスを担当される業者が英語以外の言語を要求される場合、翻訳作業はその業者の責任で行うものとさせていただきます。</li> <li>このサービスマニュアルを熟読し理解せずに、装置のサービスを行わないでください。</li> <li>この警告に従わない場合、サービスを担当される方、操作員あるいは患者さんが、感電や機械的又はその他の危険により負傷する可能性があります。</li> </ul>   |

|                            |  |
|----------------------------|--|
| <b>경고</b><br>(KO)          | <p>본 서비스 매뉴얼은 영어로만 이용하실 수 있습니다 .</p> <ul style="list-style-type: none"> <li>• 고객의 서비스 제공자가 영어 이외의 언어를 요구할 경우, 번역 서비스를 제공하는 것은 고객의 책임입니다 .</li> <li>• 본 서비스 매뉴얼을 참조하여 숙지하지 않은 이상 해당 장비를 수리하려고 시도하지 마십시오 .</li> <li>• 본 경고 사항에 유의하지 않으면 전기 쇼크 , 기계적 위험 , 또는 기타 위험으로 인해 서비스 제공자 , 사용자 또는 환자에게 부상을 입힐 수 있습니다 .</li> </ul>  |
| <b>BRDINJUMS</b><br>(LV)   | <p>Šī apkopes rokasgrāmata ir pieejama tikai angļu valodā.</p> <ul style="list-style-type: none"> <li>• Ja klienta apkopes sniedzējam nepieciešama informācija citā valodā, klienta pienākums ir nodrošināt tulkojumu.</li> <li>• Neveiciet aprīkojuma apkopi bez apkopes rokasgrāmatas izlasīšanas un saprašanas.</li> <li>• Šī brīdinājuma neievērošanas rezultātā var rasties elektriskās strāvas trieciena, mehānisku vai citu faktoru izraisītu traumu risks apkopes sniedzējam, operatoram vai pacientam.</li> </ul>   |
| <b>ĮSPĖJIMAS</b><br>(LT)   | <p>Šis eksploatavimo vadovas yra tik anglų kalba.</p> <ul style="list-style-type: none"> <li>• Jei kliento paslaugų tiekėjas reikalauja vadovo kita kalba – ne anglų, suteikti vertimo paslaugas privalo klientas.</li> <li>• Nemėginkite atlikti įrangos techninės priežiūros, jei neperskaitėte ar nesupratote šio eksploatavimo vadovo.</li> <li>• Jei nepaisysite šio įspėjimo, galimi paslaugų tiekėjo, operatoriaus ar paciento sužalojimai dėl elektros šoko, mechaninių ar kitų pavojų.</li> </ul>   |
| <b>ADVARSEL</b><br>(NO)    | <p>Denne servicehåndboken finnes bare på engelsk.</p> <ul style="list-style-type: none"> <li>• Hvis kundens serviceleverandør har bruk for et annet språk, er det kundens ansvar å sørge for oversettelse.</li> <li>• Ikke forsøk å reparere utstyret uten at denne servicehåndboken er lest og forstått.</li> <li>• Manglende hensyn til denne advarselen kan føre til at serviceleverandøren, operatøren eller pasienten skades på grunn av elektrisk støt, mekaniske eller andre farer.</li> </ul>  |
| <b>OSTRZEŻENIE</b><br>(PL) | <p>Niniejszy podręcznik serwisowy dostępny jest jedynie w języku angielskim.</p> <ul style="list-style-type: none"> <li>• Jeśli serwisant klienta wymaga języka innego niż angielski, zapewnienie usługi tłumaczenia jest obowiązkiem klienta.</li> <li>• Nie próbować serwisować urządzenia bez zapoznania się z niniejszym podręcznikiem serwisowym i zrozumienia go.</li> <li>• Niezastosowanie się do tego ostrzeżenia może doprowadzić do obrażeń serwisanta, operatora lub pacjenta w wyniku porażenia prądem elektrycznym, zagrożenia mechanicznego bądź innego.</li> </ul> |
| <b>ATENÇÃO</b><br>(PT-BR)  | <p>Este manual de assistência técnica encontra-se disponível unicamente em inglês.</p> <ul style="list-style-type: none"> <li>• Se outro serviço de assistência técnica solicitar a tradução deste manual, caberá ao cliente fornecer os serviços de tradução.</li> <li>• Não tente reparar o equipamento sem ter consultado e compreendido este manual de assistência técnica.</li> <li>• A não observância deste aviso pode ocasionar ferimentos no técnico, operador ou paciente decorrentes de choques elétricos, mecânicos ou outros.</li> </ul>                              |

|                                    |  |
|------------------------------------|--|
| <p><b>ATENÇÃO</b><br/>(PT-PT)</p>  | <p>Este manual de assistência técnica só se encontra disponível em inglês.</p> <ul style="list-style-type: none"> <li>• Se qualquer outro serviço de assistência técnica solicitar este manual noutro idioma, é da responsabilidade do cliente fornecer os serviços de tradução.</li> <li>• Não tente reparar o equipamento sem ter consultado e compreendido este manual de assistência técnica.</li> <li>• O não cumprimento deste aviso pode colocar em perigo a segurança do técnico, do operador ou do paciente devido a choques eléctricos, mecânicos ou outros.</li> </ul>  |
| <p><b>ATENȚIE</b><br/>(RO)</p>     | <p>Acest manual de service este disponibil doar în limba engleză.</p> <ul style="list-style-type: none"> <li>• Dacă un furnizor de servicii pentru clienți necesită o altă limbă decât cea engleză, este de datoria clientului să furnizeze o traducere.</li> <li>• Nu încercați să reparați echipamentul decât ulterior consultării și înțelegerii acestui manual de service.</li> <li>• Ignorarea acestui avertisment ar putea duce la rănirea depanatorului, operatorului sau pacientului în urma pericolelor de electrocutare, mecanice sau de altă natură.</li> </ul>   |
| <p><b>ОСТОРОЖНО!</b><br/>(RU)</p>  | <p>Данное руководство по техническому обслуживанию представлено только на английском языке.</p> <ul style="list-style-type: none"> <li>• Если сервисному персоналу клиента необходимо руководство не на английском, а на каком-то другом языке, клиенту следует самостоятельно обеспечить перевод.</li> <li>• Перед техническим обслуживанием оборудования обязательно обратитесь к данному руководству и поймите изложенные в нем сведения.</li> <li>• Несоблюдение требований данного предупреждения может привести к тому, что специалист по техобслуживанию, оператор или пациент получит удар электрическим током, механическую травму или другое повреждение.</li> </ul> |
| <p><b>UPOZORENJE</b><br/>(SR)</p>  | <p>Ovo servisno uputstvo je dostupno samo na engleskom jeziku.</p> <ul style="list-style-type: none"> <li>• Ako klijentov serviser zahteva neki drugi jezik, klijent je dužan da obezbedi prevodilačke usluge.</li> <li>• Ne pokušavajte da opravite uređaj ako niste pročitali i razumeli ovo servisno uputstvo.</li> <li>• Zanemarivanje ovog upozorenja može dovesti do povređivanja serviser, rukovaoca ili pacijenta usled strujnog udara ili mehaničkih i drugih opasnosti.</li> </ul>   |
| <p><b>UPOZORNENIE</b><br/>(SK)</p> | <p>Tento návod na obsluhu je k dispozícii len v angličtine.</p> <ul style="list-style-type: none"> <li>• Ak zákazníkovi poskytovateľ služieb vyžaduje iný jazyk ako angličtinu, poskytnutie prekladateľských služieb je zodpovednosťou zákazníka.</li> <li>• Nepokúšajte sa o obsluhu zariadenia, kým si neprečítate návod na obsluhu a neporozumiete mu.</li> <li>• Zanedbanie tohto upozornenia môže spôsobiť zranenie poskytovateľa služieb, obsluhujúcej osoby alebo pacienta elektrickým prúdom, mechanické alebo iné ohrozenie.</li> </ul>   |

|                          |  |
|--------------------------|--|
| <b>ATENCION</b><br>(ES)  | <p>Este manual de servicio sólo existe en inglés.</p> <ul style="list-style-type: none"><li>• Si el encargado de mantenimiento de un cliente necesita un idioma que no sea el inglés, el cliente deberá encargarse de la traducción del manual.</li><li>• No se deberá dar servicio técnico al equipo, sin haber consultado y comprendido este manual de servicio.</li><li>• La no observancia del presente aviso puede dar lugar a que el proveedor de servicios, el operador o el paciente sufran lesiones provocadas por causas eléctricas, mecánicas o de otra naturaleza.</li></ul> |
| <b>VARNING</b><br>(SV)   | <p>Den här servicehandboken finns bara tillgänglig på engelska. .</p> <ul style="list-style-type: none"><li>• Om en kunds servicetekniker har behov av ett annat språk än engelska, ansvarar kunden för att tillhandahålla översättningstjänster.</li><li>• Försök inte utföra service på utrustningen om du inte har läst och förstår den här servicehandboken.</li><li>• Om du inte tar hänsyn till den här varningen kan det resultera i skador på serviceteknikern, operatören eller patienten till följd av elektriska stötar, mekaniska faror eller andra faror.</li></ul>         |
| <b>OPOZORILO</b><br>(SL) | <p>Ta servisni priročnik je na voljo samo v angleškem jeziku.</p> <ul style="list-style-type: none"><li>• Če ponudnik storitve stranke potrebuje priročnik v drugem jeziku, mora stranka zagotoviti prevod.</li><li>• Ne poskušajte servisirati opreme, če tega priročnika niste v celoti prebrali in razumeli.</li><li>• Če tega opozorila ne upoštevate, se lahko zaradi električnega udara, mehanskih ali drugih nevarnosti poškoduje ponudnik storitev, operater ali bolnik.</li></ul>   |
| <b>DIKKAT</b><br>(TR)    | <p>Bu servis kılavuzunun sadece ingilizcesi mevcuttur.</p> <ul style="list-style-type: none"><li>• Eğer müşteri teknisyeni bu kılavuzu ingilizce dışında bir başka lisandan talep ederse, bunu tercüme ettirmek müşteriye düşer.</li><li>• Servis kılavuzunu okuyup anlamadan ekipmanlara müdahale etmeyiniz.</li><li>• Bu uyarıya uyulmaması, elektrik, mekanik veya diğer tehlikelerden dolayı teknisyen, operatör veya hastanın yaralanmasına yol açabilir.</li></ul>   |



## DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. June 13, 2006

## CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Healthcare, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Healthcare, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS



**CAUTION**  
**Risk of**  
**Explosion**

**Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**



**ATTENTION**  
**Danger**  
**d'Explosion**

**Il y a danger d'explosion s'il y a remplacement incorrect de la batterie. Remplacer uniquement avec une batterie du même type ou d'un type recommandé par le constructeur. Mettre au rebut les batteries usagées conformément aux instructions du fabricant.**

## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

# Revision History

| Revision | Date     | Reason for change  |
|----------|----------|--|
| 18       | 08/15/17 | <b>Chapter 5:</b><br>Section 1.0 Updated OC and Table Net Weight.<br><b>Chapter 6:</b><br>Section 2.0 Add subsystem packaging dimension.   |
| 17       | 03/02/17 | <b>Chapter 1:</b><br>Section 3.0 Updated Pre-Installation Checklist according to DOC1809666_r2.<br><b>Chapter 3:</b><br>Section 1.0 Update Figure 3-1 Regulatory Clearance Requirement for gantry left side.<br>Section 3.0 Add note for minimum room requirement.<br>Section 4.0 Updated OC minimum clearances to wall from 6 in. to 4 in.<br>Section 5.0 Update Floor Strength to delete 28 days requirement.<br>Section 6.0 Added host computer rear ports figure for Z840.<br><b>Chapter 4:</b><br>Section 5.0 Added Operator Console Noise level information.<br><b>Chapter 5:</b><br>Section 1.0 Updated Table 5-1 System Floor Loading.<br><b>Chapter 6:</b><br>Section 7.1 Updated Storage relative humidity according to SRS (DOC1019695_r10).<br>Section 8.0 Updated system to temperatures or humidity outside of the following specifications according to SRS (DOC1019695_r10).<br><b>Chapter 8:</b><br>Section 3.0 Updated NIO16 Console Cables. |
| 16       | 08/20/15 | Add Product Name for Discovery RT.<br><b>Chapter 5:</b><br>Update anchor information according to DOC1687529.  |
| 15       | 06/02/15 | <b>Chapter 5:</b><br>Section 1.1 Added new anchor (5487992-2) information on this section<br><b>Chapter 6:</b><br>Section 4.2 Updated door opening value   |
| 14       | 10/10/14 | <b>Chapter 5:</b><br>Section 1.0 Updated Table Weight Table 5-1<br><b>Chapter 6:</b><br>Section 1.0 Added rigging statement per global requirement<br><b>Chapter 8:</b><br>Section 3.0 Added Long and Short Cable Kit for new power pan phase in   |

| Revision | Date     | Reason for change  |
|----------|----------|--|
| 13       | 11/04/13 | <b>Chapter 7:</b><br>Section 3.0 Updated purchasable option transformer BCAT# on Table 7-2   |
| 12       | 09/16/13 | <b>Chapter 3:</b><br>Section 2.2.1 Added Gantry Service Clearance per global service requirement<br>Section 4.2.7 Added Smart Workspace Desk information<br><b>Chapter 5:</b><br>Section 1.0 Added table (SWS) weight and size on Table 5-1<br><b>Chapter 8:</b><br>Section 6.0 Corrected PDU information on Table 8-13  |
| 11       | 06/18/13 | <b>Chapter 3:</b><br>Section 4.2.2 Updated table dimension on Figure 3-10 and 3-11<br><b>Chapter 4:</b><br>Section 5.0: Added Section 5.0 for System Component Noise Levels per program requirement<br><b>Chapter 5:</b><br>Section 1.0: Updated GT1700 Table weight to 453kg on Table 5-1 per Uehara's mail<br>Section 2.0: Updated GT1700 table dimension on Figure 5-5<br><b>Chapter 6:</b><br>Section 8.0: Updated GT1700 Table weight to 453kg on Table 6-2 and Table 6-3 per Uehara's mail |
| 10       | 03/15/13 | Change table name from " <b>GT650lbs to High capacity table</b> " on whole manual according to RA requirements<br><b>Chapter 3:</b><br>Section 4.2: Updated GT1700 and High capacity Table values<br>Section 7.0: Updated head and body scatter plot<br><b>Chapter 4:</b><br>Section 2.0: Updated Table 4-1 cooling requirement<br><b>Chapter 6:</b><br>Section 7.0: Updated door opening values per service feedback  |
| 9        | 11/19/12 | <b>Chapter 3:</b><br>Section 2.2.4: Added the Note for storage cabinet as option<br><b>Chapter 8:</b><br>Section 3.3: Updated Table 8-7 for NIO console cables information   |
| 8        | 02/23/12 | <b>Chapter 1:</b><br>Section 4.0: Updated 4.1 General Scope for IEC3 requirements  |
| 7        | 11/09/11 | The following changes are from SPR HCSDM00107061<br><b>Chapter 8:</b><br>Section 3.0: Added Table 8-3 GE Supplied Cables (Standard Run) (2281840-14)- UL Information and Table 8-5 GE Supplied Cables (Optional, Long Run) (2281840-13)- UL Information  |

| Revision | Date     | Reason for change  |
|----------|----------|--|
| 6        | 07/14/11 | <b>Chapter 2:</b><br>Deleted minimum concrete floor thickness for High Capacity Table is 127 mm (5") note<br><b>Chapter 3:</b><br>Deleted this chapter<br><b>Chapter 5:</b><br>Updated Figure 5-2 Typical Floor Anchor, Table (High Capacity Table)<br>Added Figure 5-10 NIO16 Console Center-of-Gravity<br><b>Chapter 6:</b><br>Section 10.0: Added NIO16 Console Considerations<br><b>Chapter 8:</b><br>Section 3.0: Added System Interconnect Diagram with NIO16 Console<br>Section 5.0: Added 2218570-4 BOM Fuse Kit |
| 5        | 12/02/10 | <b>Chapter 8:</b><br>Section 2.2: Updated (average) power demand at maximum duty cycle to 30kVA per LSD feedback   |
| 4        | 08/23/10 | <b>Chapter 5:</b><br>Section 2.0: Updated Total system heat consumption on Table 5-1   |
| 3        | 04/13/10 | Update the Manual's title due to 8 slice phase in  |
| 2        | 12/18/09 | <b>Chapter 4:</b><br>Updated Figure 4-1 Regulatory Clearance Requirements and Figure 4-3 Minimum Service Clearances<br><b>Chapter 5:</b><br>Updated Figure 5-1 Sample Room Layout, showing approximate EMI requirements  |
| 1        | 10/10/09 | Initial Release.   |



# Table of Contents

|   |           |
|---|-----------|
| <b>Chapter 1</b>  |           |
| <b>Introduction.....</b>  | <b>21</b> |
| <b>Section 1.0</b>  |           |
| <b>Site Readiness.....</b>  | <b>21</b> |
| <b>Section 2.0</b>  |           |
| <b>Responsibility of Purchaser.....</b>   | <b>22</b> |
| 2.1 Customer Room Prep Items .....  | 22        |
| 2.2 Purchaser Site Preparation Work .....   | 23        |
| 2.3 Manufacturer's System Level Siting Requirements .....   | 23        |
| 2.3.1 Meeting Site Ready Requirements .....   | 24        |
| 2.3.2 Quick Installs .....  | 24        |
| 2.3.3 "Two-Step" and Upgrade Installs .....   | 24        |
| 2.3.4 Site Ready Inspection Visit .....   | 25        |
| <b>Section 3.0</b>  |           |
| <b>Pre-Installation Checklist. ....</b>   | <b>26</b> |
| <b>Section 4.0</b>  |           |
| <b>Medical Electrical Equipment for EMC.....</b>  | <b>28</b> |
| 4.1 General Scope.....  | 28        |
| 4.2 Electromagnetic Emission .....  | 28        |
| 4.3 Electromagnetic Immunity .....  | 29        |
| 4.3.1 Limitations Management :.....   | 31        |
| 4.4 Use Limitation: .....   | 31        |
| 4.4.1 External components .....   | 31        |
| 4.5 Installation Requirements & Environment Control :.....  | 32        |
| 4.5.1 Cable Shielding & Grounding .....   | 32        |
| 4.5.2 This product complies the radiated emission as per CISPR11 Group 1 Class A standard limits 32 |           |
| 4.5.3 Subsystem & Accessories Power Supply Distribution .....                                       | 32        |
| 4.5.4 Stacked Components & Equipment.....   | 32        |
| 4.5.5 Low Frequency Magnetic Field.....   | 32        |
| 4.5.6 Static Magnetic Field Limits .....  | 32        |
| 4.5.7 Electrostatic Discharge Environment & Recommendations .....                                   | 32        |
| <b>Chapter 2</b>  |           |
| <b>Pre-Installation Overview.....</b>   | <b>33</b> |
| <b>Section 1.0</b>  |           |
| <b>Dust/Dirt Contamination.....</b>   | <b>33</b> |
| <b>Section 2.0</b>  |           |
| <b>Chemical Contamination .....</b>   | <b>33</b> |
| <b>Section 3.0</b>  |           |
| <b>Walls, Ceiling, and Floor.....</b>   | <b>33</b> |

|   |           |
|---|-----------|
| <b>Section 4.0</b>  |           |
| <b>Broad-band</b> .....   | <b>34</b> |
| <b>Section 5.0</b>  |           |
| <b>Phone Line (for optional modem)</b> .....                          | <b>34</b> |
| <b>Section 6.0</b>  |           |
| <b>Review</b> .....   | <b>34</b> |
| <br>  |           |
| <b>Chapter 3</b>  |           |
| <b>Room Planning</b> .....  | <b>37</b> |
| <br>  |           |
| <b>Section 1.0</b>  |           |
| <b>Regulatory Clearances</b> .....                                    | <b>37</b> |
| 1.1 Regulations .....   | 37        |
| 1.2 Clearance Requirements .....                                      | 37        |
| 1.2.1 Minimum Regulatory Workspace Clearances by Major Subsystem..... | 39        |
| 1.3 How to Measure .....  | 42        |
| 1.4 NEC Conduit and Duct Fill Rate.....                               | 43        |
| 1.5 Terms and Definitions .....                                       | 44        |
| <br>  |           |
| <b>Section 2.0</b>  |           |
| <b>Service Clearances Requirements</b> .....                          | <b>46</b> |
| 2.1 Measuring Service Clearances .....                                | 46        |
| 2.2 Special Service Clearances .....                                  | 47        |
| 2.2.1 Gantry Service Clearance .....                                  | 47        |
| 2.2.2 Cover Removal.....  | 47        |
| 2.2.3 Power Distribution Unit (NGPDU) .....                           | 48        |
| 2.2.4 Console .....   | 48        |
| 2.2.5 Storage Cabinet .....   | 48        |
| <br>  |           |
| <b>Section 3.0</b>  |           |
| <b>Room Size</b> .....  | <b>50</b> |
| 3.1 Room Dimensions.....  | 50        |
| 3.1.1 Suggested Room Size.....  | 50        |
| 3.1.2 Typical Room Size .....   | 50        |
| 3.1.3 Minimum Room Size .....   | 50        |
| 3.2 Suggested and Typical Room Layouts .....                          | 52        |
| 3.3 Minimum Room Layouts .....  | 53        |
| 3.4 Control Room Considerations .....                                 | 55        |
| 3.4.1 Typical Control Room Layout.....                                | 55        |
| 3.5 CT Simulation Laser Alignment Systems .....                       | 56        |
| <br>  |           |
| <b>Section 4.0</b>  |           |
| <b>System Component Dimensions</b> .....                              | <b>57</b> |
| 4.1 Minimum Operating Clearances.....                                 | 57        |
| 4.1.1 Ceiling Pedestal Mount Installation .....                       | 57        |
| 4.1.2 Injector Control Installation.....                              | 57        |
| 4.1.3 System Operational Clearances.....                              | 58        |
| 4.2 Component Dimensions .....  | 59        |
| 4.2.1 System Dimensions.....  | 59        |
| 4.2.2 Table and Gantry Dimensions.....                                | 61        |



|       |   |    |
|-------|---|----|
| 4.2.3 | Power Distribution Unit Dimensions..... | 63 |
| 4.2.4 | Operator Console Dimensions.....        | 64 |
| 4.2.5 | Freedom Workspace Table.....            | 66 |
| 4.2.6 | Optima Table .....                      | 67 |
| 4.2.7 | Smart Workspace Desk .....              | 68 |

## **Section 5.0**

### **Structural Requirements ..... 69**

|       |   |    |
|-------|---|----|
| 5.1   | Table and Gantry Mounting Requirements..... | 69 |
| 5.2   | Floor Anchors .....                         | 69 |
| 5.3   | Floor Strength .....                        | 70 |
| 5.4   | Floor Levelness .....                       | 70 |
| 5.5   | Floor Vibration .....                       | 70 |
| 5.5.1 | Steady State Vibration .....                | 70 |
| 5.5.2 | Transient Vibration.....                    | 70 |
| 5.5.3 | Equipment Location .....                    | 70 |
| 5.6   | Walls: Scan Window .....                    | 71 |

## **Section 6.0**

### **Network Connections ..... 72**

|     |  |    |
|-----|--|----|
| 6.1 | US Broad Band Process Overview.....        | 75 |
| 6.2 | Customer Broad-Band Responsibilities ..... | 75 |

## **Section 7.0**

### **Radiation Protection ..... 76**

## **Chapter 4**

### **Environmental Conditions ..... 79**

#### **Section 1.0**

#### **Temperature and Humidity Specifications ..... 79**

|     |  |    |
|-----|--|----|
| 1.1 | Temperature (Scan Room & Control Room) ..... | 79 |
| 1.2 | Humidity (Scan Room & Control Room) .....    | 80 |
| 1.3 | Other Guideline.....                         | 80 |

#### **Section 2.0**

#### **Cooling Requirement ..... 81**

#### **Section 3.0**

#### **Altitude ..... 81**

#### **Section 4.0**

#### **Electro-Magnetic Interference (EMI) ..... 82**

|     |                                   |    |
|-----|-----------------------------------|----|
| 4.1 | Gantry .....                      | 82 |
| 4.2 | Console / Computer Equipment..... | 82 |
| 4.3 | Magnetic Media .....              | 82 |
| 4.4 | PDU .....                         | 82 |
| 4.5 | EMI Reduction .....               | 82 |
| 4.6 | UPS .....                         | 83 |
| 4.7 | Equipment EMI "Envelopes" .....   | 83 |

#### **Section 5.0**

#### **System Component Noise Levels..... 85**

## Chapter 5

### Floor Loading and Weights ..... 87

#### Section 1.0

#### Floor Loading..... 87

##### 1.1 Floor Loading and Anchoring Guidelines ..... 89

##### 1.2 GE Supplied Anchoring ..... 89

##### 1.2.1 Specifications of GE-supplied Anchors ..... 89

#### Section 2.0

#### Mounting Data, Including Seismic ..... 91

## Chapter 6

### Delivery and Storage Requirements..... 103

#### Section 1.0

#### Delivery to the Facility..... 103

##### 1.1 Loading Dock Deliveries (Preferred methods) ..... 103

##### 1.2 Ground (Non-Loading Dock) Deliveries ..... 103

##### 1.2.1 Lift-Gate Truck..... 103

##### 1.2.2 Tilt Bed Truck ..... 103

##### 1.2.3 Forklift Truck..... 104

##### 1.2.4 Rigging ..... 104

#### Section 2.0

#### Delivery to the Scan Suite..... 105

##### 2.1 Lifting..... 105

##### 2.1.1 Stairway Deliveries..... 105

##### 2.2 Floor Protection ..... 105

##### 2.3 Un-Loading and Un-packing the System ..... 105

#### Section 3.0

#### Dollies ..... 106

##### 3.1 Installations within the United States ..... 106

##### 3.2 Zero Clearance Dollies..... 106

##### 3.3 Tilting Table Dollies..... 106

##### 3.4 Installation Outside of the United States ..... 106

#### Section 4.0

#### Gantry Delivery Considerations ..... 107

##### 4.1 Gantry Shipping State ..... 107

##### 4.2 Door Openings ..... 107

##### 4.3 Elevator Requirements ..... 108

#### Section 5.0

#### Table Delivery Considerations ..... 108

#### Section 6.0

#### Console Delivery Considerations ..... 109

#### Section 7.0

#### Storage Requirements ..... 110

##### 7.1 Short-term Storage (Less than Six Months)..... 110

|                     |   |            |
|---------------------|---|------------|
| 7.2                 | Construction-Site Storage.....                                | 110        |
| <b>Section 8.0</b>  |   |            |
|                     | <b>Extreme Temperature Transportation and Deliveries.....</b> | <b>110</b> |
| <b>Section 9.0</b>  |   |            |
|                     | <b>Site Environmental Considerations.....</b>                 | <b>111</b> |
| 9.1                 | Dust/Dirt Contamination.....                                  | 111        |
| 9.2                 | Chemical Contamination.....                                   | 111        |
| <b>Section 10.0</b> |   |            |
|                     | <b>Handling Requirements .....</b>                            | <b>112</b> |
| 10.1                | Transportation.....   | 112        |
| 10.2                | Handling Requirements .....                                   | 112        |
|                     | 10.2.1 Avoid Dropping .....                                   | 112        |
|                     | 10.2.2 Avoid Shocks and Vibrations .....                      | 112        |
|                     | 10.2.3 Avoid Tipping .....                                    | 112        |
| 10.3                | Inclines and Flat-bed Truck Removal .....                     | 113        |

## Chapter 7

### Power Requirements ..... 115

|                    |  |            |
|--------------------|--|------------|
| <b>Section 1.0</b> |  |            |
|                    | <b>Introduction .....</b>                          | <b>115</b> |
| <b>Section 2.0</b> |  |            |
|                    | <b>System Input Power .....</b>                    | <b>116</b> |
| 2.1                | Power Source Configuration .....                   | 116        |
| 2.2                | Rating .....                                       | 116        |
| 2.3                | Regulation.....                                    | 117        |
| 2.4                | Phase Imbalance .....                              | 117        |
| 2.5                | Sags, Surges & Transients .....                    | 117        |
| 2.6                | Grounding .....                                    | 117        |
| <b>Section 3.0</b> |  |            |
|                    | <b>Recommended Power Distribution System .....</b> | <b>118</b> |
| <b>Section 4.0</b> |  |            |
|                    | <b>Ground System.....</b>                          | <b>120</b> |

## Chapter 8

### Interconnection Data ..... 121

|                    |   |            |
|--------------------|---|------------|
| <b>Section 1.0</b> |   |            |
|                    | <b>Introduction .....</b>                         | <b>121</b> |
| <b>Section 2.0</b> |   |            |
|                    | <b>Component Designators.....</b>                 | <b>121</b> |
| <b>Section 3.0</b> |   |            |
|                    | <b>Interconnect Runs, Wiring and Cables .....</b> | <b>122</b> |
| 3.1                | GE Supplied (Standard Length).....                | 122        |
| 3.2                | GE Supplied (Optional, Long Run) .....            | 128        |

|   |  |            |
|---|--|------------|
| 3.3   | GE Supplied (Console Cables) .....             | 131        |
| 3.4   | GE Supplied (Cables of Options) .....          | 132        |
| 3.5   | Contractor (Customer) Supplied .....           | 134        |
| <b>Section 4.0</b>                            |  |            |
| <b>Contractor Supplied Components .....</b>   |  | <b>136</b> |
| <b>Section 5.0</b>                            |  |            |
| <b>Fuse .....</b>                             |  | <b>137</b> |
| <b>Section 6.0</b>                            |  |            |
| <b>UPS Interconnect.....</b>                  |  | <b>138</b> |
| <b>Section 7.0</b>                            |  |            |
| <b>Typical Customer Supplied Wiring .....</b> |  | <b>139</b> |
| 7.1   | Primary Power Disconnect - U.S.....            | 139        |
| 7.2   | Scan Room Warning Light & Door Interlock ..... | 140        |

# Chapter 1

## Introduction

This direction contains physical and electrical data necessary for planning and preparing a site. Pre-installation work is defined as site preparation for installation of the GE CT scanner. It is the responsibility of the purchaser to arrange and pay for this work. Pre-installation work includes:

- Installation of electrical conduit, junction boxes, ducts, outlets, and line safety switches.
- Installation of interconnection wiring that is AWG stranded copper. The electrical contractor shall ring out and tag all wires at both ends. Color-coded wires are recommended for easier identification. Wires shall be continuous without splices. Ground wires must conform to local codes.
- Any site renovation.
- Alterations and modifications to products not specifically included in the sales contract.

All work must conform to local building and safety codes. Unless specifically mentioned, GE does not provide or install wires, conduits, junction boxes, and ducts as illustrated in this publication.

All CT site plans, preliminary concepts and final working drawings must be reviewed by General Electric Headquarters Architectural Planning prior to construction or approval.

Contact your local GE sales representative for complete information regarding your site-specific room layout.

## Section 1.0

### Site Readiness

Site ready is a requirement that must be achieved to install a CT product. For your convenience, a site ready visit inspection shall be performed at least three (3) days prior to the installation date. The site inspection must conclude with a minimum of a conditional pass status to be ready on the requested installation delivery date. Site ready inspections on the delivery date will not be acceptable unless prior arrangements have been made.

Pre-Installation and Site Ready Tools:

- Floor template
- Pre-Installation check List
- Pre-Installation Block Diagram
- Site Room Layouts
- Power and Grounding Inspection
- Pre-Installation Support

## Section 2.0

# Responsibility of Purchaser

### 2.1 Customer Room Prep Items

The CT air intake is near the bottom of the gantry and draws air in through a filter in the gantry heater assembly. Fine dust as listed below will clog the filter and be deposited throughout the gantry, table, console and PDU electronics. This fine dust cannot be completely removed and can be damaging to electronic components.

For these reasons, the scanner should be the last item installed in your CT suite area.

“Pre-installation” is work necessary to plan and prepare a site for installation of equipment.

Pre-installation work helps the user (customer) avoid:

- Application delay and scheduling
- Surprise siting discoveries
- Installation confusion
- Waste of manpower

The following **MUST** be completed before installation work can begin for a GE CT scanner:

- Completely finished:
    - Wall painted or have final wall covering
    - Ceiling tiles installed and no remaining ceiling work is required
    - Final floors covering installed with no remaining dust causing floor work required
    - All room millwork installed as shown on the site print
    - All plumbing work in the CT suite is completed
    - No construction in or around the scan suite AREA that will produce:
      - \* Concrete dust
      - \* Drywall dust
      - \* Ceiling tile dust
      - \* Wood sawdust or shaving
      - \* Dust tracked into the CT suite area
  - Active Broad Band connection
    - A completed network connection is required for ALL CT installations.
    - A GE network specialist may be required to complete the VPN connection. This may take a week or longer to schedule.
  - Power available to A1, with provision for Lockout/Tagout at the A1 disconnect
- If a UPS is required, a GE A1 breaker\* will be needed to complete this installation. Refer to the electrical section for more details.



#### NOTICE

**SERVICE NOTICE:** An improperly prepared site (i.e., one that is in a state of construction) can result in increased installation time.

**A CT scanner installed in a dirty environment is more prone to contamination, which can result in decreased reliability and increased scanner downtime.**

## 2.2 Purchaser Site Preparation Work

This list below will describe many of the items to consider when planning for a system replacement or designing a room for new equipment

- Determine room dimensions and verify that doorways are large enough for the scanner system.
- Install appropriate conduits and duct work for system cables. If additional components are required in the CT suite, their connection consideration must be determined and completed.
- Install junction boxes of correct size with covers at locations shown in installation plan.
- A1 main disconnect installation
- Install power supply of correct voltage output and adequate KVA rating.
- Install local disconnects, including proper over-current protection.
- Install “steelwork” or other suitable support work for mounting equipment on walls or from ceiling.
- Camera should be on-site at the time of installation.
- Complete all suite and room alterations and modifications.
- Verify that room shielding is adequate for the system being installed
- Review structural requirement - including floor vibration, levelness, and thickness
- Review HVAC requirements including system regulation and patient comfort.
- Review operational clearances to see if your daily used items fit, such as beds and carts.
- Emergency medical equipment should also be considered
- Storage cabinets and sink (if required) must be shown on the site print
- These contractors and others may be required to help confirm that the site meet all installation requirements:
  - Structural Engineer and /or Architect
  - HVAC contractor
  - Electrical contractor
  - Qualified radiological health physicist

The above items can be found in chapter 2 through 9 in this manual.

It is suggested that this work be completed at least three days prior to delivery

## 2.3 Manufacturer’s System Level Siting Requirements

These siting requirements are the minimum that must be met in order to install a new or replacement system.

- Network Communication in place and active
- Meets all scan room regulatory and service requirements
- Meets all minimum scan room structural requirements
- Meets minimum scan room HVAC requirements
- Meets minimum scan room electrical requirements
- Reviewed radiation protection section in the Pre-Installation manual
- All in room items shown on the final GE site print and the final print is on site
- No construction in the scan room or neighboring suite areas

It is suggested that this work be completed at least three days prior to delivery

### 2.3.1 Meeting Site Ready Requirements

The site ready visit will take place at least three days prior to the delivery date. The site ready visit is intended to verify that all of the siting requirement are met and the site is ready for installation.

The site ready visit will result in a report to the project manager indicating one of the following:

**Pass** - All required items are present, completed and the site is ready for installation.

**Conditional Pass** - is issued when 80% of all of the task are completed and all parties agree that the 20% will be completed by the installation delivery date.

If a "Conditional Pass" is granted on the inspection date, the project manager must present conclusive evidence that unfinished tasks are completed and that the site is ready for delivery one business day prior to delivery.

**Fail** - is issued when less than 80% of the task are completed and all parties cannot agree that the remaining work will be completed by the requested installation delivery date. Failed sites will be rescheduled when all items are completed.

### 2.3.2 Quick Installs

Quick installations are described as sites with minimum room improvements required. These include, but are not limited to the following items:

- Existing electrical disconnect device, wire size and grounds meet all of the above requirements.
- Existing structural items including floor thickness meet all of the above requirements
- Existing HVAC capacity and regulation meet all of the above requirements
- Existing CT suite meets all of the above regulatory and minimum size requirements
- Existing facility can accommodate the delivery and meet all of the above delivery requirements

Quick Installs are subject to the following restrictions:

- Quick installs must have a new room print that accurately reflects the rooms to be upgraded.
- New floor anchors must be a minimum of 4" from any existing floor penetrations.

Quick Installs typically involve a weekend de-install and room prep completion, with a next business day delivery and install.

### 2.3.3 "Two-Step" and Upgrade Installs

A "Two-Step" installation is the practice of temporarily installing one CT system in a site with the intention of upgrading the site to a different CT system at a later date.

- For a "two-step" installation to be considered, the room must meet the minimum room requirements for the project being upgraded.
- As with any upgrade installation, "two-steps" are subject to ALL of the siting requirements imposed by the upgrade/final system. This includes the recommended room size as well as electrical, structural and HVAC requirements.
- Two-steps and other upgrades may be done as "Quick" Installs. In this case, all requirements described in Section 2.3.2 (above) also apply.
- It is the customer's responsibility to check that all requirements are met.
- Rooms that do not meet the minimum requirements for the final product must either upgrade (or enlarge) their room, or consider the "Left-Side Limited Access" option.



### 2.3.4 Site Ready Inspection Visit

Must meet all in section 1 and 2 plus these additions reviewed at the site ready visit.

The GE project manager will review the site delivery process with you to determine how to best transfer the equipment from the transportation truck to your room.

This site ready inspection will review and check these items:

#### Delivery information

- Determine delivery route into the scan room
- Determine if tilt dollies or riggers are required
- Determine if elevators, doorways and hallways are adequate for delivery
- Determine if floor protection is required
- Determine if a tilt bed truck is required for ground delivery and ordered.

#### Regulatory Requirement

- Room size meets the minimum requirements
- Site print is present and accurately reflects the room size and layout.
- No grounded walls are present in the regulatory clearance areas
- All regulatory clearances space is met
- Room meets all local codes

Manufacturer Requirements: As listed in section 2 - all requirements are met

Purchaser's Site Preparation Work: As listed in section 1 - all actions are completed

## Section 3.0 Pre-Installation Checklist

| Global Site Readiness Checklist  |                               |   |
|--|-------------------------------|---|
| Customer Name: Customer Name   | PMI Name:                     |   |
| GON Number:  | FSE Name:                     |   |
| Equipment:   | Country / City:               |   |
| Required site assessment milestones  | Date of completion (dd/mm/yy) |   |
| 1) Check site before Delivery/Storage  | Storage Site Ready Date       |   |
| 2) Check site before installation start  | Actual Site Ready Date (SRD)  |   |
| Place an "X" in either Y or N column   |                               |   |
| Site Ready Checks before Delivery/Storage  | Y                             | N |
| Sufficient & secured storage space is planned with the customer.   |                               |   |
| Environmental requirements for storage place are met per GE requirements.  |                               |   |
| All permits, plans and permissions received for rigging and/or delivery.   |                               |   |
| Rooms that will contain equipment, including staging areas if applicable, are dust free. Precautions must be taken to prevent dust from entering rooms containing equipment.   |                               |   |
| Adequate delivery route from truck to final place of installation has been reviewed with all stakeholders, all communications/notifications have occurred, arrangements have been made for special handling (rigging, elevator, fork lift, etc.)   |                               |   |
| All floors along delivery route will support weight of the equipment, temporary reinforcements arranged if needed.   |                               |   |
| Site Ready Checks at Installation  | Y                             | N |
| General Site Planning  |                               |   |
| Room dimensions, including ceiling height, for all Exam, Equipment/Technical & Control rooms meets GE specifications.  |                               |   |
| Ceiling support structure, if indicated on the GE drawing, is in the correct location, at the correct height, levelness and spacing has been measured.   |                               |   |
| Rooms that will contain equipment, including staging areas if applicable, are construction debris free. Precautions must be taken to prevent debris from entering rooms containing equipment.  |                               |   |
| Finished ceiling is installed. If applicable ceiling tiles installed per PMI discretion.   |                               |   |
| Adequate delivery route from truck to final place of installation has been reviewed with all stakeholders, all communications/notifications have occurred, arrangements have been made for special handling (rigging, elevator, fork lift, etc.). All floors along delivery route will support weight of the equipment, temporary reinforcements arranged if needed. |                               |   |
| System power & grounding (PDB/MDP) is available as per GE specifications, installed at point of final connection and ready to use. Lock Out Tag Out is available.  |                               |   |
| Adequate room illumination installed and working.  |                               |   |
| Cable ways (floor/wall/ceiling/Access Flooring) are available for installation of GE cables and are of correct length and diameter. Surface floor duct can be installed at time of system installation.  |                               |   |
| HVAC systems Installed, and the site meets minimum environmental storage requirements.   |                               |   |
| Network outlets installed and computer network available and working.  |                               |   |
| Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications.   |                               |   |

|  |                       |                |  |
|--|-----------------------|----------------|--|
| Customer supplied countertops where GE equipment will be installed are in place.   |                       |                |  |
| <b>Specific</b>  |                       |                |  |
| Doors and windows complete or scheduled to be installed. If applicable, radiation protection (shielding) finished & radioprotection regulatory approval for installation obtained. |                       |                |  |
|  | <b>STATUS OF WORK</b> |                |  |
| <b>GENERAL COMMENTS</b>  |                       |                |  |
| <b>SYSTEM CAN BE DELIVERED</b>   |                       | PMI Name:      |  |
| <b>SITE READY FOR INSTALLATION</b>   |                       | PMI Signature: |  |

## Section 4.0

# Medical Electrical Equipment for EMC

### 4.1 General Scope

The Discovery CT590 RT / Optima CT580 / Discovery RT scanner complies with IEC 60601-1-2: 2004 (for the scanner with the gantry model number of 2374681-8, 2374681-9, 2374681-10 and 2374681-11) and IEC 60601-1-2: 2007 (for the scanner with the gantry model number of 2374681-12, 2374681-13, 2374681-14 and 2374681-15).

The system is suitable to be used in the electromagnetic environment, as per the limits & recommendations described in the tables hereafter:

- Emission Compliance level & limits (see [Table 1-1](#))
- Immunity Compliance level & recommendations to maintain equipment clinical utility (see [Table 1-2](#), [Table 1-3](#), and [Table 1-4](#)).

Note: This system complies with above mentioned EMC standard when used with supplied cables. If different cable lengths are required, contact a qualified GE service representative for advice.

### 4.2 Electromagnetic Emission

**Table 1-1 Emission Declaration**


| System EMC Emissions Guidance & Declaration  |                |  |
|--|----------------|--|
| The system is intended for use in the electromagnetic environment specified below. The customer or the user of RT system should assure that it is used in such an environment. |                |  |
| Emissions Test   | Compliance     | Electromagnetic Environment Guidance   |
| RF emissions<br>CISPR 11   | Group 1        | The system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                          |
| RF emissions<br>CISPR 11   | Class A        |  |
| Harmonic emissions<br>IEC 61000-3-2  | Not applicable | The system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-2  | Not applicable |  |

## 4.3 Electromagnetic Immunity

Table 1-2 Immunity Declaration

| System EMC Immunity Guidance & Declaration   |  |  |   |
|--|--|--|---|
| System is intended for use in the electromagnetic environment specified below. The customer or the user of RT system should assure that it is used in such an environment. |  |  |   |
| Immunity Test  | IEC 60601-1-2 Test Level                                       | Compliance Level   | Electromagnetic Environment Guidance  |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ± 6 kV contact<br>± 8 kV air                                   | ± 6 kV contact<br>± 8 kV air                                   | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Electrical fast transient/burst<br>IEC 61000-4-4   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines | ± 2 kV for power supply lines<br>± 1 kV for input/output lines | Mains power quality should be that of typical commercial or hospital environment.   |
| Surge<br>IEC 61000-4-5   | ± 1 kV line-line<br>± 2 kV line-earth                          | ± 1 kV line-line<br>± 2 kV line-earth                          | Mains power quality should be that of a typical commercial or hospital environment.   |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11   | < 5 % $U_T$<br>(> 95% dip in $U_T$ )<br>for 5 sec              | < 5 % $U_T$<br>(> 95% dip in $U_T$ )<br>for 5 sec              | Mains power quality should be that of a typical commercial or hospital environment. If the user of RT system requires continued operation during power mains interruptions, it is recommended that RTsystem be powered from an uninterruptible power supply or a battery. |
| Powerfrequency (50/60 Hz) magnetic field<br>IEC 61000-4-8  | 3 A/m  | 3 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.   |
| Note : $U_T$ is the a.c. mains voltage prior to application of the test level.   |  |  |   |

Table 1-3 Immunity Declaration con't

| System EMC Immunity Guidance & Declaration  |  |  |  |
|---|--|--|--|
| The system is intended for use in the electromagnetic environment specified below. The customer or the user of RT system should assure that it is used in such an environment.  |  |  |  |
| Immunity Test   | IEC 60601-1-2 Test Level   | Compliance Level   | Electromagnetic Environment Guidance   |
| <p>Conducted RF<br/>IEC 61000-4-6</p> <p>Radiated RF<br/>IEC 61000-4-3<br/>(alternative method: IEC 61000-4-21)</p>   | <p><math>3 \text{ V}_{\text{RMS}}</math><br/>150kHz to 80 MHz</p> <p>3 V/m<br/>80 MHz to 2.5 GHz</p> | <p>3 V<br/>150kHz to 80 MHz</p> <p>3 V/m<br/>80 MHz to 2.5 GHz</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the system, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p><b>Recommended Separation Distance</b></p> $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$ <p>(see Table 1-4)</p> $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$ <p>80 MHz to 800 MHz (see Table 1-4)</p> $d = \left[ \frac{7}{3} \right] \sqrt{P}$ <p>800 MHz to 2.5GHz (see Table 1-4)</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which RT System is used exceeds the applicable RF compliance level above, RT System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the system.</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> <p>Note : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p> |  |  |  |

**Table 1-4 Separation Distances**

| Recommended separation distances between portable and mobile RF communications equipment and the system.  |   |   |   |
|---|---|---|---|
| The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of RT System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment. |   |   |   |
| Rated Maximum Output Power (P) of Transmitter Watts (W)   | Separation distance according to frequency of transmitter |   |   |
|   | 150 kHz to 80 MHz   | 80 MHz to 800MHz                            | 800 MHz to 2.5 GHz                        |
|   | $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$               | $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$ | $d = \left[ \frac{7}{3} \right] \sqrt{P}$ |
|   | Separation Distance meters                                | Separation Distance meters                  | Separation Distance meters                |
| 0.01  | 0.12  | 0.12  | 0.23                                      |
| 0.1   | 0.37  | 0.37  | 0.74                                      |
| 1   | 1.17  | 1.17  | 2.33                                      |
| 10  | 3.69  | 3.69  | 7.38                                      |
| 100   | 11.7  | 11.7  | 23.3                                      |
| For transmitters rated at a maximum output power not listed above, the separation distance can be estimated using the equation in the corresponding column, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.   |   |   |   |
| Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.  |   |   |   |
| Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.  |   |   |   |

### 4.3.1 Limitations Management :

Adhering to the distance separation recommended in [Table 1-4](#), between 150kHz & 2.5GHz, will reduce disturbances recorded at the image level but may not eliminate all disturbances. However, when installed and operated as specified herein, the system will maintain its essential performance by continuing to acquire, display, and store diagnostic quality images safely.

(\*) For example, a 1W mobile phone (800MHz to 2.5GHz carrier frequency) shall be put 2.3 meters apart from the system (in order to avoid image interference risks).

## 4.4 Use Limitation:

### 4.4.1 External components

The use of accessories, transducers, and cables other than those specified may result in degraded ELECTROMATHNETIC COMPATIBILITY of the system.

## **4.5 Installation Requirements & Environment Control :**

In order to minimize interference risks, the following requirements shall apply.

### **4.5.1 Cable Shielding & Grounding**

All interconnect cables to peripheral devices must be shielded and properly grounded. Use of cables not properly shielded and grounded may result in the equipment causing radio frequency interference.

### **4.5.2 This product complies the radiated emission as per CISPR11 Group 1 Class A standard limits**

System is predominantly intended for use, in non-domestic environments, and not directly connected to the Public Mains Network. System is predominantly intended for use (e.g. in hospitals) with a dedicated supply system, and with a X-ray shielded room. In case of using in a domestic environment (e.g. doctor's offices), in order to avoid interferences, it is recommended to use a separated AC power distribution panel & line with a X-ray shielded room.

### **4.5.3 Subsystem & Accessories Power Supply Distribution**

All components, accessories subsystems, systems which are electrically connected to the system, must have all AC power supplied by the same power distribution panel & line.

### **4.5.4 Stacked Components & Equipment**

System should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, System should be observed in order to verify normal operation in the configuration in which it will be used.

### **4.5.5 Low Frequency Magnetic Field**

In case of a digital System, the Gantry (digital detector) shall be apart 1 meter from the generator cabinet, and 1 meter apart from the analog (CRT) monitors. These distance specifications will minimize the low frequency magnetic field interference risk.

### **4.5.6 Static Magnetic Field Limits**

In order to avoid interference on system, static field limits from the surrounding environment are specified.

Static field is specified less than <1 Gauss in Examination room, and in the Control Area.

### **4.5.7 Electrostatic Discharge Environment & Recommendations**

In order to reduce electrostatic discharge interference, install a charge dissipative floor material to avoid electrostatic charge buildup.

The relative humidity shall be at least 30 percent.

The dissipative material shall be connected to the system ground reference, if applicable.



# Chapter 2

## Pre-Installation Overview

Before a system can be installed, all pre-installation requirements must be complete.

- Chapter 4, Section 6.0 and 7.0 Structural Requirements
- Chapter 4, Section 8.0 Radiation Protection
- Broad-band standard
- Site Ready Visit
- [Chapter 4, Environmental Conditions](#) Section 1.0 & 2.0 HVAC Requirements
- [Chapter 7, Power Requirements](#) (Site Power Audit Required)
- [Chapter 8, Interconnection Data](#)

Site-specific items must be verified before the installation can begin.

### Section 1.0 Dust/Dirt Contamination

The systems (consisting of: Console, PDU, Table and Gantry) are highly susceptible to airborne contaminants, especially concrete and drywall dust. Due to the possibility of contamination, these systems should NEVER be installed in a construction site.



**NOTICE** Any site with unfinished floors, walls or ceilings is considered a construction site, and is not suitable for system installation.

### Section 2.0 Chemical Contamination

Wet film processors must never be installed in the same room as the scanner, due to the possibility of chemical contamination of the components. Such chemicals can contribute to increased equipment failures, increased system downtime, and decreased reliability. Film processor equipment installation must meet the manufacturer's requirements (e.g. ventilation specifications) and all applicable national and local codes. Also, consideration's should be given to the location of this equipment and chemical fumes relative to human contact as it relates to locating this equipment and chemicals in the control room.

### Section 3.0 Walls, Ceiling, and Floor

All walls, ceiling, and flooring must be completed before installation can begin.  
System scanners can be installed on a 102 mm (4") concrete floor for GT Table, 110 (4-1/3 in) for High Capacity Table.

## Section 4.0 Broad-band

For information on Broad-band requirements, refer to [Chapter 3, Room Planning, Section 6.0](#).

## Section 5.0 Phone Line (for optional modem)

Two phone lines must be installed at or near the console and be operational prior to installation.

- 1-Analog line (for modem use)
- 1-Voice line

## Section 6.0 Review

The systems use adjustable leveling pads to support the gantry and table. The gantry has four (4) primary leveling pads. Refer to [Table 2-1](#) below, to determine the number of leveling pads, by table style.

**Table 2-1 Number of Leveling Pads and Floor Template P/N, by Table Type**

| TABLE               | # OF PADS | TEMPLATE  |
|---------------------|-----------|---|
| GT1700              | 4         | 5164728   |
| High Capacity Table | 4         | 5164728<br>(part number on Rating Plate: <b>5272966-2</b> ) |


**Note:** This Room Layout Template (5164728) is shipped with the system. It may also be ordered via the Web, from Coakley-Tech.

Using the GE print to establish the room layout, make sure all the operating and service clearances shown on the print are observed. Using the template (see [Table 2-1](#)) shipped with the system, locate the anchor holes. Make sure they clear structural interferences in the floor.

Clean the area. Free the mounting surface of any material that may interfere with the positioning and leveling of the system.

- 1.) Lay out the 2 floor templates.
- 2.) Start with the Table template-align per the GE print.
- 3.) Place the Gantry template over the top of the Table template. Align the scan and table center-lines and secure the templates to the floor. Make sure there are no potential clearance issues.
- 4.) Check the level of the floor (See [Figure 2-1](#)) across the templates.

**Note:** Tiles (or other resilient flooring) around all holes will be cut during the installation process.



**Figure 2-1 Hole Locations**

## FLOOR LEVELNESS SPECIFICATION

- 6 mm (1/4 in.) over 3 m (10 ft.)

This should be measured on the template over the table/gantry area, as shown in [Figure 2-1](#), above.



# Chapter 3

## Room Planning

### Section 1.0 Regulatory Clearances

#### 1.1 Regulations

Federal regulations and national standards that determine the minimum clearances for U.S. installations include:

- 29 CFR 1910 (OSHA)
- NFPA 70E (STANDARD FOR ELECTRICAL SAFETY IN THE WORKPLACE)
- NFPA 101 (LIFE SAFETY CODE)
- AMERICANS WITH DISABILITIES ACT



**NOTICE**

All systems installed within the United States and United States territories, and within U.S. government facilities, regardless of country, must comply with all United States Federal and local regulations. All systems installed outside the United States must comply with either the national, state, or local regulatory clearance requirements for the country in which the installation occurs, or U.S. Federal regulations, whichever is greater.

#### 1.2 Clearance Requirements

A map of clearance requirements necessary for proper operation and servicing of the system is provided in [Figure 3-3, on page 46](#). This is for a standard layout in the suggested room size.



**NOTICE**

The maps and dimensions shown in this manual depict the required clearances for proper equipment operation and service only. The customer/purchaser is responsible for federal, state and/or local codes regarding facility egress and related facility requirements.



**NOTICE**

The use of alternative layouts from the appendix puts severe limitations on space for patient care and work flow. Customer approval of site drawings signifies customer agreement to these limitations.




Figure 3-1 Regulatory Clearance Requirements for System Configuration

## 1.2.1 Minimum Regulatory Workspace Clearances by Major Subsystem

Note the following when referring to the tables below:

- These requirements apply to equipment operating at 600 V or less, where examination, adjustment, servicing, or maintenance is likely to be performed with live parts exposed.
- The customer MUST maintain the required regulatory clearance distances and may NOT use these areas for storage. This applies during normal system operation as well as during service inspection and maintenance.
- Direction of Service Access refers to a direction perpendicular to the surface of the equipment serviced.

| Work Space Requirement                        | Minimum Clear Space | Additional Conditions   |
|---|---------------------|---|
| Direction of Service Access: front of console | 914 mm (36 in.)     | There are no exposed live part hazards with the cover in place. If the console is placed under a counter, the front edge of the console must be even with the vertical edge of the console workspace.<br><br>Note: This component is typically serviced from the front with access to the rear. |
| Service access width: Front of console        | 762 mm (30 in.)     | This is the width of the workspace in front of the equipment. A minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater, is required.  |
| Head clearance                                | 1981.2 mm (78 in.)  | This is the height of the workspace measured from the floor at the front edge of the equipment to the ceiling or overhead obstruction(s).<br><br>A minimum of 1981.2 mm (78 in.) or the height of the equipment, whichever is greater, is required.   |

Table 3-1 Console Subsystem

| Work Space Requirement                         | Minimum Clear Space | Additional Conditions  |
|--|---------------------|--|
| Direction of Service Access (Front of NGPDU)   | 914.4 mm (36 in.)*  | There are no exposed live part hazards with the cover in place.<br><br>This component is typically serviced from the front with access to the rear.<br><br>*If exposed live parts of 151 - 600 volts are present, 1219 mm (48 in.) is required on both sides of the workspace with the operator between.<br><br>*If the opposite wall is grounded and exposed live parts of 151 - 600 volts are present, 1067 mm (42 in.) is required. |
| Service Access Width (Left-Right of workspace) | 762 mm (30 in.)     | This is the width of the working space in front of the equipment. A minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater, is required.   |

Table 3-2 NGPDU Subsystem

| Work Space Requirement | Minimum Clear Space | Additional Conditions  |
|------------------------|---------------------|--|
| Head Clearance         | 1981 mm (78 in.)    | This is the height of the workspace measured from the floor at the front edge of the equipment to the ceiling or overhead obstruction(s). A minimum of 1981 mm (78 in.) or the height of the equipment, whichever is greater, is required. |

**Table 3-2 NGPDU Subsystem**

- For the gantry and table, distances are measured from the enclosure, not the finish covers.

| Work Space Requirement                         | Minimum Clear Space | Additional Conditions   |
|--|---------------------|---|
| Direction of Service Access (All Sides)        | 914 mm (36 in.)     | If exposed live parts of 151 - 600 volts are present, 1219 mm (48 in.) on both sides of workspace with the operator between is required.<br><br>If the opposite wall is grounded and exposed live parts of 151 - 600 volts are present, 1067 mm (42 in.) is required. |
| Service Access Width (Left-Right of workspace) | 762 mm (30 in.)     | This is the width of the working space in front of the equipment. A minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater, is required.  |

**Table 3-3 Gantry Subsystem**

| Work Space Requirement                         | Minimum Clear Space | Additional Conditions   |
|--|---------------------|---|
| Direction of Service Access (Table Head)       | N/A                 |   |
| Direction of Service Access (Table Sides)      | 914 mm (36 in.)*    | * Can reduce to 711 mm (28 in.) provided the local team obtains written and signed approval from the local AHJ (Authority Having Jurisdiction). GE must have the signed document on file.                                   |
| Front Cover Removal (Table Foot)               | 686 mm (27 in.)     | 457 mm (18 in.) minimum for Front Gantry Cover removal. Refer to the appendix for alternate front cover removal options.  |
| Service Access Width (Left-Right of workspace) | 762 mm (30 in.)     | Refer to the width of working space in front of the equipment. 762 mm (30 in.) minimum or the equipment width, whichever is greater.  |
| Head Clearance                                 | 1981 mm (78 in.)    | Refers to the height of the workspace measured from the floor at the front edge of the equipment to the ceiling or any overhead obstructions.<br><br>1981 mm (78 in.) or the height of the equipment, whichever is greater. |

**Table 3-4 Table Subsystem**



| Work Space Requirement                              | Minimum Clear Space | Additional Conditions   |
|---|---------------------|---|
| Direction of Service Access (Front of UPS)          | 914.4 mm (36 in.)*  | <p>There are no exposed live part hazards with the cover in place.</p> <p>This component is typically serviced from the front with access to the rear.</p> <p>*If exposed live parts of 151 - 600 volts are present, 1219 mm (48 in.) is required on both sides of the workspace with the operator between.</p> <p>*If the opposite wall is grounded and exposed live parts of 151 - 600 volts are present, 1067 mm (42 in.) is required.</p> |
| Service Access Width (Right side and length of UPS) | 762 mm (30 in.)     | This is the width of the working space in front of the equipment. A minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater, is required.  |
| Head Clearance                                      | 1981 mm (78 in.)    | This is the height of the workspace measured from the floor at the front edge of the equipment to the ceiling or overhead obstruction(s). A minimum of 1981 mm (78 in.) or the height of the equipment, whichever is greater, is required.  |

Table 3-5 UPS Subsystem

| Work Space Requirement  | Minimum Clear Space | Additional Conditions   |
|---|---------------------|---|
| Direction of Service Access (Front of A1 Disconnect)          | 914.4 mm (36 in.)*  | <p>There are no exposed live part hazards with the cover in place.</p> <p>This component is typically serviced from the front with access to the rear.</p> <p>*If exposed live parts of 151 - 600 volts are present, 1219 mm (48 in.) is required on both sides of the workspace with the operator between.</p> <p>*If the opposite wall is grounded and exposed live parts of 151 - 600 volts are present, 1067 mm (42 in.) is required.</p> |
| Service Access Width (Right side and length of A1 Disconnect) | 762 mm (30 in.)     | This is the width of the working space in front of the equipment. A minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater, is required.  |
| Head Clearance  | 1981 mm (78 in.)    | This is the height of the workspace measured from the floor at the front edge of the equipment to the ceiling or overhead obstruction(s). A minimum of 1981 mm (78 in.) or the height of the equipment, whichever is greater, is required.  |

Table 3-6 A1 Disconnect Subsystem



## 1.4 NEC Conduit and Duct Fill Rate

Full operation, service, and safety of the system requires the maintenance of sufficient regulatory and service clearances around equipment.

Cable length is an important consideration in room layout. The system ships with standard (short) length cables, with a set of longer cables available as an option. Refer to the electrical page of your GE site print for the specific requirements of your site. The following rules govern cable usage for the system:

- When possible, use the rear cable cover assembly to let cables enter the gantry from the rear.
- Do not cut or otherwise shorten long cables.
- Do not store excess cable length behind the operator console, gantry, or PDU.
- Store excess cable in wall or floor ducts, if desired, provided that sufficient space exists. Refer to NEC code to determine cable fill rates for conduits and ducts.
- All installed systems shall comply with NEC 70-E Electrical Regulations governing conduit or duct fill.

## 1.5 Terms and Definitions

### CLEARANCES

Clearances are the clear space or distance between or around objects and equipment, governed by all applicable safety, service, and regulatory requirements and representing the lowest margin of freedom permissible for equipment siting.

### DIMENSIONS

Dimensions are the length, width, depth, and height of equipment.

### EGRESS



#### NOTICE

The maps and dimensions shown in this manual depict the required clearances for proper equipment operation and service only. The customer/purchaser is responsible for federal, state and/or local codes regarding facility egress and related facility requirements.

### PRE-INSTALLATION ESCALATION

Pre-installation escalation is the process used to consult CT Engineering, the Design Center, or EHS to resolve pre-installation issues related to siting concerns and requirements.

### GROUNDING WALL

A grounded wall is any wall with electrical conductivity to earth. Conductive materials generally found in walls include masonry, concrete, and tile. Treat as grounded additional elements commonly found in walls, including but not limited to:

- Medical gas ports and plates
- Metal door and window frames
- Water sources and metallic sink structures
- Metallic wall mounted cabinets
- A1 disconnect panel
- Equipment Emergency Off panels
- Industrial equipment such as air conditioners and vents
- Expansion joints
- Surface raceway
- Exposed wall conduits
- Floor outlet boxes
- Floor HVAC boxes
- Floor medical gas

Common wall components NOT constituting grounded elements include:

- Standard wall outlet
- Light switches
- Telephones
- Communication wall jacks
- Ceiling tile grids

## HEAD CLEARANCE

*Head clearance* represents the height dimension of the workspace, measured from the floor at the front edge of the equipment to the ceiling or any overhead obstructions. It requires a minimum of 1981 mm (78 in.), or the height of equipment, whichever is greater.

## MINIMUM

Minimum indicates the lowest limit permitted by law or other authority.

## SERVICE ACCESS WIDTH

*Service access width* refers to the width of the working space in front of the equipment, and requires a minimum of 762 mm (30 in.) or the width of the equipment, whichever is greater.

## WORK SPACE


*The workspace* represents a three-dimensional box of space required for safe inspection or service of energized equipment. It consists of depth, width, and height, with the depth dimension measured perpendicular to the direction of access. U.S. regulation requires a minimum depth of 914 mm (36 in.). Additional conditions can increase the minimum requirement. For example, FCT defines workspace as the envelope of the component superstructure, measured for the NGPDU with the front panel removed, and measured for the gantry and table with the external covers removed.

## Section 2.0 Service Clearances Requirements

### 2.1 Measuring Service Clearances

The service envelopes defined in [Figure 3-3](#) provide enough space to safely allow CT System servicing. System Servicing requires:

- Sufficient space to remove the covers from the system.
- One service engineer to be able to accomplish all service component replacement tasks without needing special tools or equipment.
- ALL room layouts to provide service space and access around the table to the gantry right side. This is needed for replacement procedures that require components that ship in large boxes, such as the tube, detector, and HV tank.



#### Minimum Clearances

| Item                            | Dimension  |
|---------------------------------|--|
| 1 GT1700<br>High Capacity Table | 6096 mm (240 in.)<br>6891 mm (271 in.)<br>(travel distance of table) |
| 2                               | 2548 (100 in.)<br>(cover with dolly)                                 |
| 3                               | 3234 mm (127 in.)  |
| 4                               | 1617 mm (64 in.)   |
| 5                               | 2545 mm (100 in.)  |
| 6                               | 914 mm (36 in.)<br>minimum   |
| 7                               | 356 mm (14 in.)  |
| 8                               | 711 mm (28 in.)  |
| 9                               | 914 mm (36 in.)  |
| 10                              | 914 mm (36 in.)  |
| 11                              | 2794 mm (110 in.)  |
| 12                              | 851 mm (34 in.)  |
| 13                              | 680 mm (26.7 in.)  |

#### CT Components

| Item | Description   |
|------|---|
| A    | Gantry with covers installed                                    |
| B    | Table cradle footprint, coverage as extended in both directions |
| C    | Front gantry cover removed with dolly                           |
| D    | Back gantry cover removed with dolly                            |

Figure 3-3 Minimum Service Clearances

## 2.2 Special Service Clearances

Note: When calculating service clearances, refer to [Figure 3-3](#) for all service clearance needs.

### 2.2.1 Gantry Service Clearance

Specification for Boom Assembly clearance are defined in [Figure 3-4](#). The boom assembly is used during tube and detector replacement. The minimum ceiling height within the clearance radius is 2286 mm (90 in).




Figure 3-4 Boom Assembly Clearance

### 2.2.2 Cover Removal

- Gantry front cover removal requires the use of the Tilting Cover Dollies and a minimum clearance space of 3234 mm (127 in.) to maneuver the cover. The dollies allow the service engineer to separate the cover from the gantry, tilt it 90 degrees, roll it to the foot end of the table, and then tilt it an additional 90 degrees, so that it is upside-down relative to its normal system-mounted condition. After removal, the service engineer must then move the gantry front cover to a position that satisfies the minimum regulatory clearances.
- The gantry rear cover, with service dollies installed, requires a clearance width of 2545 mm (100 in.) and a depth of 914 mm (36 in.) for removal. Sufficient space to allow the service engineer to move the cover either straight back or to one side of the table to satisfy the minimum service clearances shown in [Figure 3-3](#) must be maintained. The rear cover with dollies cannot extend past the allowable clearance space within the room. If the system is not sited straight (it is positioned diagonally), full service space is still required. The PMI and customer should discuss this consideration and make the necessary provisions.
- The scan room must offer sufficient space to allow adequate egress during service operations that require both front and rear cover removal. If the customer and PMI have any concern that site will not provide adequate space for egress under these conditions, they should discuss

these requirements and make the necessary provisions to accommodate this event.

- A single service engineer can safely perform servicing of the table. Ensure sufficient clear space to maintain egress clearances with the table covers or cradle removed.
- A single service engineer can safely perform servicing of the system. Ensure sufficient clear space to maintain egress clearances with covers or cradle removed.
- A tube change box is 737 mm x 1524 mm (29 in. x 5 ft.), with the handles extended. The box rolls like a wheelbarrow and must have access to the right side of the gantry. It is the PMI's responsibility to demonstrate that the tube change box can be positioned in the tube change area next to the gantry and that the front and rear covers can be removed.

### 2.2.3 Power Distribution Unit (NGPDU)

When positioning the Power Distribution Unit (NGPDU), consider regulatory compliance, as defined in [Section 1.0 - Regulatory Clearances](#). See [Table 3-2](#) in that section.

### 2.2.4 Console

The operator console does not present an exposed live parts hazard. However, the site shall maintain a working space at all times with a minimum depth of 1219 mm (48 in.), extending the full width of the operator console for service activity.

The console is on wheels. As some service activities require access to the rear of the console, be sure to maintain sufficient space for moving the console to allow rear service access.

See [Figure 3-7](#) for a typical control room layout.

### 2.2.5 Storage Cabinet

GE provides a storage cabinet for storing all supplied service equipment (see [Table 3-7](#)). Situate this storage cabinet within the scan room suite area to allow easy service access. The dimensions of the cabinet measure 457 mm D x 914 mm W x 1067 mm H; ~ 41 kg (18 in. D x 36 in. W x 42 in. H; ~90 lbs).

Note: A storage cabinet is provided as option (B77292CA).

| ITEM                      | SIZE   | WEIGHT (TOTAL) |       |
|---------------------------|--|----------------|-------|
| QA Phantom (water filled) | 20 cm x 15 cm<br>(7.9 in. x 5.9 in.)   | 5.5 kg         | 12 lb |
| 35CM Phantom              | 35 cm x 7 cm<br>(13.8 in. x 2.8 in.)   | 8.2 kg         | 18 lb |
| 48CM Phantom              | 48 x 7 cm<br>(18.9" x 2.8)   | 11.4 kg        | 25lb  |
| Phantom Holder            | 25 cm x 25 cm<br>(9.8 in. x 9.8 in.)   | 3.6 kg         | 8 lb  |
| FE Box (Purple)           | 30 cm x 38 cm x 30 cm<br>(11.8 in. x 15 in. x 11.8 in.)                        | 6.8 kg         | 15 lb |
| Rear Cover Dollies        | 158 cm x 82 cm<br>(62.2 in. x 32.3 in.)  | 11.4 kg        | 25 lb |
| Front Cover Dollies       | 85 cm x 20 cm and 85 cm x 15 cm<br>(33.5 in. x 7.9 in. and 33.5 in. x 5.9 in.) | 15.9 kg        | 35 lb |

**Table 3-7 Equipment Stored in Storage Cabinet**



| ITEM                      | SIZE  | WEIGHT (TOTAL) |       |
|---------------------------|---|----------------|-------|
| Install Support Kit (box) | 30 cm x 30 cm x 38 cm<br>(11.8 in. x 11.8 in. x 15 in.)                     | 9.1 kg         | 20 lb |
| Tube Hoist Kit            | 77 cm x 8 cm and 38 cm x 15 cm<br>(30.3 in. x 3.1 in. and 15 in. x 5.9 in.) | 9.1 kg         | 20 lb |
| Balance Weight Kit        |   | 33 kg          | 73 lb |

**Table 3-7 Equipment Stored in Storage Cabinet (Continued)**

## Section 3.0 Room Size

### 3.1 Room Dimensions

| System Configuration | Suggested Scan Room Size                           | Typical scan Room Size                              | Minimum Scan Room Size                                    | Recommended Ctrl Room Size                        |
|----------------------|--|---|---|---|
| GT1700 mm Table      | 4420 mm x 6706 mm<br>(14 ft. 6 in. x 22 ft.)       | 4267 mm x 6096 mm <sup>1</sup><br>(14 ft. x 20 ft.) | 3708 mm x 6096 mm <sup>2</sup><br>(12 ft. 2 in. x 20 ft.) | 1700 mm x 4420 mm<br>(5 ft. 7 in. x 14 ft. 6 in.) |
| High Capacity Table  | 4420 mm x 6915 mm<br>(14 ft. 6 in. x 22 ft. 8 in.) | 4267 mm x 6706mm<br>(14 ft. x 22 ft.)               | 3708 mm x 6706 mm<br>(12 ft. 2 in. x 22 ft.)              | 1700 mm x 4420 mm<br>(5 ft. 7 in. x 14 ft. 6 in.) |

**Table 3-8 Suggested Room Size Dimensions**

<sup>1</sup>All service/regulatory requirements apply.

<sup>2</sup>All service requirements apply, with the addition of no energized left-side service.

#### 3.1.1 Suggested Room Size

The suggested room configuration offers the most flexibility for future upgrades. It provides both ample workspace and space to add millwork and still meet all regulatory requirements. When local regulations require a sink in the scan room, this room size also provides sufficient space for a sink. This room size accommodates the needs of larger hospitals and medical teaching facilities, where patients may require transportation to the scan area in beds, gurneys, and larger wheelchairs and where they may require the assistance of larger medical care teams. Likewise, this room offers adequate access for crash carts and other emergency medical equipment on both sides of the table.

The suggested size supports all service activities, including tube change, and accommodates all future two-step installations.

#### 3.1.2 Typical Room Size

The typical room configuration represents what is most commonly found at installation sites. It offers adequate workspace and may provide adequate space for a sink, but allows only LIMITED space to add millwork and still meet all regulatory requirements. The size of this room accommodates the needs of larger clinics and medium-sized hospitals, where patients may require transportation into the scan area using gurneys and wheelchairs, and where they may require the assistance of small medical care teams. It provides access for crash carts and other emergency medical equipment on only one side of the table.

The typical room size supports all service activities, including tube change, and may offer compatibility with SOME future upgrades and two-step installations.

#### 3.1.3 Minimum Room Size

The minimum room configuration represents the smallest functionally acceptable space for this product and represents the type of room often found at doctor's offices and smaller clinics and outpatient facilities. Due to its limited size, and to functional and regulatory requirements, this room usually provides only LIMITED workspace, and leaves to NO space to add in-room millwork and sinks and still meet the necessary regulatory and service requirements. This room can accommodate the transportation of patients into the scan area using wheelchairs, and provides access for crash carts and other emergency medical equipment on only one side of the table. Sites

considering a minimum room size may not have been designed with the structural requirements necessary to support the system and consequently may require upgrading prior to installation.

Customers considering a minimum room size should discuss their workspace requirements and future upgrade plans with their PMI, as the size and layout of these rooms often eliminates them from any future upgrade considerations and offers NO compatibility with future two-step installations.

If using the square meters (square footage) to determine regulatory compliance, please note that the front and rear cover clearances are wider than the regulatory clearance along the table length, and that the cover park position is behind the table in the home position.

Note: Sites must provide sufficient space to allow the removal of the rear cover, which is on wheels, from behind the gantry during service operations.



**CAUTION**

**Operational Caution:** In a minimum room layout, the customer should consider workflow, customer access for patient care, and critical-care operations space requirements. Additionally, this room provides only limited equipment access on the gantry left side when loading patients or when positioning patient equipment in the room between the gantry and the wall.

## 3.2 Suggested and Typical Room Layouts

Figure 3-5 and Figure 3-6 show the suggested and typical room layouts, both with and without in-room cabinets. You need to know the locations for medical gas, surface duct work, or other items that make a grounded wall.

Note: Your room layout may meet the suggested or typical room requirements but appear different than Figure 3-5 or Figure 3-6. Your salesperson can provide a detailed room layout for your site.




Figure 3-5 Room Layout with cabinets





Figure 3-6 Room Layout without cabinets

### 3.3 Minimum Room Layouts


**Note:** If the room size does not meet the typical room layouts, but meet the minimum room layouts, you don't have to submit the concession.

**Room A** -Less than 711 mm (28 in.), but greater than 356 mm (14 in.), measured from the covers to the left sidewall. In this configuration service, egress and workspace are compromised around the gantry's left side.



**Room B** - Less than 914 mm (36 in.) but greater than 711 mm (28 in.) measured from the covers to the left sidewall.

In this configuration service, egress, and workspace are acceptable around the gantry.



## 3.4 Control Room Considerations

- The control room must provide a suitable operating environment for the console electronic, and operator working comfort.
- The operator console must remain in the same configuration as shipped. Do not dismantle the console, or remove or rearrange its components except the peripherals.
- TIO should be placed side by side with FWS table considering cable length. Refer to [Chapter 8 - Interconnection Data Section 3.3](#) for cables information. The cables are shipped with TIO Console.
- If operationally possible, the monitor desktop and user desktop components may be removed and placed on a counter-top, providing the cable lengths shipped are not altered or changed. The console cabinet then can be remote mounted, provided the cooling and venting requirements are met.
- A suitable work area, which is within reach of the operator's console, should be provided for placement of the injector control. Injector controls differ in dimensions depending on the brand selected.
- A PACS, workstation, image printer, or filming device are often placed in the console control room area, and sometimes may be directly linked to the console.
- Additional components although linked via network or ethernet cable, are not powered from the CT console.
- Additional room power and network connection should be considered when reviewing the console work space.

### 3.4.1 Typical Control Room Layout




Figure 3-7 Typical Control Room Layout

### 3.5 CT Simulation Laser Alignment Systems

CT simulation laser alignment systems are used to reproduce the patient position for simulation and subsequent treatment on a linear accelerator. These lasers can be wall mounted or can be free standing. Work with your installation planner to determine the correct installation solution for your site.

Laser Options: [www.gammex.com](http://www.gammex.com) and [www.lap-laser.com](http://www.lap-laser.com)

Custom freestanding mounting structures are available for rooms that have no suitable wall or ceiling mounting space. Stanchion systems can provide flexibility in room design for placing alignment devices in rooms that may not be suitable for wall-mounted lasers.




Figure 3-8 Typical Room Layout with CT Simulation Laser Alignment System



## Section 4.0

# System Component Dimensions

### 4.1 Minimum Operating Clearances

The sections in this chapter provide the minimum dimension and operating clearance information for each category of components listed. Be sure that the site conforms to each of these specifications.

#### 4.1.1 Ceiling Pedestal Mount Installation

The distance from the floor to the lowest point of the ceiling pedestal mount for the Injector or Monitor CANNOT measure LESS than 2286 mm (90 in.). Refer to the installation guides of those components for the length of the mounting post.



#### NOTICE

Failure to maintain a distance of at least 2286 mm (90 in.) from the floor to the lowest point of the Injector or Monitor ceiling pedestal mount may pose a safety hazard. For installations with a finished ceiling height that is less than suggested, consideration should be given to utilizing floor mounted components, or attaching the mounting plate in the overhead (for example, above dropped ceiling tiles).

#### 4.1.2 Injector Control Installation

Minimum dimensions and clearances include the following requirements for the injector control:

- Provision of a suitable work area for placement of the injector control, within reach of the operator console.
- Wall mounted, ceiling mounted, and pedestal units require routing of cables from the gantry area to the operator console area. The supplied cable measures 15.2 m (50 ft).
- Injectors require an AC power source that is powered from the console. The IEC power cord is supplied with the injector.
- Available mounts come in several different lengths and configurations. Refer to the injector documentation for detailed installation instructions.

### 4.1.3 System Operational Clearances

The clearances listed in [Table 3-9](#) govern system operation; be sure that the site maintains each of these clearances.

| SYSTEM OPERATION   | MM      | INCHES  |
|--|---------|---------|
| Ceiling Pedestal mount (optional)<br>(lowest point to floor injector or monitor) | 2286 mm | 90 in.  |
| Finished ceiling to floor (recommended)  | 2743 mm | 108 in. |
| Finished ceiling to floor (minimum)  | 2438 mm | 96 in.  |
| Table extension head end w/extender to<br>obstruction                            | 152 mm  | 6 in.   |
| Back of Console to wall  | 96 mm   | 4 in.   |
| Back of PDU to wall  | 152 mm  | 6 in.   |

**Table 3-9 Minimum Dimensions and Clearances**

## 4.2 Component Dimensions

### 4.2.1 System Dimensions




Figure 3-9 Dimensions without Covers




Figure 3-10 Gantry Front Cover with Service Dolly Dimensions

| DESCRIPTION                | WIDTH |       | DEPTH |      | HEIGHT |      |
|----------------------------|-------|-------|-------|------|--------|------|
|                            | MM    | INCH  | MM    | INCH | MM     | INCH |
| Remote Color Monitor (LCD) | 413   | 16.25 | 203   | 8    | 406    | 16   |
| Color printer              | 584   | 23    | 457   | 18   | 178    | 7    |

Table 3-10 Dimensions of Accessories

## 4.2.2 Table and Gantry Dimensions




Figure 3-11 Table and Gantry (Side View)




Figure 3-12 Gantry shown tilted +30° (top) and -30° (bottom) - Table Options

### 4.2.3 Power Distribution Unit Dimensions




Figure 3-13 Power Distribution Unit (NGPDU)




Figure 3-14 Power Distribution Unit (NGPDU)

## 4.2.4 Operator Console Dimensions




Figure 3-15 True-In-One Operator's Console







Figure 3-16 NIO16 Operator's Console

## 4.2.5 Freedom Workspace Table



NOTES:  
DIMENSIONS ARE IN MILLIMETERS (INCHES).

Figure 3-17 Freedom Workspace Table - 5168666-2 (Wide FWS for console)

## 4.2.6 Optima Table




Figure 3-18 Console with Optima Table

## 4.2.7 Smart Workspace Desk




Figure 3-19 Smart Workspace Desk (SWS) (part 5449758-2)

## Section 5.0 Structural Requirements

### 5.1 Table and Gantry Mounting Requirements



**WARNING** **POTENTIAL FOR PATIENT INJURY.**  
**IMPROPERLY SECURED TABLE MAY TIP, DISLODGING PATIENT.**  
**PROPER ANCHORING IS KEY TO MAINTAINING PATIENT SAFETY DURING SYSTEM OPERATION.**



**NOTICE** It is the purchaser's responsibility to provide an approved support structure and mounting method for all floor types other than those listed. General Electric is not responsible for any failure of the support structure or method of anchoring, including seismic requirements.  
**GE is not responsible for methods other than those listed.**

Table and gantry mounting dimensions are shown in [Figure 5-4](#), [Figure 5-5](#), [Figure 5-6](#). Refer to [Chapter 5](#) for additional details of floor loadings, component weights, and Gantry and Table installation and anchoring.

Anchor gantry and table to floor by a means that will maintain their relative alignment and meet applicable building and other local codes, including seismic structural mounting requirements.

Floor structure must be capable of withstanding the occupied weight of table and gantry, and the individual contact area loading of these components. Refer to [Section 4.2](#) for each of the three (3) major components of the system.

Support areas of the patient table and gantry must rest on solid concrete, not resilient tile or carpeting which will slowly yield over a period of time and disturb alignment of table to gantry.

**Factors that could cause misalignment between gantry and table due to floor sag should be considered. The cradle can potentially carry a 227 kg (500 lb) patient or 295 kg (650lb) for High capacity table. Center of gravity changes as cradle cantilevers.**

Take into consideration all other moving weights such as gurneys or personal equipment. Refer to [Chapter 5](#) for gantry and table mounting details.

No part of floor surface within table and gantry, nor the two interface areas between table and gantry, should be higher than the support area for table and gantry.

### 5.2 Floor Anchors

Provided floor anchors are designed for use ONLY on concrete floors that meet the 102 mm (4 in.) concrete floor requirements. At the customer's expense, all other anchoring methods (on floor types other than the 102 mm (4 in.) concrete minimum) must be determined by their structural contractor to meet the stated GE minimum load requirements. The customer's contractor is responsible for the installation of all anchors other than those shipped with the system.

## 5.3 Floor Strength

Concrete floors must have a minimum strength of  $f'_c = 2000$  psi ( $1.4 \times 10^7$  Pa) for mounting floor anchors. It is the responsibility of each customer to have appropriate tests performed to determine and measure concrete strength.

Note: If installing the GE LS scanner on a floor type thinner than a 102 mm (4 in.) concrete floor, the purchaser, at their expense, shall provide acceptable anchoring and mounting methods that meet all structural specifications provided in sections 5.1 through 5.5 of this Pre-Installation Manual.

## 5.4 Floor Levelness

The CT Room floor levelness requirement is important for accurate patient positioning. Floor levelness in the Scan Room must not be greater than 6 mm (1/4 in.) between depression and high spots over any 3048 mm (120 in.) distance within the area of the gantry/table template (see the envelope shown in [Figure 2-1](#), on page 39).

Note: The floor must meet levelness specification to properly align the table/gantry. Minimum gantry height at this specification is 15 mm (1/2 in.) to prevent cable crushing.

Table level may not be achievable if the overall floor levelness is greater than the specification. The overall floor level must be 0 in. to use under gantry cable entrance. The minimum gantry height is 20 mm (3/4 in.) with this option to prevent cable crushing.

The use of floor shims is not suitable to achieve floor levelness. It is recommended that the concrete be leveled to meet this requirement.

## 5.5 Floor Vibration

The CT equipment may be sensitive to vibration in the frequency range of 0.5 to 20 Hz depending on the amplitude of the vibration. It is the customers responsibility to contract a vibration consultant or qualified engineer to implement design modifications to meet the specific limits, However, it is ultimately the customer/architect/engineer responsibility to design the site solution.

### 5.5.1 Steady State Vibration

The maximum steady state vibration transmitted through the floor should not exceed  $10^{-3}$  m/s<sup>2</sup> rms maximum single frequency above ambient baseline from 0.5 to 80 Hz (measured in any 1 hour during a normal operating period).

### 5.5.2 Transient Vibration

The behavioral characteristics must be such that any measurable transient disturbance must also be minimized to less than 0.01 m/s<sup>2</sup> peak-to-peak.

### 5.5.3 Equipment Location

To minimize the interference, the CT equipment should be placed on a solid floor, located as far as possible from the following vibration sources:

- Parking lots
- Roadways
- Subways
- Trains
- Hallways
- Elevators
- Heliports
- Hospital power plants containing pumps, motors, air handling equipment and air conditioning units

## 5.6 Walls: Scan Window

The recommended patient viewing window dimensions are 1219 mm wide x 1067 mm high (48 in. x 42 in.). The location of the window is dependent on the position of Operator Workspace position. Consult [Section 7.0](#) of this chapter and a **qualified radiological health physicist** for radiation protection requirement of glass (lead content and thickness).

Note: The operator at the Operator Workspace must be able to view the patient during a scan.

## Section 6.0 Network Connections

Broad-band is considered the standard network connection for the systems. (A dial-up modem is optional.) A 1000 baseT high-speed network is desired, with 100 baseT network service acceptable. Broad-band connections should use one of the following Category 5 patch cables:

| CAT Num | GE Part Num | Length |
|---------|-------------|--------|
| K9000WB | 2215028-10  | 20 m   |
| K9000KP | 2215028-5   | 10 m   |
| K9000JR | 2215028-4   | 5 m    |
| K9000WA | 2215028-9   | 3 m    |

The CT system is connected to the network through the Console.

- A patch cable (not to exceed 3.05 m (10 ft.)) should be provided by the customer, and it is used to connect the console to a wall box. (See Notes on [Figure 8-4](#))
- Some customer-site units may require cable duct-work or conduit to route connecting network cables to the workstation, camera and console.
- The run from the hospital switch to the CT wall outlet must not exceed 88 m (290 ft.). Bandwidth performance is degraded when the length reaches 91 m (300 ft.) or greater.
- For the optional modem: **Two phone lines should be provided by the facility.** One line is for use with a modem and must be an analog line. The second line is a voice only line.




Figure 3-20 True-In-One Console Host Computer Rear






Figure 3-21 NIO16 Console with Z800 Host Computer Rear




Figure 3-22 NIO16 Console with Z840 Host Computer Rear

## 6.1 US Broad Band Process Overview

The United States network connectivity requirement for this product is broad-band. The US process relies on the GE PMI to select a Customer Champion and identify an IT contact for the site. Together, those individuals then complete a site assessment to gauge what tasks are needed to fulfill the connection.

Anyone can contact the GE Connectivity team at 800.321.7937, Option #3, with questions.

## 6.2 Customer Broad-Band Responsibilities

Provide GE Installation Project Manager with an accurate site address, telephone number, contact name, and email address for the:

- Customer Champion
  - Coordinate VPN activities between Radiology/Cardiology and the Information Technology (IT) departments
  - Act as a focal point in assuring site broad-band infrastructure meets GE requirements for connection as determined by a mutual assessment with the GE Connectivity team.
- IT Contact
  - Complete an equipment assessment with GE Connectivity team to determine site readiness for broad-band
- Contact your GE PMI, for the name of the zone broad-band specialist to:
  - Work with the Customer Champion to complete any identified infrastructure changes
  - Provide IP addresses for new CT equipment
  - Provide a VPN compatible appliance that will support the IPSec tunneling protocol and 3DES data encryption
  - Utilize an Internet Service Provider that supports static routing

## Section 7.0 Radiation Protection



**NOTICE** Scanner-room shielding requirements should be reviewed by a qualified radiological health physicist taking into consideration:

- Scatter radiation levels within the scanning room (see [Figure 3-23](#))
- Equipment placement
- Weekly projected work-loads (# patient/day technique (kvp\*ma))
- Materials used for construction of walls, floors, ceiling, doors, and windows
- Access to surrounding scan room areas
- Equipment in surrounding scan room areas (e.g., film developer, film storage)

[Figure 3-23](#) depicts measurable radiation levels within the scanning room while scanning a 32 cm CTDI phantom (body) and a 20 cm water phantom (head) with the technique shown.

Note:  
Actual  
measurements  
can vary.

All measurements have an accuracy of  $\pm 20\%$  because of measurement equipment, technique, and system-to-system variation.

Use the correction factors shown in [Table 3-11](#) to adjust exposure levels to the usual scan technique at your site.

| CHANGED PARAMETER                         | MULTIPLICATION FACTOR |
|---|-----------------------|
| mAs                                       | new mAs/100           |
| 80 kV                                     | 0.21                  |
| 120 kV                                    | 0.71                  |
| 140 kV                                    | 1.0                   |
| 4 x 3.75mm images                         | 0.82                  |
| 16 x 0.625 LD<br>8 x 1.25 LD<br>Fluro 5mm | 0.59                  |
| 4 x 1.25 LD<br>5mm (1i)<br>Fluro 2.5 mm   | 0.40                  |
| 1 x 1.25mm images                         | 0.20                  |
| 4 x 0.625 LD<br>1 x 1.25                  | 0.10                  |

**Table 3-11 Shielding Requirements Scaling**



**NOTICE** This publication uses  $\mu\text{Gy}$  (micrograys) to measure radiation levels. The conversion factor from mR to  $\mu\text{Gy}$  (micrograys) is:  $1 \text{ mR} = 8.76 \mu\text{Gy}$




Figure 3-23 Typical Scatter Survey (Head Filter)




Figure 3-24 Typical Scatter Survey (Body Filter)

# Chapter 4

## Environmental Conditions

Ratings and duty cycles of CT subsystems apply if site environment meets the standards of this section. Maintain environmental conditions listed below at all times – including, for example, overnight, weekends and holidays. Shut down the CT system if air conditioning is not working. When system is shut down for major repair, air conditioning may be shut down also.



**NOTICE** Silver, copper, gold films used in the CT system are especially sensitive to the presence of sulfide, chloride and nitrate contaminants. Sulfur is the most damaging. If high levels of contaminants exist site may want to consider appropriate air filtration systems.

### Section 1.0

#### Temperature and Humidity Specifications

Environmental specifications apply to the table, gantry, power distribution unit, and operating console.



**NOTICE** System Operation and Image Quality may be affected, if environmental specifications are exceeded.

#### 1.1 Temperature (Scan Room & Control Room)

|  |                               |
|--|-------------------------------|
| Maximum allowable ambient room temperature:                | 26° C (79° F)                 |
| Recommended ambient room temperature:                      | 22° C (72° F)                 |
| Minimum allowable ambient room temperature:                | 18° C (64° F)                 |
| Maximum allowable ambient room temperature rate of change: | 3° C per hour (5° F per hour) |

**Note:** Any cooling equipment cycle control range must be taken into account, such that the maximum and minimum ambient room temperatures shown above are not exceeded, during room thermal cycling. For example, if the HVAC is capable of  $\pm 2^{\circ}$  C control, then the limits would 20° C - 24° C, to maintain absolute limits.

## 1.2 Humidity (Scan Room & Control Room)

|   |             |
|---|-------------|
| Maximum allowable non-condensing relative humidity: | 60%         |
| Minimum allowable non-condensing relative humidity: | 30%         |
| Maximum allowable relative humidity rate of change: | 5% per hour |

## 1.3 Other Guideline

- To help determine the hospital room environmental conditions, a temperature and humidity recorder may be temporarily installed close to where the system will be installed. Note readings before and after installation, to verify the true temperature and humidity in your environment.
- Consider heating, ventilating, air conditioning (HVAC) needs and redundancy. An air conditioner with two compressor units, rather than one, may prevent system downtime. A back-up (redundant) air conditioner permits CT system operation during an extended repair of the primary air conditioner.



## Section 2.0 Cooling Requirement

Use [Table 4-1](#) to assist in cooling requirements planning. Over half the cooling used by your scanner is required for gantry operation. Locate a wall air-conditioning vent at floor level beside and behind gantry to meet both gantry cooling needs and provide patient comfort. Do not locate any cooling vents directly above the gantry. Air returns above the gantry are recommended.

| SYSTEM COMPONENT  | BTU/HR        | WATT          |
|---|---------------|---------------|
| Gantry maximum  | 30,570        | 8,964         |
| Table   | 1030          | 300           |
| Power Distribution Unit   | 3400          | 1000          |
| <b>Scan Room Subtotal</b>   | <b>35,000</b> | <b>10,264</b> |
| <b>Recommended Control Room Subtotal</b><br>(w/o option):   | <b>8,189</b>  | <b>2,400</b>  |
| Operator's console w/ 2 monitor & PMT   |               |               |
| <b>System Total (Recommended)</b><br><b>(See NOTE 1)</b>  | <b>43,189</b> | <b>12,664</b> |
| <b>Option:</b> Remote Color Monitor   | 170           | 50            |
| <b>Option:</b> UPS  | 2900          | 850           |
| <b>ROOM TOTAL (SEE NOTE 2)</b>  | <b>46,259</b> | <b>13,564</b> |
| <b>NOTE 1:</b> With 75 scan rotations per patient:<br>Recommended BTU/hr. provides for up to six patients per hour. It is also needed during calibration of the system. |               |               |
| <b>NOTE 2:</b> Cooling requirements do not include cooling for room lighting, personnel or non-CT equipment.  |               |               |

Table 4-1 Cooling Requirements (Worksheet)

## Section 3.0 Altitude

The system shall meet all functional and performance specification when placed in a room that is at an elevation of -150 m to 2,400m (-492 ft to 7,875 ft) above sea level.

Note: For site with altitudes 2,400 m to 3,048 m (7,875 ft to 10,000 ft), you need a deviation to site a product at this altitude. Altitudes above 2,400 m (7,875 ft) require engineering approval.

## Section 4.0

# Electro-Magnetic Interference (EMI)

Note: If power sub-stations exist under or above the scan room, or near the control room, consider EMI testing to determine if your proposed room meets the published acceptable EMI room limits. This also includes high voltage lines under the scan or control room floor.

### 4.1 Gantry

Locate gantry in ambient static magnetic fields of less than  $10^{-4}$  tesla (1,000 milligauss) to guarantee specified imaging performance. Ambient AC magnetic fields must be below  $10^{-6}$  tesla (10 milligauss) peak.

### 4.2 Console / Computer Equipment

Locate computer equipment in ambient static magnetic fields of less than  $10^{-3}$  tesla (10,000 milligauss) to guarantee data integrity. See [Figure 4-1](#).

### 4.3 Magnetic Media

Locate magnetic media in ambient static magnetic fields of less than  $10^{-3}$  tesla (10,000 milligauss).

### 4.4 PDU

The PDU produces an electromagnetic field that radiates outward from its cabinet in all directions. Do not place the gantry or patient table within 0.3 meters (12 in.) of the edge of the Power Distribution Unit. Do NOT place other sensitive electronics (e.g. the computer equipment) within 1 m (39 in.) of the edge of the Power Distribution Unit in any direction, including above or below it. See [Figure 4-1](#).

Note: The UPS is not classified as sensitive electronics.

### 4.5 EMI Reduction

If fields of excessive EMI are known or suspected to be present, consult GE Sales & Service for recommendations. Consider the following if you attempt to reduce EMI:


- External field strength decreases rapidly with distance from source of magnetic field.
- External leakage magnetic field of a three-phase transformer is much less than that of a bank of three single phase transformers of equivalent power rating.
- Large electric motors are a source of substantial EMI.
- High-powered radio signals are a source of EMI.

Maintain good screening of cables and cabinets.

## 4.6 UPS

The Uninterruptable Power Supply (UPS) provides a consistent power supply to various electrical components of the system. Also, it continues to provide electrical power to components during a site-wide power outage so components can be safely shut down. The UPS should be kept at least one meter (1m) away from sensitive electronics (the PDU does not include sensitive electronics). For UPS interconnect information, please refer to [Chapter 8 - Interconnection Data, Section 6.0 on page 138](#).

## 4.7 Equipment EMI “Envelopes”



**Figure 4-1 Sample Room Layout, showing approximate EMI requirements**



## Section 5.0

# System Component Noise Levels

**Maximum Gantry Audible Noise Level** The maximum ambient noise level is produced by the gantry during a CT scan acquisition. It is less than 75dBA when measured at a distance of one meter from the nearest gantry surface, in any direction.

**Maximum Console Audible Noise Level** The maximum audible noise level is less than 56dBA when measured one meter up and one meter away from the console at an ambient temperature of 26°C.



# Chapter 5

## Floor Loading and Weights

### Section 1.0 Floor Loading

The systems have a total floor load of approximately 3119kg (6869lbs). About 2561kg (5642lbs), including patient, is concentrated in the table-gantry assembly.

Table 5-1 lists CT components with weight, size, floor loading and normal mounting requirements.

| ITEM                                    | NET WEIGHT<br>KG / LB / N              | OVERALL<br>W X D<br>MM (INCH) | LOAD PATTERN<br>MM (IN.)  | NORMAL METHOD OF<br>MOUNTING MM (IN.)<br>(GE SUPPLIED) <sup>1</sup>   |
|---|--|-------------------------------|---|---|
| Gantry (w/ covers)                      | ~1816 / ~4003 / ~17808                 | 2439 X 1006<br>(96.0 X 39.6)  | Rectangular base plate 700 X 2165.7 (27.6 X 85.26) with four round pads, each 63.5 (2.5) in contact with floor.<br><br>Individual pad loadings are: 499 KG (1100 lb), 518 KG (1140 lb), 558 KG (1230 lb), and 593 KG (1306 lb) (see Figure 5-5) | Hilti Kwik-Bolt II 12.7mm (1/2 in) diameter by 203mm (8 in) long per P/N 2106573 at four leveling pads into concrete floor.         |
| Dollies (each)                          | 114 / 250 / 1118                       |                               |   |   |
| Top Cover (each)                        | 11.2 / 24.6 / 110                      |                               |   |   |
| Side Cover (each)                       | 11.3 / 25 / 111                        |                               |   |   |
| Front Cover                             | 43 / 95 / 422                          |                               |   |   |
| Rear Cover                              | 45 / 100 / 441                         |                               |   |   |
| Patient Table: GT 1700 w/ 227kg patient | 446 / 982 / 4370<br>673 / 1481 / 6595  | 650 x 2370<br>(25.6 x 93.3)   | Four round 63.5 mm (2.5 in.) pads: 19.7 x 40.3  | Hilti Kwik-Bolt II 12.7mm (1/2 in) diameter by 203mm (8 in) long per P/N 2106573 at leveling pads into concrete floor.              |
| High capacity table w/295kg patient     | 507 / 1117 / 4970<br>802 / 1767 / 7860 | 650 x 2910<br>(25.6 x 114.5)  | Four round 63.5 mm (2.5 in.) pads: 19.7 x 40.3  | Hilti Kwik-Bolt 5/8 in.diameter by 8 1/2 in.long per P/N 5314211 at the leveling pads into concrete floor.                          |
| Footswitch (GT)                         | 15 / 33 / 147                          | -                             | -   | -   |
| Power Distribution Unit                 | ~370 / ~815 / ~3626                    | 711 X 559<br>(28 X 22)        | Four Casters support area of 711 X 559 (28 X 22) .  | Casters are for positioning and service. Set on floor. May be anchored to floor using angle brackets <sup>2</sup> in seismic zones. |
| True-In-One Console w/o monitors        | 87 / 192 / 853                         | 470 x 740<br>(18.5 x 29.2)    |   |   |
| NIO16 Console w/o monitors              | 80 / 176 / 784                         | 470 x 740<br>(18.5 x 29.2)    |   |   |
| Monitor - LCD (ea)                      | 10 / 22 / 98                           | 415 x 453<br>(16.3 x 17.8)    |   |   |
| Freedom Workspace (5168666-3)           | 44 / 97 / 431                          | 1300 x 620<br>(51 x 24)       |   |   |

Table 5-1 RT System Floor Loading

| ITEM                               | NET WEIGHT<br>KG / LB / N | OVERALL<br>W X D<br>MM (INCH) | LOAD PATTERN<br>MM (IN.) | NORMAL METHOD OF<br>MOUNTING MM (IN.)<br>(GE SUPPLIED) <sup>1</sup> |
|------------------------------------|---------------------------|-------------------------------|--------------------------|---|
| Optima table<br>(5371587)          | 57 / 126 / 559            | 1300 x 895<br>(51 x 35)       |                          |   |
| Aurora<br>table_SWS<br>(5449758-2) | 40 / 88 / 392             | 1300 x 850<br>(51 x 33)       |                          |   |

**Notes:**

- 1.) Use the GE Supplied mounting hardware ONLY IF APPROVED by qualified personnel.  
[See statements in [5.1 - Table and Gantry Mounting Requirements.](#)]
- 2.) Seismic angle brackets are included and shipped with the PDU.
- 3.) Seismic angle brackets are included on the shipping skid.

**Table 5-1 RT System Floor Loading**



## 1.1 Floor Loading and Anchoring Guidelines

Follow the floor loading and anchoring guidelines below when preparing a site for system installation:

- The table and gantry require secure anchoring to the scan room floor. The power distribution unit and the console sit on the floor with casters; anchoring of these components to the floor is optional, unless required because of seismic considerations.
- For total floor load of a Discovery CT590RT with a GT1700 / High capacity table and no UPS refer to [Table 5-1](#).
- Anchors mount through the table and gantry supports. Use the floor template or its dimensions to locate the table and gantry support positions within the scan room, making sure that any anchors that pass through the supports clear all structural beams and interferences in the floor.
- If a loading analysis determines that the gantry and table position should change relative to their position on the GE site print, be sure to take into account the clearance requirements when determining an appropriate location for the system.
- Hospitals and scanning facilities throughout the world may utilize a variety of floor types, and the disposition of different floor types may necessitate additional planning to adequately accommodate the system:
  - Wood floors often require substantial reinforcement. GE does not recommend using wood floors.
  - Temperature variation in blacktop or marble floors may allow anchor movement and pullout. GE does not recommend using these floors.
  - GE recommends using concrete floors with a minimum thickness of 102 mm (4 in.) for Gantry, GT1700 Table and 110 mm (4-1/3 in.) for High capacity Table, when using GE-supplied anchoring or any other equivalent anchoring method.



### NOTICE

Responsibility for providing an approved support structure and mounting method for all floor types other than the GE-recommended floor rests with the purchaser. General Electric accepts no responsibility for any failure of the support structure or anchoring method, including those used for seismic mounting. GE accepts no responsibility for methods other than those listed.

## 1.2 GE Supplied Anchoring

GE supplies anchors for mounting the table and gantry. The console and power distribution unit do not require anchoring to the floor. It is the responsibility of the customer to have a structural engineer and trained contractor use either the GE-supplied anchoring method or to provide an equivalent anchoring method to mount the table and gantry to the floor.

Consult your architect, structural engineer, contractor, or PMI to resolve any questions.



### WARNING

#### POTENTIAL FOR PATIENT INJURY!

AN IMPROPERLY SECURED TABLE MAY TIP, DISLODGING THE PATIENT. PATIENT SAFETY DURING SYSTEM OPERATION REQUIRES PROPER ANCHORING OF SYSTEM COMPONENTS.

### 1.2.1 Specifications of GE-supplied Anchors

[Table 5-2](#) lists the specifications of GE-supplied anchors for the system. There are three types of anchors used in this product depending on manufacturing date and table type. For a detailed view, including dimensions and additional specifications, see [Figure 5-1](#)~ [Figure 5-3](#) of this section.

| PART NUMBER | 2106573           | 5487992-2         | 5314211 (FOR HIGH CAPACITY TABLE |
|-------------|-------------------|-------------------|----------------------------------|
| Description | Hilti Kwik Bolt 2 | Hilti Kwik Bolt 3 | Hilti Kwik Bolt 3                |
| Diameter    | 12.7 mm (0.5 in.) | 12.7 mm (0.5 in.) | 15.9 mm (0.625 in.)              |
| Length      | 203 mm (8 in.)    | 178 mm (7 in.)    | 216 mm (8-1/2 in.)               |

Table 5-2 GE-Supplied Anchor Specifications

## Section 2.0

### Mounting Data, Including Seismic



#### WARNING

**POTENTIAL FOR PATIENT INJURY.**

**IMPROPERLY SECURED TABLE MAY TIP, DISLODGING PATIENT.**

**PROPER ANCHORING IS KEY TO MAINTAINING PATIENT SAFETY DURING SYSTEM OPERATION.**

The following pages show center-of-gravity information for system components:

- Gantry: [Figure 5-4](#)
- Table: [Figure 5-6](#)
- Power Distribution Unit: [Figure 5-8](#)
- Operator's Console/Computer: [Figure 5-10](#)

Floor mounting hole locations for components that don't have templates are also in this section.

Customer is responsible for seismic mounting. Refer to all applicable codes in your area.

GE provided floor anchors ([Figure 5-1](#)) are designed to be used ONLY on concrete floors that meet the concrete floor requirement. Supplied floor anchors must be installed by a trained contractor, and shall be set to a minimum depth of 3 inches at each anchor point. Any anchors having more than 1 inch of thread showing above the nut, when torqued to 55 ft.-lbs, (60 ft.-lbs for High capacity table) shall have a second anchor installed in the closest adjacent hole. the second anchor shall be installed to the standard depth and torque specifications.

| MOUNTING REQUIREMENTS                 | ANCHOR P/N<br>2106573             | ANCHOR P/N<br>5487992-2           |
|---------------------------------------|-----------------------------------|-----------------------------------|
| Minimum Floor Thickness:              | 102 mm (4 in.)                    | 102 mm (4 in.)                    |
| Recommended Drilling Depth:           | 91 mm (3.58 in.)                  | 85 mm (3-3/8 in.)                 |
| Minimum Anchor Embedment:             | 78 mm (3.07 in.)                  | 75 mm (2.95 in.)                  |
| Available Alternate Anchor Locations: | Yes                               | Yes                               |
| Shipped Anchor Size:                  | 203 mm (8 in.)                    | 178 mm (7 in.)                    |
| Alternate Anchoring Methods:          | Yes (see notes, below)            | Yes (see notes, below)            |
| Floor Levelness Requirement:          | 6 mm (0.24 in.) over 3 m (10 ft.) | 6 mm (0.24 in.) over 3 m (10 ft.) |

**Table 5-3 Gantry and GT1700 Table Mounting Requirements**

| MOUNTING REQUIREMENTS                 | ANCHOR P/N (GANTRY)               |                                   | ANCHOR P/N (TABLE)                |
|---------------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                                       | 2106573                           | 5487992-2                         | 5314211                           |
| Minimum Floor Thickness:              | 102 mm (4 in.)                    | 102 mm (4 in.)                    | 110 mm (4.33 in.)                 |
| Recommended Drilling Depth:           | 91 mm (3.58 in.)                  | 85 mm (3.35 in.)                  | 100 mm (3.94 in.)                 |
| Minimum Anchor Embedment:             | 78 mm (3.07 in.)                  | 75 mm (2.95 in.)                  | 90 mm (3.54 in.)                  |
| Available Alternate Anchor Locations: | Yes                               | Yes                               | Yes                               |
| Shipped Anchor Size:                  | 203 mm (8 in.)                    | 178 mm (7 in.)                    | 216 mm (8.5 in.)                  |
| Alternate Anchoring Methods:          | Yes (see notes, above)            | Yes (see notes, above)            | Yes (see notes, below)            |
| Floor Levelness Requirement:          | 6 mm (0.24 in.) over 3 m (10 ft.) | 6 mm (0.24 in.) over 3 m (10 ft.) | 6 mm (0.24 in.) over 3 m (10 ft.) |

**Table 5-4 Gantry and High Capacity Table Mounting Requirements**

If the Installers cannot set at four anchors for the GT1700 or High capacity table, the installer must inform the customer that the minimum anchoring cannot be met, and a structural engineering contractor is required to determine the anchoring method and certify that their anchoring meets the stated GE minimum load requirement and torque specifications.

- Note: All other anchoring methods on floor types other than the concrete minimum must be determined at the customer's expense by a structural engineering contractor, and anchoring method must be certified to meet the stated GE minimum load requirement and torque specification.
- Note: If installing the GE scanner on a floor type other than 102 mm (4 in.) concrete floor, all structural specifications in this document must be reviewed and met. For High Capacity Table, Minimum Floor Thickness is 110 mm (4.33 in.)




Figure 5-1 Gantry and Table Anchoring with 2106573 (8 in.) Anchor Bolt (GT 1700)




Figure 5-2 Gantry and Table Anchoring with 5487992-2 (7 in.) Anchor Bolt (GT 1700)




Figure 5-3 Typical Floor Anchor, Table (High capacity table)




Figure 5-4 Gantry Center of Gravity




Figure 5-5 Gantry Anchor Locations







Figure 5-6 Patient Table (CT1) - GT 1700 (Cradle *Mechanical* Limits shown)



CRADLE : INMAX POSITION  
IMS : INMAX POSITION  
TABLE HEIGHT : 991 mm (39)

NOTES:  
DIMENSIONS ARE IN MILLIMETERS (INCHES)

**Figure 5-7 Patient Table (CT1) - High Capacity Table**




Figure 5-8 Power Distribution Unit (NGPDU)




Figure 5-9 NGPDU Seismic Mounting Bracket (dimensions are in mm)




Figure 5-10 True-In-One Operator's Console




Figure 5-11 NIO16 Console Center-of-Gravity

# **Chapter 6**

## ***Delivery and Storage Requirements***

---

This chapter provides information necessary for planning a safe and successful delivery of the system from GE Healthcare to the receiving area of the installation site, and from the receiving area of that facility to the scan suite.

### **Section 1.0**

#### **Delivery to the Facility**

Your Project Manager of Installation will determine the most appropriate means of transporting the system to your facility. However, the type of receiving area at the facility where the installation will occur determines, to a large extent, the method used to transport the system to that facility. When planning for delivery, facilities fall into two general categories: those with a loading dock, and those without a loading dock.

#### **1.1 Loading Dock Deliveries (Preferred methods)**

Facilities with a loading dock in the receiving areas can generally accommodate delivery of the system by semi-tractor trailer. This is the preferred method for system delivery. Dock-to-dock shipment minimizes the possibility of dropping the gantry or damaging other subsystems during the transition from the trailer to the facility. This method also allows for the most efficient packing and unpacking of the system.

#### **1.2 Ground (Non-Loading Dock) Deliveries**

Facilities without a loading dock require a Lift Gate or Tilt Bed truck. Such deliveries require unloading the system components from the truck bed to ground level and then transported to the facility over a smooth surface such as a concrete sidewalk or driveway or paved area. These paved surfaces must be able to support the weight of the subsystems. It may be necessary to protect these surfaces as well.

##### **1.2.1 Lift-Gate Truck**

If a truck equipped with a lift-gate is used, the delivery truck requires a lift gate rated for a 2722.0 kg (3.0 Tons) capacity. When the gantry or table is lowered to ground level they should be lowered at a steady rate using the slowest speed as possible to minimize G-loads when the lift gate reaches the ground. Keep gantry or table level during movement to avoid flipping. Failure to smoothly transition the table and gantry to ground level may cause serious damage to the table, gantry or their transport dollies.

##### **1.2.2 Tilt Bed Truck**

If a tilt bed delivery truck is used, a GE representative shall supervise the delivery of the CT scanner to ensure the system is safely delivered without damage. To avoid damaging the table and gantry, the representative shall direct the driver to attach strapping to the lowest point (not the wheels) of each dolly. When the table or gantry is moved from the back of the delivery truck to ground level, both shall be lowered at the slowest reasonable steady rate until wheel contact is made at ground level. Movement should be temporarily halted when the dolly wheels come in contact with the

ground. Further movement should resume minimizing any G-loads as the final wheels meet ground level. Failure to smoothly transition the table and gantry to ground level may cause serious damage to the table, gantry or their transport dollies.

### 1.2.3 Forklift Truck

A forklift can be used to unload the gantry, provided that the lifting option is ordered and delivered. The system will arrive with a lifting skid attached to the gantry and table. This option cannot be added later as an on-site addition.

### 1.2.4 Rigging

The CT gantry assemblies shall not be lifted by their dollies. The CT gantry assemblies shall not be transported across any surface by any means other than the dollies provided by GE. The CT gantry assemblies have no lifting points on them and are not designed to be lifted by any special rigging attached to the gantry assemblies themselves.



**POSSIBLE SEVERE PERSONAL INJURY OR DEATH.  
THE DOLLIES ARE NOT DESIGNED TO BE USED AS AN ATTACHMENT POINT FOR  
ANY METHOD OF LIFTING THE SUBSYSTEMS.  
ATTACHING LIFTING STRAPS, CABLES OR MECHANISMS TO THE DOLLY  
HANDLES OR ANY OTHER PART OF THE DOLLY IS STRICTLY PROHIBITED.**



#### NOTICE

**If it is determined that the subsystems must be lifted by crane or other lifting method the PM or person responsible for local siting of the system shall NOT proceed with the installation without consulting directly with GE Engineering.**

Lifting the subsystems by crane or other lifting method should always be avoided. All alternate methods of delivery should be evaluated including the removal of any obstructions, doorways, walls, and windows.

If lifting is still required:

- 1.) The entire gantry assembly and both gantry transport side dollies must be placed on a lifting platform. GE does not provide a lifting platform.
- 2.) The entire patient table must be on its dollies and lifted while sitting on a lifting platform.  
The patient table on its dolly shall be lowered to its transport position so the table base is in contact with the platform.
- 3.) The platform must be designed so no lifting straps or cables come in contact with any part of the gantry or table subsystems or its side dollies.



- 4.) The lifting platform shall bear the entire load. No part of the subsystem shall bear any load during the lift.

Note: If delivery requires vertical or horizontal lifting, the PM needs to add the necessary identifier to the order.

## Section 2.0

### Delivery to the Scan Suite

Once at the installation site, conveyance of the system into the scan suite may involve special considerations, such as vertical lifting, or transportation through stairwells, which involves additional planning by the Project Manager of Installation.

#### 2.1 Packaging Dimension

| SUBSYSTEM | LENGTH               | WIDTH                | HEIGHT                |
|-----------|----------------------|----------------------|-----------------------|
| Gantry    | 2620 mm<br>(103 in.) | 1370 mm<br>(54 in.)  | 2270 mm<br>(89.4 in.) |
| Table     | 3200 mm<br>(126 in.) | 900 mm<br>(35 in.)   | 1360 mm<br>(53.5 in.) |
| PDU       | 900 mm<br>(35 in.)   | 700 mm<br>(27.6 in.) | 1230 mm<br>(48.4 in.) |
| OC        | 880 mm<br>(34.6 in.) | 590 mm<br>(23.2 in.) | 840 mm<br>(33 in.)    |

Table 6-1 Packaging Dimension

#### 2.2 Lifting

Both vertical and horizontal lifting require professional riggers. The PMI should always notify CT engineering before attempting either lifting procedure and should make sure that the order includes the necessary lifting fixtures, as both vertical and horizontal fixtures must appear on the order for them to ship with the system.

If delivery requires vertical lifting, the PMI adds the appropriate identifier to the order. The gantry ships in a vertical lifting crate with lifting instructions for riggers.

If delivery requires horizontal lifting, the PMI adds the corresponding identifier to the order. The gantry ships in a horizontal lifting crate with lifting instructions for riggers.

##### 2.2.1 Stairway Deliveries

Stairways with angles at or less than 45 degrees can accommodate delivery of system components. If the site requires delivery through stairwells, the PMI adds the appropriate identifier to the order to ensure proper packaging of the system, and notifies CT engineering before attempting the procedure. The components ship attached to special lifting skids with lifting instruction for riggers.

## 2.3 Floor Protection

GE recommends floor protection along the delivery path from the dock/receiving area to scan room.

## 2.4 Un-Loading and Un-packing the System

Retain the packaging surrounding the following components:

- Console-Shipped on a shock resistant skid. Do not remove the skid.
- UPS-Shipped on a shock resistant skid. Do not remove the skid.

# Section 3.0 Dollies

## 3.1 Installations within the United States

Typically, domestic shipments (shipments within the United States) involve the use of dollies for moving the gantry, table, and console. After completing installation, return the dollies to GE using the shipping document found in Box #1.

## 3.2 Zero Clearance Dollies

Deliveries involving small elevators with a depth of at least 2692 mm (106 in.) require zero clearance dollies. Zero clearance dollies allow movement of the gantry in tight areas. Avoid using them for normal dock or van deliveries. To order zero clearance dollies, go to <http://www.umi-dollyshop.com>.

## 3.3 Tilting Table Dollies

Deliveries involving small elevators with a depth of at least 2438 mm (96 in.) require tilting table dollies. If storing the system prior to installation, do not order tilt dollies. If you are unable to obtain tilt dollies for delivery, substitute riggers in their place. A limited number of tilt dollies exist for U.S. deliveries. To order tilt dollies, go to <http://www.umi-dollyshop.com>.

## 3.4 Installation Outside of the United States

Customers may purchase dollies (B7850LD) for shipments outside of the United States. After removing the system from the crates, DO NOT return dollies shipped outside of the US to GE Healthcare in Milwaukee, WI, USA. Instead, forward them to the local GE office or warehouse. Zero Clearance and Tilting Table dollies can be purchased through UMI, To buy tilt dollies, go to: <http://www.umi-dollyshop.com>.

## Section 4.0

# Gantry Delivery Considerations

### 4.1 Gantry Shipping State

The gantry is shipped with most covers installed. The assembly is mounted between two dollies. See [Figure 6-1](#). Two side rails are bolted to the dollies to stabilize dollies and protect gantry. Use dolly elevating casters to lift gantry off its base and roll it into position.



Figure 6-1 Gantry with Shipping Dollies and Side Rails

### 4.2 Door Openings

Unobstructed door openings, for moving equipment into building, must measure 1400 mm X 2150 mm (55 in. X 85 in.) minimum. Routing through corridors with a width of 2800 mm (9 ft.) also prove helpful.

There are 4 options of Gantry considerations for delivery RT system, all of them are fit with the minimum requirement. please consider the situation in fact:

Note: Option **1) and 2)**, we define the gantry as the primary point, so the corridor must be fit with gantry; Option **3) and 4)**, we define the corridor as the primary point, use a reference and fixed data of corridor.

- 1.) 42" door with 120" corridor (dollies on, side rails removed)
- 2.) 52" door with 125" corridor (dollies on, side rails on)
- 3.) 55" door with 96" corridor (dollies on, side rails removed)
- 4.) 68" door with 96" corridor (dollies on, side rails on)

If you cannot meet the above specifications on side, you can order the Zero clearance dollies from <http://www.umi-dollyshop.com>.

## 4.3 Elevator Requirements

When moving the gantry from the receiving location to the scanning room, pay special attention to elevator size and capacity. Removing side rails and one dolly after placing the gantry in the elevator reduces the gantry width/length and elevator requirements.

Due to gantry component weight differences all weight listed here are averages. This difference can measure  $\pm 18.14$  kg ( $\pm 40$  lb). Contact the elevator manufacturer if the gantry weight exceeds elevator capacity (see [Table 6-2](#)).

| CONFIGURATION                  | LENGTH               | WIDTH                 | HEIGHT              | WEIGHT               |
|--------------------------------|----------------------|-----------------------|---------------------|----------------------|
| Dollies On, Side Rails On      | 3100 mm<br>(122 in.) | 1290 mm<br>(51 in.)   | 2030 mm<br>(80 in.) | 2016 kg<br>(4445 lb) |
| Dollies On, Side Rails Removed | 3100 mm<br>(122 in.) | 1010 mm<br>(39.7 in.) | 2030 mm<br>(80 in.) | 1993 kg<br>(4395 lb) |
| Dollies Off, Covers Off        | 2160 mm<br>(85 in.)  | 860 mm<br>(34 in.)    | 1870 mm<br>(73 in.) | 1655 kg<br>(3648 lb) |

**Table 6-2 Size of Gantry & Dollies, with and without Side Rails**

Minimum hallway and door size for gantry with covers and dollies attached, but side rails removed, is 1016 mm (40 in.). For alternative lifting arrangements and instructions, contact GE Installation Support Services.

## Section 5.0 Table Delivery Considerations

The table is shipped without side covers installed. Covers are shipped in four separate boxes. The table is mounted between two dollies.

Table shipping dimensions are 3200 mm (126 in.) long, 650 mm (25.6 in.) wide, and 1200 mm (47.2 in.) high.

|                                 | LENGTH |     | WIDTH |    | HEIGHT |    | WEIGHT |      |
|---------------------------------|--------|-----|-------|----|--------|----|--------|------|
|                                 | MM     | IN  | MM    | IN | MM     | IN | KG     | LB   |
| High Capacity Table w/Dollies   | 2997   | 118 | 762   | 30 | 1143   | 45 | 639    | 1407 |
| GT1700 w/Dollies                | 2489   | 98  | 762   | 30 | 1143   | 45 | 607    | 1337 |
| High Capacity Table w/o Dollies | 2997   | 118 | 660   | 26 | 889    | 35 | 507    | 1117 |
| GT1700 w/o Dollies              | 2489   | 98  | 660   | 26 | 889    | 35 | 475    | 1046 |

**Table 6-3 GT Table Dimensions - with and without dollies**

|                             | LENGTH |    | WIDTH |    | HEIGHT |    | WEIGHT |     |
|-----------------------------|--------|----|-------|----|--------|----|--------|-----|
|                             | MM     | IN | MM    | IN | MM     | IN | KG     | LB  |
| High Capacity Table Dollies | 2210   | 87 | 762   | 30 | 127    | 5  | 132    | 291 |
| GT1700 Dollies              | 2210   | 87 | 762   | 30 | 127    | 5  | 132    | 291 |
| High Capacity Table Tilting | 762    | 30 | 762   | 30 | 889    | 35 | 136    | 300 |
| GT1700 Tilting              | 762    | 30 | 762   | 30 | 889    | 35 | 132    | 291 |

Table 6-4 GT Table Dolly Dimensions

|  | LENGTH      |          | WIDTH |    | HEIGHT      |       | WEIGHT |      |
|--|-------------|----------|-------|----|-------------|-------|--------|------|
|  | MM          | IN       | MM    | IN | MM          | IN    | KG     | LB   |
| GT1700   | 2489        | 98       | 762   | 30 | 1143        | 45    | 602    | 1325 |
| High capacity table                                | 2997        | 118      | 762   | 30 | 1143        | 45    | 632    | 1392 |
| High capacity table - Tilting (approx. dimensions) | 2489 - 2921 | 98 - 115 | 660   | 26 | 1778 - 2032 | 70-80 | 636    | 1400 |

Table 6-5 GT Table Elevator Delivery Dimensions

## Section 6.0 Console Delivery Considerations

**For TIO Console:** True-In-One console has no table top. FWS table is shipped with the True-In-One console.

The dimensions of the console alone (as shipped) are 740mm (29 in.) deep, 470 mm (19 in.) wide, and 640mm (25 in.) high.

**For NIO16 Console:** The console is shipped without covers installed. The covers are delivered in the Product Grade Collector.

The dimensions of the console alone (as shipped) measure 830 mm (33 in.) deep, 470 mm (19 in.) wide, and 656 mm (26 in.) high.

## Section 7.0

# Storage Requirements



**NOTICE** Failure to adhere to Storage Requirements will likely result in equipment damage.

### 7.1 Short-term Storage (Less than Six Months)

If storing the CT system before installation for less than six months, store in a temperature and humidity controlled warehouse. Protect from weather, dirt and dust. Meeting the following requirements prevents rust and corrosion from forming on bearing surfaces due to condensation:

- Storage temperature should not exceed 0° to 30° C (40° to 80° F).
- Storage relative humidity up to 70% , non-condensing.
- Maximum relative humidity rate of change is 5%/hr.
- The maximum temperature rate of change is 3° C/hr. (5° F/hr.)
- Storage longer than 6 months is not recommended



**NOTICE** Between delivery is considered short-term storage. Van storage must meet the same specifications as above.

### 7.2 Construction-Site Storage

When storing the CT system at a construction site be sure to adhere to the following storage requirements:

- Do not damage or puncture the shipping crate.
- Do not remove packaging until all construction is completed at the site and all dust created by the construction is removed.
- Maintain a storage temperature within the range of 10° to 32° C (50° to 90° F).
- Maintain a relative humidity (non-condensing) between 20% and 70%.

## Section 8.0

# Extreme Temperature Transportation and Deliveries



**NOTICE** Failure to adhere to Extreme Temperature Transportation and Delivery requirements will likely result in equipment damage.

Extreme temperatures should be avoided during system transportation and delivery.

Extreme temperatures consist of temperatures below -18° C (0° F), or above 49° C (120° F), without humidity control.

When transporting the CT system, prevent extended exposure of the system to temperatures or humidity outside of the following specifications:

- Time: Up to two weeks duration
- Temperature: -40° to +70° C (-40° to +158° F)
- Humidity: 10% to 100%, including condensing

- Altitude: -1,800ft to 18,000ft



**NOTICE**

**Component Freezing occurs if CT system is exposed to temperatures below -18° C (0° F) for a period longer than two days.**

**Allow a minimum of 12 hours for the CT system to adjust to ambient room temperature, prior to installation.**

## Section 9.0 Site Environmental Considerations

### 9.1 Dust/Dirt Contamination

The system (consisting of: Console, PDU, Table and Gantry) are highly susceptible to airborne contaminants, especially concrete and drywall dust. Due to the possibility of contamination, these systems should NEVER be installed in a construction site. Any site with unfinished floors, walls or ceilings is considered a construction site, and is not suitable for system installation.



**NOTICE**

**The act of installing a GE CT scanner in a construction (i.e., unfinished) site will likely result in the following adverse effects:**

- **Increased installation time**
- **Decreased installation quality**
- **Increased scanner downtime, due to increased service calls**

### 9.2 Chemical Contamination

Wet film processors must never be installed in the same room as the scanner, due to the possibility of chemical contamination of system CT scanner components. Such chemicals can contribute to increased equipment failures, increased system downtime, and decreased reliability. Film processor equipment installation must meet the manufacturer's requirements (e.g. ventilation specifications) and all applicable national and local codes. Also, consideration's should be given to the location of this equipment and chemical fumes relative to human contact as it relates to locating this equipment and chemicals in the control room.

## Section 10.0

# Handling Requirements

Communicate the information in this chapter to any personnel who will transport, move, or otherwise handle the system components during transportation and delivery of the system.

### 10.1 Transportation

To avoid dropping the gantry, it is recommended that the system is transported from GE Healthcare to the facility of the installation site, shipping dock-to-dock in a van. However, facilities without a loading dock may transport the system using lift-gate or flatbed trucks, provided that no dropping or mis-handling of the system occurs. These methods involve unloading system components from the truck and then rolling them across SMOOTH sidewalks or other paved surfaces.

### 10.2 Handling Requirements

The design of the system does not tolerate dropping, shock, vibration, tipping, or hoisting. Be sure to communicate these handling requirements to all parties involved in transporting, moving, and handling system components.

#### 10.2.1 Avoid Dropping

Never drop the gantry, console, table, or PDU. A drop from a height greater than 13 mm (0.5 in.) may cause structural damage to the frame or other major components. Damage resulting from a drop (e.g., bent frame, misalignment) may not become apparent until after the system is installed.

#### 10.2.2 Avoid Shocks and Vibrations

The design of the system, including the gantry, console, table, and PDU, does not tolerate excessive shock or vibration, which may occur during unloading. For example, rolling the console across a "washboard" style ramp may vibrate components, causing loose or broken connections. Damage resulting from shock or vibration (e.g., monitor, CD-ROM, hard-drive, or console failure) may not become evident until after the system is installed.

#### 10.2.3 Avoid Tipping

All system components must remain upright at all times; avoid tipping them. Move the gantry by rolling it on its dollies ONLY, do NOT hoist it. Avoid tipping or lifting the gantry when moving it through hallways, doorways, elevators, etc.



#### NOTICE

**Never lift the gantry with a forklift. Lifting the gantry requires engineering approval for each occurrence. Your GE PMI should contact CT Engineering for all special lifting requirements, as unauthorized gantry lifting can cause gantry bearing damage.**



## 10.3 Inclines and Flat-bed Truck Removal

Inclines and Flat-bed Truck Removal wrecker, attach the straps to the LOWEST possible point on the dolly, and lower the gantry at the SLOWEST reasonable rate, (see [Figure 6-2](#)).




Figure 6-2 Proper Gantry Strap Location



**WARNING** SOME ASSEMBLIES ARE TOP-HEAVY. BE CAREFUL NOT TO TIP.



# ***Chapter 7***

## ***Power Requirements***

---

### **Section 1.0**

#### **Introduction**

The power distribution unit (PDU) supplied with the system transforms and distributes power to all system components. The PDU is the only power entry point required to operate system.

To minimize voltage regulation effects, power wiring between the facility main distribution panel and the PDU should be kept as short as possible.

When routing the power wiring all three phase wires and ground must be run in the same conduit or raceway duct. Power wires should be routed separately from system control and signal cables, using a separate conduit or trough in raceway duct.

Metallic conduit, floor duct or surface raceway may be used for running cables, depending upon local codes and practices. However, cable passageways should be large enough to install any cable with all other cables already installed. Use of non-metallic conduit is not recommended.

## Section 2.0

# System Input Power

### 2.1 Power Source Configuration

The CT Scanner is designed to operate on a three-phase, four-wire wye power source. A solidly grounded wye source is preferred. The neutral wire does not need to be run to the system, i.e., four-wire connection. If a neutral wire is run, then it should be terminated in the A1 box.

A dedicated feeder from the nearest Main Distribution Panel (MDP) should supply power to the scanner. In accordance with the National Electric Code (U.S.) and similar applicable national and local codes, a protective disconnect device must be provided in the power line supplying the PDU. It must be located within 10 m (32 ft.) of the PDU, visible to PDU service personnel, and must have "lockout / tagout" provisions. This disconnect device is identified as "A1" in the interconnection schematic diagrams.

### 2.2 Rating

The system operates on three-phase power meeting the following specifications.

|           |                      |
|-----------|----------------------|
| Voltage   | 380 to 480 VAC       |
| Capacity  | 150 kVA              |
| Frequency | 50 or 60 Hz +/- 3 Hz |

- Maximum power demand = 150 kVA @ 0.85 PF at a selected technique of 140 kV, 715 mA.
- Average (continuous) power demand at maximum duty cycle = 30 kVA.

The "A1" disconnect device referenced above must provide overcurrent protection for the system and facilitate multi-point remote "Emergency Off" control of system power. A disconnect utilizing undervoltage release control is preferred over shunt trip devices. The rating of the "A1" disconnect device depends on the nominal line voltage at the site. Refer to [Section 3.0: Recommended Power Distribution System](#), for minimum rating requirements and suggested disconnect devices.



**TO PREVENT POWER LOSS TO OTHER LOADS IN CASE OF AN UNEXPECTED CT OR PET SYSTEM FAULT, THE POWER FEEDER MUST HAVE OVERCURRENT PROTECTION SUCH THAT THE DOWN-STREAM OVERCURRENT PROTECTION DEVICES (E.G. GE A1 PANEL) CLEAR THE FAULT BEFORE ANY UP-STREAM OVERCURRENT PROTECTION DEVICE OPENS.**

## 2.3 Regulation

Total load regulation as measured at the PDU input terminals must not exceed 6%. The capacity of the facility transformer and size & length of feeder wires directly affect the load regulation presented to the system. Refer to [Section 3.0: Recommended Power Distribution System](#), for recommended single-unit installation specifics.

## 2.4 Phase Imbalance

The difference between the highest line-to-line voltage and lowest line-to-line voltage must not exceed 2% of the lowest line-to-line voltage.

## 2.5 Sags, Surges & Transients

Sags and surges of the power line must not exceed the absolute range limits shown in [Table 7-1](#). Maximum transient voltages should be limited to 1500V peak.

## 2.6 Grounding

The customer's electrician needs to perform the following tasks:

- Bond metal conduit, raceway, or the armor of armored cable used to power the system to the PDU cabinet and to the A1 Disconnect
- Run a dedicated 1/0 (55 mm<sup>2</sup>) or larger insulated copper ground wire from the main distribution panel to the PDU with the phase wires.
- Run the ground wire with the three-phase wires from the power source to the A1 Disconnect and from A1 Disconnect to the PDU. Grounding does not require a neutral wire.

Note: The shield or armor of armored cable ALONE does NOT provide sufficient grounding.

Bond the ground wire to the intermediate distribution panels through which it passes in accordance with local codes. The resistance between the PDU ground and the facility earth ground must not exceed 0.5 ohm. In addition, the total resistance between the PDU ground and earth must not exceed 2 ohms.

## Section 3.0

# Recommended Power Distribution System

A dedicated feeder run from the facility main isolation transformer is recommended to power the CT scanner. If the scanner must be powered from an existing distribution transformer and secondary feeder, such as the equipment distribution panel of an X-ray department, installation with other X-Ray equipment that use rapid film changers should be avoided. These changers use a large number of high powered, closely spaced exposures, which may coincide with the CT scan and produce image artifacts.



**IF THE POWER FEED FOR THE A1/PDB PANEL IS NOT ON A DEDICATED POWER TRANSFORMER ANY DEVICE THAT SHARES POWER FROM THAT TRANSFORMER MAY BE IMPACTED BY INADVERTENT POWER INTERRUPTION CAUSED BY AN A1/PDB POWER PANEL FAULT. CONVERSLY, THE OPERATION OF OTHER DEVICES SHARING THE POWER TRANSFORMER MAY ALSO IMPACT THE OPERATION OF THE CT/PET SCANNER.**

If a dedicated distribution transformer is provided for the scanner, the minimum recommended transformer size is 225 kVA, rated 2.4% regulation at unity power factor. For this configuration, the minimum recommended feeder size and overcurrent protection device based on line voltage is shown in [Table 7-2](#) Minimum Feeder Wire Size.

In all cases, qualified personnel must verify that the transformer and feeder, at point of take-off, plus the run to the CT scanner meet all the requirements stated in this document.

### SYSTEM CHARACTERISTICS

- Maximum power demand = 150kVA @ 0.85 PF: at a Selected Technique of 140 kV, 715 mA.
- Continuous (average) power demand at maximum duty cycle = 30kVA.
- Maximum allowable total source regulation is 6%.
- Minimum recommended transformer size: 225 kVA, with 2.4% rated regulation at unity power factor. Resultant maximum allowable feeder regulation is 3.4%

### NOMINAL LINE VOLTAGE

The nominal line voltage must fall within one of the ranges listed below

| Nominal Line Voltage                          | 380 | 400 | 420 | 440 | 460 | 480 |
|---|-----|-----|-----|-----|-----|-----|
| Hi-Line Limit, +10%                           | 418 | 440 | 462 | 484 | 506 | 528 |
| Lo-Line Limit, -10%                           | 342 | 360 | 378 | 396 | 414 | 432 |
| Continuous Line Current                       | 38  | 36  | 34  | 33  | 31  | 30  |
| Momentary Line Current                        | 228 | 217 | 206 | 197 | 188 | 180 |
| Maximum Line Current                          | 253 | 241 | 229 | 219 | 209 | 200 |
| Minimum Recommended Circuit Protection Rating | 150 | 150 | 150 | 125 | 125 | 125 |

Table 7-1 Nominal Line Voltage

## MINIMUM FEEDER WIRE SIZE

| FEEDER LENGTH<br>(MDA TO A1)<br>FEET (METERS) | MINIMUM FEEDER WIRE SIZE, AWG OR MCM<br>(SQ. MM)/ VAC |           |          |          |          |          |
|---|---|-----------|----------|----------|----------|----------|
|   | 380 VAC   | 400 VAC   | 420 VAC  | 440VAC   | 460VAC   | 480VAC   |
| 15 m (50 ft)                                  | 1/0 (55)  | 1/0 (55)  | 1/0 (55) | 1 (45)   | 1 (45)   | 1 (45)   |
| 30 m (100 ft)                                 | 1/0 (55)  | 1/0 (55)  | 1/0 (55) | 1 (45)   | 1 (45)   | 1 (45)   |
| 46 m (150 ft)                                 | 1/0 (55)  | 1/0 (55)  | 1/0 (55) | 1 (45)   | 1 (45)   | 1 (45)   |
| 61 m (200 ft)                                 | 1/0 (55)  | 1/0 (55)  | 1/0 (55) | 1 (45)   | 1 (45)   | 1 (45)   |
| 76 m (250 ft)                                 | 2/0 (70)  | 2/0 (70)  | 1/0 (55) | 1/0 (55) | 1 (45)   | 1 (45)   |
| 91 m (300 ft)                                 | 3/0 (85)  | 3/0 (85)  | 2/0 (70) | 2/0 (70) | 1/0 (55) | 1/0 (55) |
| 107 m (350 ft)                                | 4/0 (100)   | 3/0 (85)  | 3/0 (85) | 2/0 (70) | 2/0 (70) | 1/0 (55) |
| 122 m (400 ft)                                | 250 (125)   | 4/0 (100) | 3/0 (85) | 3/0 (85) | 3/0 (85) | 2/0 (70) |

**Table 7-2 Minimum Feeder Wire Size**

Note: In all cases the recommended ground wire is a 55 sq. mm (1/0) ground wire.

## MINIMUM SUB-FEEDER WIRE SIZE

| SUB-FEEDER LENGTH<br>(A1 TO PM) FEET<br>(METERS) | MINIMUM SUB-FEEDER WIRE, AWG OR MCM<br>(SQ. MM) |          |          |        |        |        |
|--|---|----------|----------|--------|--------|--------|
|  | 380 VAC   | 400 VAC  | 420 VAC  | 440VAC | 460VAC | 480VAC |
| 9.7536 m (32 ft)                                 | 1/0 (55)  | 1/0 (55) | 1/0 (55) | 1 (45) | 1 (45) | 1 (45) |

**Table 7-3 Minimum Sub-Feeder Wire Size**

- 1.) [Table 7-1](#), [Table 7-2](#), and [Table 7-3](#) above are based on the use of copper wire, rated 75 C and run in steel conduit. Ampacity is determined in accordance with the National Electrical Code (NFPA 70), Table 310-16 (2002)
- 2.) The minimum feeder size is determined by the ampacity of the circuit protection device listed above, except where a larger size is necessary to meet total source regulation limits.

## Section 4.0 Ground System

The CT Scanner has been designed to use an equal potential grounding system. The required ground system is shown in Figure 7-1. There are three primary grounding points:

- A system power ground point located in the PDU.
- A reference ground point located between gantry and table base.
- A patient ground point located at the front of the table base.

All exposed metal surfaces in the patient vicinity are grounded to the reference ground point.

For additional information, refer to Electrical Safety Equipment, Direction 46-014505.




Figure 7-1 System Ground Map



# Chapter 8

## Interconnection Data

### Section 1.0 Introduction

Figure 8-4 shows interconnection runs for a 50/60 Hz system.

Table 8-1 shows component designators for supplied equipment and options and wall power outlets.

Table 8-11 lists customer-installed wiring and supplied cables. Actual length of each run is less than the length of supplied cables to allow for routing inside equipment. Cable diameters and sizes of connectors are provided to aid in sizing conduit and access plates.

Table 8-2 and Table 8-5 list details for connection to the equipment, using standard (short) length and non-standard (long) length cables, respectively. Details are listed for the following types of runs as appropriate:

- Flush-floor duct
- Computer floor
- Through-wall bushing
- Junction box
- Surface floor duct
- Through-floor duct
- Wall duct
- Conduit

Need for additional junction boxes is minimized by use of either a cable raceway system or a raised computer floor. Systems use prefabricated cables with large plugs. Therefore, conduit or pipe is not recommended for cable runs.

Note: Use dry cleaning for electro components.

### Section 2.0 Component Designators

| DESIGNATOR | APPLIES TO                    | SOURCE              |
|------------|-------------------------------|---------------------|
| A1         | Primary power disconnect      | Contractor supplied |
| CT1        | Patient table                 | System              |
| CT2        | Gantry                        | System              |
| ITL        | InSite telephone lines        | Contractor supplied |
| LP         | Line printer                  | Option              |
| OC1        | Operator's console/computer   | System              |
| PDU        | Power distribution unit       | System              |
| SEO        | System emergency off          | Contractor supplied |
| SM         | Slave monitor                 | Option              |
| WL         | "X-ray on" warning light      | Contractor supplied |
| DS         | Door Interlock Switch         | Contractor supplied |
| XCVR       | Ethernet transceiver          | System              |
| BBNC       | Broad-Band Network Connection | Contractor supplied |

Table 8-1 Component Designators

## Section 3.0 Interconnect Runs, Wiring and Cables

### 3.1 GE Supplied (Standard Length)

| RUN # | LENGTH, ACTUAL (USABLE) |             | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |
|-------|-------------------------|-------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|
|       | ft                      | m           |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |
| 050   | 28 (20)                 | 8.5 (6.1)   | 2343529-2 | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |
| 051   | 28 (20)                 | 8.5 (6.1)   | 2343530-2 | HVAC, Axial Drive Power PDU to Gantry    | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |
| 052   | 28 (20)                 | 8.5 (6.1)   | 2343528-2 | LVAC, PDU to Gantry 120VAC Power         | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |
| 053   | 65 (60)                 | 19.8 (18.3) | 2343531-2 | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |
| 054   |                         |             | 2343596   | LVAC, Gantry to Table                    | 1015                 |              | 600            | 120VAC         |                  |                | 3         | 14          |  |
| 055   | 28 (20)                 | 8.5 (6.0)   | 2371450-2 | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |
| 056   | 68 (57)                 | 20.8 (17.4) | 2371450-4 | Ground, Raceway to Console               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |
| 100   | 32.5 (20)               | 9.9 (6.1)   | 5120646-2 | Signal, PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 101   | 71 (60)                 | 21.7 (18.3) | 5419981-2 | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 102   | 71 (63)                 | 21.7 (19.3) | 2373436-3 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |
| 103   | 68 (60)                 | 20.7 (18.3) | 2117848-7 | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |
| 104   |                         |             | 2333151   | Signal, Gantry to Table                  |                      | FT-4         | 300            |                | 80               |                | 25        | 22          |  |

Table 8-2 GE Supplied Cables for GOC or TIO (Standard Run) (2281840-5)- UL Information

| RUN # | LENGTH, ACTUAL (USABLE) |             | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |
|-------|-------------------------|-------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|
|       | ft                      | m           |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |
| 050   | 28 (20)                 | 8.5 (6.1)   | 2343529-2 | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |
| 051   | 28 (20)                 | 8.5 (6.1)   | 2343530-2 | HVAC, PDU to Gantry                      | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |
| 052   | 28 (20)                 | 8.5 (6.1)   | 2343528-2 | LVAC, PDU to Gantry                      | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |
| 053   | 65 (60)                 | 19.8 (18.3) | 2343531-2 | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |
| 054   |                         |             | 2343596   | LVAC, Gantry to Table                    | 1015                 |              | 600            | 120VAC         |                  |                | 3         | 14          |  |
| 055   | 28 (20)                 | 8.5 (6.0)   | 2371450-2 | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |
| 056   | 68 (57)                 | 20.8 (17.4) | 2371450-4 | Ground, Raceway to Console               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |
| 100   | 32.5 (20)               | 9.9 (6.1)   | 5419992-2 | Signal, PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 101   | 71 (60)                 | 21.7 (18.3) | 5419981-2 | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 102   | 71 (63)                 | 21.7 (19.3) | 2373436-3 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |
| 103   | 68 (60)                 | 20.7 (18.3) | 5432019   | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |
| 104   |                         |             | 2333151   | Signal, Gantry to Table                  |                      | FT-4         | 300            |                | 80               |                | 25        | 22          |  |

Table 8-3 GE Supplied Cables for NIO16 Console (Standard Run) (2281840-14)- UL Information

Note: For short cable kits GEHW Cat# is B71022RT, WSO Cat# is B7580JY

| RUN # | LENGTH, ACTUAL (USABLE) |             | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |
|-------|-------------------------|-------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|
|       | ft                      | m           |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |
| 050   | 28 (20)                 | 8.5 (6.1)   | 2343529-2 | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |
| 051   | 28 (20)                 | 8.5 (6.1)   | 2343530-2 | HVAC, PDU Axial Drive Power to Gantry    | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |
| 052   | 28 (20)                 | 8.5 (6.1)   | 2343528-4 | LVAC, PDU to Gantry 120VAC Power         | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |
| 053   | 65 (60)                 | 19.8 (18.3) | 2343531-2 | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |
| 055   | 28 (20)                 | 8.5 (6.0)   | 2371450-2 | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |
| 056   | 68 (57)                 | 20.8 (17.4) | 2371450-4 | Ground, Console to Raceway               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |
| 100   | 32.5 (20)               | 9.9 (6.1)   | 5120646-2 | Signal, PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 101   | 71 (60)                 | 21.7 (18.3) | 5419981-2 | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 102   | 71 (63)                 | 21.7 (19.3) | 2373436-3 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |
| 103   | 68 (60)                 | 20.7 (18.3) | 5432019   | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |

Table 8-4 GE Supplied Cables for NIO16 Console with Simplified Powerpan (Standard Run) (5443710)- UL Information

Note: For short cable kits GEHW Cat# is B75372CB, WSO Cat# is B7580GA






Figure 8-1 System Interconnect Diagram with TIO Console





### 3.2 GE Supplied (Optional, Long Run)

| RUN # | LENGTH, ACTUAL (USABLE) |               | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |  |
|-------|-------------------------|---------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|--|
|       | ft                      | m             |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |  |
| 050   | 63 (55)                 | 19.3 (16.76)  | 2343529   | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |  |
| 051   | 62.5 (55)               | 19 (16.76)    | 2343530   | HVAC, PDU Axial Drive Power to Gantry    | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |  |
| 052   | 63 (58)                 | 19.3 (17.56)  | 2343528   | LVAC, PDU to Gantry 120VAC Power         | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |  |
| 053   | 80 (75)                 | 24.5 (22.86)  | 2343531   | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |  |
| 055   | 63 (55)                 | 19.3 (16.76)  | 2371450   | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |  |
| 056   | 83 (75)                 | 25.5 (22.86)  | 2371450-3 | Ground, Raceway to Console               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |  |
| 100   | 70 (62)                 | 21.4 (18.86)  | 5120646   | Cable - PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |  |
| 101   | 86 (78)                 | 26.35 (23.71) | 5419981   | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |  |
| 102   | 86 (81)                 | 26.3 (24.84)  | 2373436-2 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |  |
| 103   | 80 (75)                 | 24.3 (22.86)  | 2117848-2 | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |  |

Table 8-5 GE Supplied Cables for GOC or TIO (Optional, Long Run) (2281840-4)- UL Information



| RUN # | LENGTH, ACTUAL (USABLE) |               | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |
|-------|-------------------------|---------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|
|       | ft                      | m             |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |
| 050   | 63 (55)                 | 19.3 (16.76)  | 2343529   | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |
| 051   | 63 (55)                 | 19 (16.76)    | 2343530   | HVAC, PDU Axial Drive Power to Gantry    | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |
| 052   | 63 (58)                 | 19.3 (17.56)  | 2343528   | LVAC, PDU to Gantry 120VAC Power         | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |
| 053   | 80 (75)                 | 24.5 (22.86)  | 2343531   | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |
| 055   | 63 (55)                 | 19.3 (16.76)  | 2371450   | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |
| 056   | 83 (75)                 | 25.5 (22.86)  | 2371450-3 | Ground, Raceway to Console               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |
| 100   | 70 (62)                 | 21.4 (18.86)  | 5419992   | Signal, PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 101   | 86 (78)                 | 26.35 (23.71) | 5419981   | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |
| 102   | 86 (81)                 | 26.3 (24.84)  | 2373436-2 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |
| 103   | 80 (75)                 | 24.3 (22.86)  | 5432019   | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |

Table 8-6 GE Supplied Cables for NIO16 Console (Optional, Long Run) (2281840-13)- UL Information

Note: For long cable kits GEHW Cat# is B70992RT, WSO Cat# is B7580JZ

| RUN # | LENGTH, ACTUAL (USABLE) |             | PART #    | DESCRIPTION                              | UL CABLE INFORMATION |              |                |                |                  |                |           |             | PULL SIZE MM (INCHES)                        |  |
|-------|-------------------------|-------------|-----------|--|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|-------------|--|--|
|       | ft                      | m           |           |  | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG    |  |  |
| 050   | 63 (55)                 | 19.5 (16.8) | 2343529   | HVDC, PDU to Gantry                      | 2587                 | FT4          | 600            | + & - 350VDC   | 90               | 19 (.751)      | 3         | (2) 4 (1) 8 | 22 (.87) Dia                                 |  |
| 051   | 63 (55)                 | 19.5 (16.8) | 2343530   | HVAC, PDU Axial Drive Power to Gantry    | 2587                 | FT4          | 600            | 440Y/254       | 90               | 15.3 (.604)    | 4         | 14          | 11.2 (.44) Dia                               |  |
| 052   | 63 (58)                 | 19.5 (16.8) | 2343528-3 | LVAC, PDU to Gantry 120VAC Power         | 2587                 | FT4          | 600            | 208Y/120       | 90               | 13.8 (.542)    | 5         | 8           | 56.4 (2.22) Dia                              |  |
| 053   | 80 (75)                 | 24.5 (22.9) | 2343531   | LVAC, PDU to Console                     | 2587                 | FT4          | 600            | 120VAC         | 90               | 12.3 (.483)    | 3         | 10          | 56.4 (2.22) Dia                              |  |
| 055   | 63 (55)                 | 19.5 (16.8) | 2371450   | Ground, PDU to Raceway                   | 1284                 | VW-1 (FT-1)  | 600            | 0              | 105              | 15.5 (.608)    | 1         | 1/0         | 15.8 (.62) Dia                               |  |
| 056   | 83 (75)                 | 25.5 (22.9) | 2371450-3 | Ground, Raceway to Console               | 1283                 | VW-1 (FT-1)  | 600            | 0              | 105              | 11.9 (.467)    | 1         | 2           | 12.2 (.48) Dia                               |  |
| 100   | 70 (62)                 | 21.4 (18.9) | 5120646   | Signal, PDU Interface to Gantry TGPU J11 |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |  |
| 101   | 86 (78)                 | 26.4 (23.7) | 5419981   | Signal, Console to TGPU J9               |                      | FT-4         | 300            | <30VDC         | 80               | 11.2 (.440)    | 25        | 22          | 17 x 58 (.68 x 2.30)<br>19 x 51 (.75 x 2.01) |  |
| 102   | 86 (81)                 | 26.3 (24.9) | 2373436-2 | Signal (LAN), Gantry to OC               |                      |              | 1900           | <30VDC         |                  | 5.9 (.234)     | 8         | 24          | 15 (.59) Dia                                 |  |
| 103   | 80 (75)                 | 24.5 (22.9) | 5432019   | Fiber Optic, Gantry to OC                |                      |              | N/A            | N/A            |                  |                | 1         | N/A         | 10 (.39) Dia                                 |  |

Table 8-7 GE Supplied Cables for NIO16 Console with Simplified Powerpan (Optional, Long Run) (5438124)- UL Information

Note: For long cable kits GEHW Cat# is B75262CB, WSO Cat# is B7580GB

### 3.3 GE Supplied (Console Cables)

| PART #    | DESCRIPTION                     | CONNECT TO        | QUANTITY | LENGTH     |               |
|-----------|---------------------------------|-------------------|----------|------------|---------------|
|           |                                 |                   |          | MM         | INCHES        |
| 5160575-2 | PS_2 EXTENSION CABLE            | Mouse             | 1        | 3000 ± 100 | 118.11 ± 3.94 |
| 5160574   | USB EXTENSION CABLE             | Keyboard          | 1        | 3000 ± 100 | 118.11 ± 3.94 |
| 5332107-2 | CABLE, DVI to D-SUB VIDEO CABLE | Monitor           | 2        | 3000 ± 20  | 118.11 ± 0.79 |
| 5332093-2 | CABLE, SCIM HD50-HD50           | SCIM              | 1        | 3250 ± 25  | 127.95 ± 0.98 |
| 5315370   | CABLE, USB TYPE A-B             | PMT Tower         | 3        | 2000       | 78.74         |
| 5332100-2 | Power Cord 3.5 meter            | PMT Tower/Monitor | 4        | 3500 ± 30  | 137.80 ± 1.18 |

Table 8-8 GE Supplied TIO Console Cables

| PART #    | DESCRIPTION  | CONNECT TO                               | QUANTITY | LENGTH    |               |
|-----------|--|--|----------|-----------|---------------|
|           |  |  |          | MM        | INCHES        |
| 5366514-2 | USB EXTENSION CABLE  | Keyboard                                 | 1        | 3560 ± 30 | 140.16 ± 1.18 |
| 5450275   | PS_2 EXTENSION CABLE for <b>Z800</b>                                   | Mouse                                    | 1        | 3500 ± 50 | 137.80 ± 1.97 |
| 5366514   | USB EXTENSION CABLE for <b>Z840</b>                                    | Mouse                                    | 1        | 3000 ± 30 | 118.11 ± 1.18 |
| 5332107-2 | CABLE, DVI to D-SUB VIDEO CABLE  | Monitor                                  | 1        | 3000 ± 20 | 118.11 ± 0.79 |
| 5315370   | CABLE, USB TYPE A-B  | PMT media Tower, DVD-RW/USB external HDD | 2        | 2000      | 78.74         |
| 5408703   | DP to DVI cable, 3 meter   | Monitor                                  | 1        | 3000 ± 50 | 118.11 ± 1.97 |
| 5432953-2 | Power Cable, Peripheral Tower to NIO AC Box, IEC C14 to C13 Connectors | PMT media Tower                          | 1        | 3050 ± 50 | 120 ± 1.97    |
| 5432953-3 | Power Cable, Display monitor to NIO AC Box, IEC C14 to C13 Connectors  | Monitor                                  | 1        | 3050 ± 50 | 120 ± 1.97    |
| 5432953-4 | Power Cable, Scan monitor to NIO AC Box, IEC C14 to C13 Connectors     | Monitor                                  | 1        | 3050 ± 50 | 120 ± 1.97    |

Table 8-9 GE Supplied NIO16 Console Cables

### 3.4 GE Supplied (Cables of Options)

| OPTION   | LENGTH, ACTUAL (USABLE) |      | PART #    | DESCRIPTION   | UL CABLE INFORMATION |              |                |                |                  |                |           |          | PULL SIZE MM (INCHES) |
|----------|-------------------------|------|-----------|---|----------------------|--------------|----------------|----------------|------------------|----------------|-----------|----------|-----------------------|
|          | ft                      | m    |           |   | UL Style             | Flam. Rating | Voltage Rating | Voltage Actual | Temp. Rating (C) | Dia. mm (inch) | # of Cond | Size AWG |                       |
| Fluoro   | 75                      | 22.9 | 2403438-3 | 5 BNC MALE TO HD 15 MALE 75 FEET                            | 1015                 | FT4          |                | 1Vp-p          | 75               | 9.1 (0.358)    | 5         | 26       |                       |
|          | 70                      | 21.2 | 2213219   | POWER CABLE FOR LCD-CONSOLE TO LCD                          |                      | FT1          | 120            | 120VAC         | 105              | 9.3 (0.366)    | 3         | 14       |                       |
|          | 71                      | 21.5 | 5116219   | Grounding Cable For LCD Console To LCD                      |                      | VW-1         | 600            | 0V             | 105              |                | 1         | 8        |                       |
|          | 15                      | 4.6  | 2403438-4 | HD 15 FEMALE TO HD 15 MALE 15 FEET                          |                      |              |                | 1Vp-p          | 60               | 8.0 (0.315)    | 5         | 26       |                       |
|          | 1.3                     | 0.4  | 2355108   | JUMPER CABLE FOR ADAPTING 2286150 TO WORK WITH H-POWER MSUB |                      |              |                |                |                  |                | 8         | 22       |                       |
| UPS      | 15                      | 4.6  | 2391751   | POWER CABLE, NGPDU TO UPS                                   | 2587                 | FT4          | 600            | 208VAC         | 90               | 5.8 (0.228)    | 5         | 8        |                       |
|          | 15                      | 4.6  | 2391751-3 | POWER CABLE, UPS DISCONNECT PANEL TO NGPDU                  | 2587                 | FT4          | 600            | 208VAC         | 90               | 5.8 (0.228)    | 4         | 8        |                       |
|          | 45                      | 13.6 | 5169224   | UPS CONTROL CABLE   | 2587                 | FT4          | 600            | 120VAC         | 90               | 10.3 (0.406)   | 5         | 18       |                       |
| Injector | 100                     | 30.5 | 5169456   | GANTRY TO INJECTOR  | 1007                 | VW-1         | 300            | <30VDC         | 80               | 1.57 (0.062)   | 3         | 22       | 45 (1.78) Dia         |
|          | 8.2                     | 2.5  | 5317258   | POWER CABLE INJECTOR TO CONSOLE                             | 62                   | VW-1         | 300            | 120VAC         | 60               | 9.4 (0.37)     | 3         | 14       | 36 (1.41) Dia         |

Table 8-10 GE Supplied Cables for Options - UL Information


| OPTION      | LENGTH,<br>ACTUAL<br>(USABLE) |      | PART #  | DESCRIPTION           | UL CABLE INFORMATION |                                 |                   |                   |                     |                   |           |             | PULL SIZE<br>MM (INCHES) |
|-------------|-------------------------------|------|---------|-----------------------|----------------------|---------------------------------|-------------------|-------------------|---------------------|-------------------|-----------|-------------|--------------------------|
|             | ft                            | m    |         |                       | UL Style             | Flam.<br>Rating                 | Voltage<br>Rating | Voltage<br>Actual | Temp.<br>Rating (C) | Dia.<br>mm (inch) | # of Cond | Size<br>AWG |                          |
| Cardiac     | 30                            | 9.1  | 5198566 | GANTRY TO EKG MONITOR | 2919                 | UL<br>1685<br><br>UL<br>loading | 30                | <30VDC            | 80                  | 6.45<br>(0.254)   | 6         | 24          | 37 (1.45) Dia            |
| Adv 4D Resp | 100                           | 30.5 | 5199717 | GANTRY TO RPM UNIT    | 2464                 | FT4                             | 300               | <30VDC            | 80                  | 6.6 (0.26)        | 4         | 22          | 37 (1.45) Dia            |
|             |                               |      |         |                       |                      |                                 |                   |                   |                     |                   |           |             |                          |
|             |                               |      |         |                       |                      |                                 |                   |                   |                     |                   |           |             |                          |

Table 8-10 GE Supplied Cables for Options - UL Information

### 3.5 Contractor (Customer) Supplied

| CUSTOMER INSTALLED WIRING  |  | DESCRIPTION                                | CABLES SUPPLIED |                 |              | PLUG PULLING DIMENSIONS |    | WIRE & CABLE PIGTAILS FT. (M.) |       |
|--|--|--|-----------------|-----------------|--------------|-------------------------|----|--------------------------------|-------|
| QTY  | SIZE AWG (MM <sup>2</sup> )  |  | PART NO         | LENGTH FT. (M.) | DIA. IN (MM) | FROM                    | TO | FROM                           | TO    |
| <b>RUN NO. 1</b> FROM PRIMARY POWER SOURCE TO FACILITY DISCONNECT (POWER SOURCE - A1)    |  |  |                 |                 |              |                         |    |                                |       |
| Maximum Run Length *   |  |  |                 |                 |              |                         |    |                                |       |
| 3  | *  | POWER                                      |                 |                 |              |                         |    | 3 (1)                          | 3(1)  |
| 1  | 1/0 (50)   | GROUND                                     |                 |                 |              |                         |    | 3 (1)                          | 3 (1) |
| <b>RUN NO. 2</b> FROM FACILITY DISCONNECT TO POWER MODULE (A1 - PM) MAXIMUM RUN LENGTH * |  |  |                 |                 |              |                         |    |                                |       |
| 3  | *  | POWER                                      |                 |                 |              |                         |    | 3 (1)                          | 3(1)  |
| 1  | 1/0 (50)   | GROUND                                     |                 |                 |              |                         |    | 3 (1)                          | 3 (1) |
| 1  | *  | NEUTRAL                                    |                 |                 |              |                         |    | 3 (1)                          | 3 (1) |
| <b>RUN NO. 3</b> FROM FACILITY DISCONNECT TO SYSTEM EMERGENCY OFF (A1 - SEO)             |  |  |                 |                 |              |                         |    |                                |       |
| 2  | 14 (2)   | POWER                                      |                 |                 |              |                         |    | 6 (2)                          | 6 (2) |
| 1  | 14 (2)   | GROUND                                     |                 |                 |              |                         |    | 6 (2)                          | 6 (2) |
| <b>RUN NO. 4</b> POWER MODULE TO WARNING LIGHT CONTROL (PM - WL)                         |  |  |                 |                 |              |                         |    |                                |       |
| 2  | 14 (2)   | WARNING LIGHT 24 VOLT CONTROL A3J2-1,2,3,4 |                 |                 |              |                         |    |                                |       |
| <b>RUN NO. 5</b> POWER MODULE TO SCAN ROOM DOOR INTERLOCK (PM - DOOR SWITCH)             |  |  |                 |                 |              |                         |    |                                |       |
| 2  | 14 (2)   | SCAN ROOM DOOR INTER LOCK A3J6-1,2         |                 |                 |              |                         |    |                                |       |
| *  | REFER TO <a href="#">Table 7-2 on page 119</a> FOR AWG (MM2) WIRE SIZES. |  |                 |                 |              |                         |    |                                |       |

Table 8-11 Runs 1, 2, 3, 4 and 5 Connections




NOTES :

- 1) Used for remote diagnostics - Option
- 2) Refer to the appropriate Pre-installation / Installation documents for the Laser Camera
- 3) Category 5 cable. Use one of the following patch cords:

| CAT Num | GE Part Num | Length |
|---------|-------------|--------|
| K9000WB | 2215028-10  | 20 m   |
| K9000KP | 2215028-5   | 10 m   |
| K9000JR | 2215028-4   | 5 m    |
| K9000WA | 2215028-9   | 3 m    |
- 4) In order to avoid any violation of each National Regulation (NEC in USA, CCC in China, etc.), use of the complied cable/wire is recommended. For China market, China end-user shall purchase the power supply cable that has the CCC mark.

Only one phone connection is required for the system.



KEY:  
[ Dashed Box ] OPTION

Figure 8-4 Interconnection Runs

## Section 4.0

### Contractor Supplied Components

| REFERENCE                 | ASSOCIATED EQUIPMENT                    | MATERIAL/LABOR SUPPLIED BY CUSTOMER CONTRACTOR  | VENDOR / CAT NO. GE CATALOG  |
|---------------------------|---|---|--|
| A1<br>380-480V<br>50/60Hz | Circuit Breaker with Magnetic Contactor | 3 Pole, 380V-480V, Combination breaker with magnetic contactor. Includes control transformer, optional UPS interface, On/Off controls and auto-restart feature. Must be LOTO compatible.                          | Recommend*:<br><ul style="list-style-type: none"> <li>E4502AE (125A)</li> <li>E4502AF (150A)</li> <li>GEMS_CT_330 (European PDB)</li> </ul> Optional remote operator control available from GE Supply, Cat# GESCTR0CS1 |
| BBNC (required)           | Broad-Band Network Connection           | Broad-Band network connection wall jack, located within 1m (39inches) of console location, for internal hospital networking and InSite Broad-Band connectivity.<br>Cabling to conform to facility's IT standards. |  |
| ITL (optional)            | In-suite Telephone Lines                | Supply 2 voice-grade telephone lines. One line must be a direct number from outside the facility – do not route this line through a telephone switchboard. Telephone line operating charges are paid by customer. |  |
|                           | System Components                       | Reference the system installation drawings supplied by Installation Support Services within your geographic area.   |  |

\*Refer to [Section 3.0 - Recommended Power Distribution System](#) .

**Table 8-12 Contractor-Supplied Components**



## Section 5.0

### Fuse

| MARK | NUMBER        | QTY | FRU | DESCRIPTION/NAME          |
|------|---------------|-----|-----|---------------------------|
| 1    | 2351492       | 1.0 | Yes | 200A FUSE                 |
| 2    | 2364059       | 1.0 | Yes | GLASS FUSE                |
| 3    | 46-170021P50  | 1.0 | Yes | FUSE 12 AMPS 250 VOLTS    |
| 4    | 46-170021P106 | 1.0 | Yes | FUSE 8 AMPS 250 VOLTS     |
| 5    | 2336517-2     | 1.0 | Yes | FUSE 25 AMPS 700 VOLTS 2. |
| 6    | 46-170021P52  | 5.0 | Yes | FUSE 3 AMPS 250 VOLTS     |
| 7    | 46-170021P10  | 2.0 | Yes | 006.000A 0250V 3AG        |
| 8    | 46-170021P101 | 1.0 | Yes | FUSE 20 AMPS 700 VOLTS    |
| 9    | 46-170021P11  | 1.0 | Yes | 7.0A, 250V SLO-BLO FUSE.  |
| 10   | 46-170021P31  | 2.0 | Yes | 1/2A, 250V SLO-BLO FUSE.  |
| 11   | 2106993-5     | 2.0 | Yes | FUSE 20 AMPS 500 VOLTS .3 |
| 12   | 46-327160P1   | 1.0 | Yes | 12.0A, 125VAC, DUAL ELEME |
| 13   | 2379651       | 1.0 | Yes | FUSE - 700 V, 200 A, FAST |
| 14   | 2374694       | 2.0 | Yes | 8A TIME DELAY FUSES       |


Table 8-13 Freedom RT Fuse Kit (2218570-3 BOM, rev 1)

| MARK | NUMBER       | QTY | FRU | DESCRIPTION/NAME  |
|------|--------------|-----|-----|---|
| 1    | 2364059      | 1.0 | Yes | GLASS FUSE  |
| 2    | 46-170021P50 | 1.0 | Yes | FUSE 12 AMPS 250 VOLTS                                      |
| 3    | 5325319      | 1.0 | Yes | FUSE 25 AMPS 700 VOLTS 2.5M SECONDS                         |
| 4    | 46-170021P52 | 2.0 | Yes | FUSE 3 AMPS 250 VOLTS                                       |
| 5    | 46-170021P31 | 2.0 | Yes | 1/2A, 250V SLO-BLO FUSE.                                    |
| 6    | 5329469      | 2.0 | Yes | FUSE 20 AMPS 500 VOLTS .3 SECONDS--RoHS                     |
| 7    | 2379651      | 1.0 | Yes | FUSE - 700 V, 200 A, FAST SEMI                              |
| 8    | 2374694      | 2.0 | Yes | 8A TIME DELAY FUSES   |
| 9    | 5324765      | 1.0 | Yes | 200 A FUSE  |
| 10   | 99170530     | 1.0 | Yes | FUS T 10X38 8A 500V C 10KA ULL                              |
| 11   | 5368105      | 1.0 | Yes | 1A Time_Delay fuse,10_4X38_1CLIP                            |
| 12   | 5306477-3    | 2.0 | Yes | PROTECTION FUSE, 15A, TIME DELAY, 300VDC, 600VAC, CARTRIDGE |
| 13   | 5327448-2    | 2.0 | Yes | PROTECTION FUSE, 2A, NA, FAST ACTION, 250V, HM FUSE,        |
| 14   | 5327449      | 2.0 | Yes | PROTECTION FUSE, 3.2A, NA, FAST ACTION, 250V, HM FUSE       |
| 15   | 5306477-5    | 1.0 | Yes | PROTECTION FUSE, 12A, TIME DELAY, 300VDC, 600VAC, CARTRIDGE |

Table 8-14 RT Innovation Fuse Kit (2218570-4 BOM, rev 3)

## Section 6.0

### UPS Interconnect



**Figure 8-5 Typical UPS Interconnect**

The UPS Kit catalog number is B7999ZA. See table below, for detail.



**NOTICE** UPS Kit B7999ZA REQUIRES installation of one of the A1 Panels listed below.

| PDU Type & Model # | Max. Mom. kVA Rating | Recommended Main Disconnect (A1) Cat #  |   | Optional Partial UPS Kit Cat #   |
|--------------------|----------------------|---|---|--|
|                    |                      | Europe & Asia (380-400V or 420V)  | North America (440V or 460-480V)                                    |  |
| NGPDU 2326492-61   | 150kVA               | <b>E4502AF</b> (150A) (incl. Auto Restart & Integrated UPS Control)<br>GEMS_CT_330 (European PDB) | <b>E4502AE</b> (125A) (incl. Auto Restart & Integrated UPS Control) | <b>B7999ZA</b><br><b>alt. E4502KY</b><br>(Includes PW 9155-10GE10kVa, 2ph. UPS, & hardware kit)<br><b>REQUIRES</b> one of the A1 Panels shown at left, or equal. |

**Table 8-15 Partial UPS Back-up Options**

Conduit is required between:

- A1 and UPS
- UPS and PDU
- PDU and A1

## Section 7.0

### Typical Customer Supplied Wiring


#### 7.1 Primary Power Disconnect - U.S.



**Figure 8-6 Primary Power Disconnect (A1) - Typical GE A1 shown**

Installing this system requires a Lockout/Tagout-compatible disconnect. If a UPS is required, a GE Disconnect is strongly recommended for safe operation. The GE disconnect and UPS are designed to work together.

7.2 Scan Room Warning Light & Door Interlock



Note: This is GE Brand Room Light Shown

Figure 8-7 TS6 X-Ray Warning Light Connections





Figure 8-8 TS6 Room Door Interlock Connections - With a Door Interlock



If not using a door switch, add a jumper.  
If jumper is not in place, exposures will not be made. Check this jumper if you get scan interlock errors.

Figure 8-9 TS6 Room Door Interlock Connections - Without a Door Interlock



[www.gehealthcare.com](http://www.gehealthcare.com)

