

GE Healthcare

Australia & New Zealand

# Clinical Consumables



[www.gehealthcare.com.au](http://www.gehealthcare.com.au)

# How to order

---

It's simple – find the item you want, and simply fax or email us at the contact details below.

Our friendly customer support team will look after you every step of the way.

## Ordering Information

### Australia

Email: [AUCS@ge.com](mailto:AUCS@ge.com)

Fax: 02 9846 4762

Phone: 1800 223 615

### New Zealand

Email: [NZCS@ge.com](mailto:NZCS@ge.com)

Fax: 09 262 3898

Phone: 0800 122 277

## MyServices

---

**A simpler way for biomedical engineers to access the information they need.**

MyServices lets biomedical engineers save valuable time by providing online technical information to more effectively support clinical equipment.

MyServices web portal can help get your team off the phone and back on the job.

Available 24/7, MyServices gives you the digital tools to enhance your technical maintenance & support process. We can help you search for parts and access online training and technical documentation quickly.

MyServices is available via any web-connected device and will continually evolve to more effectively support biomedical engineering departments.

<https://services.gehealthcare.com.au>



### Track my Repair

Check the repair status of equipment shipped to the national repair depot.



### Learn

Develop technical competency by accessing a range of online training videos & technical documentation, reducing repair times and improving workflow



### Finding the Right Part

Identify parts quickly and accurately using the part identification tool and reduce part ordering errors and delays



### User Access

Manage user access by assigning various roles and privileges for individual team members

NIBP Supplies

SpO<sub>2</sub> Supplies

Anaesthesia  
Supplies

Anaesthetic Breathing Circuit Spares .....	34
Ancillaries .....	36
BP Cuff Adapters .....	13
BP Cuff Adapters DINACLICK and Neo-Snap .....	16
BP Cuff Connectors .....	17
Calibration Gases & Regulators .....	33
CLASSIC-CUF .....	8
Connector Cables Adult-Paediatric .....	11
Connector Cables Neonatal .....	12
Disposable SpO <sub>2</sub> Sensors .....	27
DURA-CUF .....	6
Entropy Cables & Sensors .....	32
Flow Sensor Spares .....	34
Gas Monitoring .....	33
MASIMO Cables and Probes .....	22
Masimo LNCS Reusable & Adhesive Sensors .....	21
Masimo LNOP Adhesive Sensors .....	20
Masimo LNOP Reusable Sensors .....	19
NELLCOR Cables & Probes .....	23
Neonatal Cuffs .....	7
Neuromuscular Transmission Sensors & Cables .....	31
Radial-Cuf .....	5
Reusable SpO <sub>2</sub> Sensors .....	26
SENSA-CUF .....	10
SOFT-CUF .....	5
Spirometry .....	35
SpO <sub>2</sub> Accessories .....	29
SpO <sub>2</sub> Interconnect cables .....	26
SpO <sub>2</sub> Sensors - Datex Connector .....	29
SpO <sub>2</sub> Sensors - GE Connector .....	28
SpO <sub>2</sub> Sensors - Ohmeda Connector .....	29
SpO <sub>2</sub> Sensors - TruSat Connector .....	28
Wash Plugs and O-Rings .....	18
Yellow CLASSIC-Cuf .....	9

# NIBP Supplies

## Performance accuracy

### Quality cuffs for every patient, in every care setting

A patient's blood pressure (BP) is often the first step in determining the best clinical course of action. That's why quality data acquisition through your BP cuff is so essential.

GE CRITIKON® BP cuffs from GE Healthcare are designed and tested for performance accuracy – from the first touch to the final measurement. As part of an end-to-end system for managing Clinical Information Logistics™ across your hospital, GE CRITIKON BP cuffs help you gather precise data up front. By starting with accurate information, you can make well-informed decisions that may ultimately enhance patient outcomes.

#### SOFT-CUF

Our most popular cuff, the CRITIKON SOFT-CUF® is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer-term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It can be cleaned and reused on multiple patients, but is priced economically so it can be disposed of after a single patient use, helping you affordably prevent infections.

#### DURA-CUF

CRITIKON DURA-CUF® is gentle and designed for durability. Made with tough, flexible nylon, they are convenient and easy-to-clean with curved corners to enhance your patient's comfort. They are constructed for multiple patient use in low-acuity care settings and ideal for physicians' offices where the risk of contamination is diminished.

#### Neonatal Cuffs

Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl. The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient use.



#### RADIAL-CUF

GE Healthcare offers the first BP forearm cuff clinically validated for use on obese or overly muscular adult patients. It provides proper fit and comfort to patients, while giving reliable blood pressure readings that support quality clinical care.

#### Yellow CRITIKON CLASSIC-CUF

The yellow CRITIKON CLASSIC-CUF® is designed for single patient use within an isolation environment. These vinyl cuffs are priced as disposable-use cuffs and support your infection prevention protocols. They provide the same quality and assurance that you expect from the entire line of CRITIKON cuffs, while tailored for isolation areas.

#### CLASSIC-CUF

CRITIKON CLASSIC-CUF® BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.

#### SENSA-CUF

The CRITIKON SENSA-CUF® is a highly versatile solution for a wide range of extended-care settings. Made using strong, breathable, easy-to-clean material, the SENSA-CUF is designed to maximize patient comfort in low-acuity areas in a hospital where risk of contamination are diminished, and durable enough for multi-patient use.

Single Patient  
short-term use









Single  
Patient use

**Limited  
reuse**

Hospital  
reusable




General care  
reusable

## SOFT-CUF

SOFT-CUF			GE Healthcare				Welch Allyn	Colin Philips/HP Spacelabs Siemens Datascop	Protocol	Manual Sphyg. Wall Manometer
Size	Colour	Limb Circumference								
Infant	Orange	8-13 cm	2500	SFT-P1-2A	2401	2450 <sup>1</sup>	002331	2360	2320	
Child	Green	12-19 cm	2501	SFT-P2-2A	2402	2451	002332	2361	2321	002584
Child Long	Green	12-19 cm	2506	SFT-P2-2A-L	2400					
Small Adult	Light Blue	17-25 cm	2502	SFT-A1-2A	2403	2452	002333	2362	2322	002583
Small Adult Long	Light Blue	17-25 cm	2507	SFT-A1-2A-L	2407		002334			
Adult	Navy	23-33 cm	2503	SFT-A2-2A	2404	2453	002335	2323	2323	002580
Adult Long	Navy	23-33 cm	2604	SFT-A2-2A-L	2116	2454	002336			002582
Large Adult	Rose	31-40 cm	2504	SFT-A3-2A	2405	2455	002337	2324	2324	002581
Large Adult Long	Rose	31-40 cm	2644	SFT-A3-2A-L	2117	2456	002338			002585
Thigh	Brown	38-50 cm	2505	SFT-T1-2A	2406	2457	002339	2325	2325	

1. Infant BP Cuff for Datex-Ohmeda monitors has Male Slip Luer Connector.

## Radial-Cuf

RADIAL-CUF		GE Healthcare			
Colour	Forearm Circumference				
Purple	26-36 cm	SFT-F1-2E	SFT-F1-2A	SFT-F1-2C	SFT-F1-2D

# NIBP Supplies

Single Patient short-term use









Single Patient use

Limited reuse

Hospital reusable

General care reusable

## DURA-CUF

DURA-CUF			GE Healthcare				Non-GE Healthcare			
										
Size	Colour	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector	1-Tube Submin Connector	DURA-CUF Inflation System
Infant	Orange	8-13 cm	002783	DUR-P1-2A	002200	2750 <sup>1</sup>	002760	002274	2740	002385
Child	Green	12-19 cm	002781	DUR-P2-2A	002201	2751	002761	002275	2741	002384
Child Long	Green	12-19 cm	2785	DUR-P2-2A-L			002769			002386
Small Adult	Royal Blue	17-25 cm	002779	DUR-A1-2A	002202	2752	002762	002276	2742	002383
Small Adult Long	Royal Blue	17-25 cm		DUR-A1-2A-L			002763			002387
Adult w/ Hanger	Navy	23-33 cm	002771							002380
Adult	Navy	23-33 cm	002774	DUR-A2-2A	002203	2753	002764	002277	2743	
Adult Long	Navy	23-33 cm	002772	DUR-A2-2A-L	002206	002756	002765	002130		2388
Large Adult	Wine	31-40 cm	002791	DUR-A3-2A	002204	2754	002766	002278	2744	002381
Large Adult Long	Wine	31-40 cm	002784	DUR-A3-2A-L	002207	002757	002767	002131		2389
Thigh	Brown	38-50 cm	002796	DUR-T1-2A	002205	2755	002768	002279	2745	002382

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

DURA-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	002699	1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh (6 cuffs/box)
	002697	2 Infant, 3 Child, 1 Small Adult (6 cuffs/box)
	002698	1 Small Adult, 2 Adult, 1 Adult Long, 1 Large Adult, 1 Large Adult long (6 cuffs/box)
2-Tube Submin Connectors	2299	1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh (6 cuffs/box)
2-Tube DINACLICK Connectors	2059301-001	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)

Single Patient  
short-term use

Single  
Patient use

**Limited  
reuse**

Hospital  
reusable

General care  
reusable

## Neonatal Cuffs

Neonatal			GE Healthcare				Philips/HP Siemens SpaceLabs Datascope			
			GE CRITIKON GE DINAMAP CARESCAPE GE Marquette							
			2-Tube Luer Connector		2-Tube Neo-Snap Connector		1-Tube Luer Connector		1-Tube Neo-Snap Connector	
Size	Color	Limb Circumference	SOFT-CUF	CLASSIC- CUF	SOFT-CUF	CLASSIC- CUF	SOFT-CUF	CLASSIC- CUF	SOFT-CUF	CLASSIC- CUF
1	Orange	3-6 cm	2521	2638	SFT-N1-2B	CLA-N1-2B	2121	2138	SFT-N1-1B	CLA-N1-1B
2	Lt. Blue	4-8 cm	2422	2633	SFT-N2-2B	CLA-N2-2B	2122	2133	SFT-N2-1B	CLA-N2-1B
3	Green	6-11 cm	2523	2628	SFT-N3-2B	CLA-N3-2B	2124	2128	SFT-N3-1B	CLA-N3-1B
4	Navy	7-13 cm	2524	2623	SFT-N4-2B	CLA-N4-2B	2125	2123	SFT-N4-1B	CLA-N4-1B
5	Rose	8-15 cm	2525	2619	SFT-N5-2B	CLA-N5-2B	2126	2119	SFT-N5-1B	CLA-N5-1B









Neonatal Assortment Pack	Cuff Type	Part Number	Description
2-Tube Luer Connectors	SOFT-CUF	2694	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2693	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
1-Tube Luer Connectors	SOFT-CUF	2696	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2680	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)



# NIBP Supplies



## CLASSIC-CUF


CLASSIC-CUF			GE Healthcare				Non-GE Healthcare			
Size	Colour	Limb Circumference								
Infant	Orange	8-13 cm	2618	CLA-P1-2A	2354	2170	002160	2300	2150	
Child	Green	12-19 cm	2613	CLA-P2-2A	2355	2171	002161	2301	2151	002287
Small Adult	Light Blue	17-25 cm	2608	CLA-A1-2A	2356	2172	002162	2302	2152	002285
Adult	Navy	23-33 cm	2603	CLA-A2-2A	2357	2173	002164	2303	2153	002283
Adult Long	Navy	23-33 cm	2647	CLA-A2-2A-L	2352	2174	002165	2306		
Large Adult	Rose	31-40 cm	2643	CLA-A3-2A	2358	2175	002166	2304	2154	002281
Large Adult Long	Rose	31-40 cm	2649	CLA-A3-2A-L	2353	2176	002167	2307		
Thigh	Brown	38-50 cm	2648	CLA-T1-2A	2359	2177	002168	2305	2155	

1. Infant BP Cuff for Datex-Ohmeda monitors has Male Slip Luer Connector.

Note: 20 cuffs/box

CLASSIC-CUF Assortment Pack	Cuff Type	Description
2-Tube Screw Connectors	2056305-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2692	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
2-Tube Submin Connectors	2056305-002	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2292	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
2-Tube DINACLICK Connectors	2059304-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)





## Yellow CLASSIC-Cuf

YELLOW CLASSIC-CUF			GE Healthcare		
					
Size	Colour	Forearm Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Conector
Infant	Yellow/Blue	8 - 13 cm	2650	CLA-P1-2A-X	2670
Child	Yellow/Blue	12 - 19 cm	2650	CLA-P2-2A-X	2671
Small Adult	Yellow/Blue	17 - 25 cm	2607	CLA-A1-2A-X	2672
Adult	Yellow/Blue	23 - 33 cm	2602	CLA-A2-2A-X	2673
Large Adult	Yellow/Blue	31 - 40 cm	2642	CLA-A3-2A-X	2674
Thigh	Yellow/Blue	38 - 50 cm	2652	CLA-T1-2A-X	2675

Note: 20 cuffs/box



# NIBP Supplies

Single Patient  
short-term use







Single  
Patient use

Limited  
reuse

**Hospital  
reusable**

General care  
resuable

## SENSA-CUF

SENSA-CUF			GE Healthcare				Non-GE Healthcare	
								
Size	Colour	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	2458	SEN-P1-2A	2482	2430	2440	2470
Child	Green	12-19 cm	2460	SEN-P2-2A	2484	2431	2441	2472
Small Adult	Lt Blue	17-25 cm	2462	SEN-A1-2A	2486	2432	2442	2474
Small Adult Long	Lt Blue	17-25 cm	2463	SEN-A1-2A-L	2487	2433	2443	2475
Adult	Navy	23-33 cm	2464	SEN-A2-2A	2488	2434	2444	2476
Adult Long	Navy	23-33 cm	2465	SEN-A2-2A-L	2489	2435	2445	2477
Large Adult	Rose	31-40 cm	2466	SEN-A3-2A	2490	2436	2446	2478
Large Adult Long	Rose	31-40 cm	2467	SEN-A3-2A-L	2491	2437	2447	2479
Thigh	Brown	38-50 cm	2468	SEN-T1-2A	2492	2438	2448	2480

Note: 5 cuffs/box

SENSA-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	2494	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)
2-Tube DINACLICK Connectors	2059300-001	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)



## Connector Cables Adult-Paediatric



### Adult-Pediatric, DINACLICK Connector

Compatibility: Eagle 3000, Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

Length: 1.2 m/4 ft      Part # 2058203-001  
 Length: 3.6 m/12 ft      Part # 2058203-002  
 Length: 7.2 m/24 ft      Part # 2058203-003



### Adult-Paediatric, Screw Connectors

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

Length: 1.2 m/4 ft      Part # 2017008-004



### Adult-Paediatric, Subminiature Connectors

Compatibility: Eagle 3000, Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, Corometrics 120, 128, 129, 259, 259cx, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

Length: 3.6 m/12 ft      Part # 2017008-001  
 Length: 1.2 m/4 ft      Part # 2017008-002



### Adult, Mated Subminiature Connectors

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, Corometrics 120, 128, 129, 259, 259cx, B30 Patient Monitor, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

Length: 3.6 m/12 ft      Part # 2020980-001



### Adult-Paediatric, Screw Connectors

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, Corometrics 120, 128, 129, 259, 259cx, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

Length: 3.6 m/12 ft      Part # 2017008-003



### Adult, Mated Subminiature Connectors

Compatibility: MRI Monitor


Length: 6 m/19.7 ft      Part # 2020980-002

NOTE: The NIBP measurement in E-PSM(P), E-PRESTN, E-PRETN, E-RSTN and MRI Monitor is intended for patients weighing over 5 kg (11 lb).

# NIBP Supplies

---

## Connector Cables Neonatal



### Neo-Snap Connector

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module

---

Length: 1.2 m/4 ft                      Part # 2058204-001

Length: 2.4 m/8 ft                      Part # 2058204-002

Length: 3.6 m/12 ft                     Part # 2058204-003

Length: 6 m/19.7 ft                     Part # 2058204-004



### Neonatal, Female Slip Luer

Compatibility: MRI Monitor

---

Length: 6 m/19.7 ft                      Part # 2017009-004



### Neonatal, Female Slip Luer

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN, B30 Patient Monitor

---

Length: 2.4 m/8 ft                      Part # 2017009-001



### Air Hose Adapter

Compatibility: connect round connector air hoses 9461-203, 9461-205, 9461-214, 9461-215 to rectangular connector NIBP port on monitor

---

Length: 0.25 m/0.83 ft                      Part # 414876-001



### Neonatal, Female Slip Luer

Compatibility: Dash 2000, 2500, 3000, 4000, 5000, Tram, CARESCAPE Patient Data Module, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN

---

Length: 1.2 m/4 ft                      Part # 2017009-002

Length: 3.6 m/12 ft                     Part # 2017009-003

## BP Cuff Adapters

GE's connection system is simple and easy to implement.

GE Healthcare provides cost-effective standardization and inventory management solutions. By standardizing blood pressure (BP) connection systems, you create a universal connection to allow the same type of BP cuff connection to work on a wide variety of non-invasive blood pressure devices. By consolidating all of your BP cuff purchases to one vendor, your Purchasing Department will have fewer purchase orders and invoices. And a single point of contact for BP cuffs is created, providing you with a dedicated implementation team.

Products included in this category are not made with natural rubber latex.

### Connecting BP Cuff to 2 Male Slip Luers



**Y, Unfitted to 2 Male Slip Luers**

Compatibility: Unfitted

Part # 330062 (Qty: 10)



**Y, 1 Male Bayonet to 2 Male Slip Luers**

Compatibility: Philips/HP®, Siemens®, Datascope®, Colin®

Part # 330071 (Qty: 10)



**Y, 1 Female Locking Luer to 2 Male Slip Luers**

Compatibility: Spacelabs®, Datascope®, Colin®

Part # 330065 (Qty: 10)



**Y, 1 Female Submin to 2 Male Slip Luers**

Compatibility: Protocol®

Part # 330066 (Qty: 10)



**X, 2 Female Submins to 2 Male Slip Luers**

Compatibility: GE, GE Marquette\*

Part # 300939 (Qty: 10)



**Y, 1 Male Submin to 2 Male Slip Luers**

Compatibility: Nihon Kohden®

Part # 330069 (Qty: 10)



**X, 2 Mated Pair Submins to 2 Male Slip Luers**

Compatibility: GE Datex-Ohmeda\*

Part # 330063 (Qty: 10)



**X, 2 Male Bayonets to 2 Male Slip Luers**

Compatibility: User Modified

Part # 330067 (Qty: 10)

# NIBP Supplies

---

## Connecting BP Cuff to Female Slip Luer



### **Y, 1 Male Slip Luer to 2 Female Slip Luers - Neonatal**

Compatibility: Philips/HP,  
Siemens, Datascope, Spacelabs

Part # 330049 (Qty: 1)



### **1 Male Micro Plug to 1 Female Locking Luer - Neonatal**

Compatibility: Spacelabs

Part # 330088 (Qty: 10)



### **Y, 1 Male Micro Plug to 2 Female Locking Luers - Neonatal**

Compatibility: Spacelabs

Part # 330078 (Qty: 10)

## Connecting BP Cuff to 2 Female Screw Connectors



### **Y, 1 Unfitted to 2 Female Screw Connectors**

Compatibility: Unfitted

Part # 330048 (Qty: 1)



### **Y, 1 Male Bayonet to 2 Female Screw Connectors**

Compatibility: Philips/HP,  
Siemens, Datascope, Colin

Part # 330046 (Qty: 1)

## Connecting BP Cuff to 2 Unfitted



### **Y, 1 Unfitted to 2 Unfitted**

Compatibility: Unfitted

Part # 300829 (Qty: 1)



## Connecting BP Cuff to 2 Male Submins



**Y, 1 Unfitted to 2 Male Submins**  
Compatibility: Unfitted  
Part # 330084 (Qty: 10)



**Y, 1 Male Bayonet to 2 Male Submins**  
Compatibility: Philips/HP, Siemens, Datascope, Spacelabs  
Part # 330080 (Qty: 10)



**Y, 1 Female Locking Luer to 2 Male Submins**  
Compatibility: Legacy Spacelabs, Datascope, Colin  
Part # 330081 (Qty: 10)



**Y, 1 Female Submin to 2 Male Submins**  
Compatibility: Protocol  
Part # 330093 (Qty: 10)



**Y, 1 Male Submin to 2 Male Submins**  
Compatibility: Nihon Kohden  
Part # 330083 (Qty: 10)



**X, 2 Male Screw Connectors to 2 Male Submins**  
Compatibility: GE DINAMAP\*  
Part # 330087 (Qty: 10)



**X, 2 Male Bayonets to 2 Male Submins**  
Compatibility: MDE, Spacelabs  
Part # 330086 (Qty: 10)



**X, Mated Pair Submins to 2 Male Submins**  
Compatibility: GE Datex-Ohmeda  
Part # 330085 (Qty: 10)



# NIBP Supplies

## BP Cuff Adapters DINACLICK and Neo-Snap

Connecting BP Cuff to NEW DINACLICK and Neo-Snap Connectors



**Y, 1 Male Bayonet to DINACLICK Female Recept**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 2059291-001 (Qty: 10)



**Y, 1 Female Submin to DINACLICK Female Recept**

Compatibility: Protocol

Part # 2059291-002 (Qty: 10)



**Y, T Male Submin to DINACLICK Female Recept**

Compatibility: Nihon Kohden

Part # 2059291-003 (Qty: 10)



**Y, 1 Unfitted to DINACLICK Female Recept**

Compatibility: Unfitted

Part # 2059291-004 (Qty: 10)



**X, 2 Male Bayonet to DINACLICK Female Recept**

Compatibility: User Modified

Part # 2059292-001 (Qty: 10)



**X, 2 Female Submin to DINACLICK Female Recept**

Compatibility: GE Marquette

Part # 2059292-002 (Qty: 10)



**X, 2 Mated Submin to DINACLICK Female Recept**

Compatibility: GE Datex-Ohmeda

Part # 2059292-003 (Qty: 10)



**X, 2 Male Screw to DINACLICK Female Recept**

Compatibility: DINAMAP

Part # 2059292-004 (Qty: 10)



**Y, 1 Male Slip Luer to Neo-Snap\* Male Plug**

Compatibility: Philips/HP, Siemens, Datascope, Spacelabs

Part # 2059290-001 (Qty: 10)



**Male Micro Plug to 2 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059290-002 (Qty: 10)



**Male Micro Plug to 1 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059297-001 (Qty: 10)



## BP Cuff Connectors



**Tubing, Gray, Two Inches Long, 5/32" ID Unfitted**

Compatibility: Unfitted

Part # 300959 (Qty: 10)



**5/32" ID Barbed to 1/8" ID Barbed, Plastic**

Compatibility: Airhose Connector

Part # 300671 (Qty: 10)



**5/32" ID Barbed to 5/32" ID Barbed, Plastic**

Compatibility: Airhose Connector

Part # 300670 (Qty: 10)



**Female Submini to 5/32" ID Barbed, Plastic**

Compatibility: GE, GE Marquette, GE Datex-Ohmeda

Part # 330092 (Qty: 10)



**5/32" ID Barbed to Male Slip Luer, Plastic**

Compatibility: Female Locking Luer and Screw Connectors

Part # 300666 (Qty: 10)



**Female Locking Luer to 5/32" ID Barbed, Plastic**

Compatibility: Spacelabs, Datascope, Colin

Part # 300668 (Qty: 10)



**Male Screw to Male Slip Luer, Plastic**

Compatibility: Female Locking Luer and Screw Connectors

Part # 300662 (Qty: 10)



**Male Screw to 5/32" ID Barbed, Plastic**

Compatibility: GE DINAMAP

Part # 300664 (Qty: 10)



**5/32" ID Barbed to Female Bayonet, Metal**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 330060 (Qty: 10)



**Male Bayonet to 3/16" ID Barbed, Metal**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 330059 (Qty: 10)



**2T DINACLICK Female Recept to 1/8" ID tubing, dual lumen hoses**

Compatibility: DINACLICK Male Plug

Part # 2059294-001 (Qty: 10)



**2-1T DINACLICK Female Recept for single lumen hoses**

Compatibility: DINACLICK Male Plug

Part # 2055630-001 1/8" tube (Qty: 10)

Part # 2055630-002 5/32" tube (Qty: 10)



**5/32" ID Barbed to Female Slip Luer**

Compatibility: Universal Neonatal

Part # 330068 (Qty: 10)



**5/32" ID Barbed to Female Legacy Marquette Connector, Plastic**

Compatibility: GE Marquette

Part # 330057 (Qty: 10)



**2T Neo-Snap Male Plug**

Compatibility: Neo-Snap Female Care hose

Part # 2059293-001 (Qty: 10)

# NIBP Supplies



## Male Submin to 1/8" ID Barbed, Plastic

Compatibility: GE, GE Marquette, GE Datex-Ohmeda

Part # 330091 (Qty: 10)



## Male Submin to 5/32" ID Barbed, Plastic

Compatibility: GE, GE Marquette, GE Datex-Ohmeda

Part # 330090 (Qty: 10)



## Mated Pair Submin to 1/8" ID Barbed, Plastic

Compatibility: GE Datex-Ohmeda

Part # 330058 (Qty: 10)



## Female Submin to 1/8" ID Barbed, Plastic

Compatibility: GE, GE Marquette, Protocol

Part # 330064 (Qty: 10)



## 1/8" ID Barbed to Female Screw, Metal

Compatibility: GE DINAMAP

Part # 300619 (Qty: 2)



## 5/32" ID Barbed to Male Slip Luer, Metal

Compatibility: Female Locking Luer and Screw Connectors

Part # 300667 (Qty: 10)

## Wash Plugs and O-Rings

GE Healthcare provides wash plugs that block water from entering the BP cuff during cleaning.



## Washing Plugs for Submin Connector Cuffs

Compatibility: GE, GE Marquette, Protocol

Part # 330072 (Qty: 50)



## Washing Plugs for Screw Connector Cuffs

Compatibility: GE DINAMAP

Part # 300877 (Qty: 50)



## O-Ring Replacement Set for Air Hose Instrument Connector

Compatibility: (330073) and Air Hoses (107365, 107368, 107366)

Part # 330077 (Qty: 20)



## Washing Plugs for Neo-Snap Female Recept and DINACLICK Male Plug

Compatibility: GE

Part # 2059295-001 (Qty: 50)

## Masimo LNOP Reusable Sensors

GE




Interconnect cable with GE connector and Masimo LNOP connector

**2017002-002** 1.2 m (4 ft)

**2017002-003** 2.4 m (8 ft)

**2017002-001** 3.6 m (12 ft)



### Reusable sensors

**2002800-001** Masimo LNOP Reusable Finger sensor LNOP/DCI, Adult

**2002799-001** Masimo LNOP Reusable Finger sensor LNOP/DCIP, Pediatric

**2010463-001** Masimo LNOP Reusable Multisite sensor LNOP-YI

**2027274-001** Masimo LNOP Reusable Tip-Clip Ear sensor LNOP TC-I

**2009745-001** Masimo LNOP Reusable Finger sensor LNOP DC-195, Adult sensor connects directly to the interconnect cable

# SpO<sub>2</sub> Supplies

---

## Masimo LNOP Adhesive Sensors

GE



Interconnect cable with GE connector and Masimo LNOP connector

**2017002-002** 1.2 m (4 ft)

**2017002-003** 2.4 m (8 ft)

**2017002-001** 3.6 m (12 ft)

Sensors connect directly to the interconnect cable

### Adhesive sensors



**2010458-001** Masimo LNOP Disposable Adhesive sensor, LNOP-Adt, Adult (20/box)



**2010459-001** Masimo LNOP Disposable Adhesive sensor, LNOP-Pdt, Pediatric (20/box)



**2017089-001** Masimo LNOP Disposable Adhesive sensor, LNOP-Neo-L, Neonatal (20/box)



**2017090-001** Masimo LNOP Disposable Adhesive sensor, LNOP-NeoPt-L, Neonatal (20/box)



**2027269-001** Masimo LNOP Disposable Adhesive sensor, LNOP-Adtx, Adult, transparent tape (20/box)



**2027270-001** Masimo LNOP Disposable Adhesive sensor, LNOP-Pdtx, Pediatric, transparent tape (20/box)

## Masimo LNCS Reusable & Adhesive Sensors

GE



Masimo LNC-10 interconnect cable  
for Dash, Tram, Corometrics  
**2027263-002** 3 m (10 ft)

### Reusable sensors

**2027258-001** Masimo LNCS DCI Reusable sensor, Adult

**2027259-001** Masimo LNCS DCIP Reusable sensor, Pediatric

**2027261-001** Masimo LNCS TC-I TipClip Reusable Ear sensor

### Adhesive sensors

**2027253-001** Masimo LNCS Adhesive sensor, LNCS-Adtx, Adult (20/box)

**2027254-001** Masimo LNCS Adhesive sensor, LNCS-Pdtx, Pediatric (20/box)

**2051068-001** Masimo LNCS Adhesive sensor, LNCS Inf, Infant (20/box)

**2051069-001** Masimo LNCS Adhesive sensor, LNCS Neo Neonatal (20/box)

**2051070-001** Masimo LNCS Adhesive sensor, LNCS-Neo-Pt, Neonatal (20/box)

## Nellcor OxiMax Sensors and Cables

### Monitors



#### OxiMax

- CARESCAPE Monitor B850 with CARESCAPE Patient Data Module with Nellcor OxiMax™ option with Tram module 451N5/851N5
- CARESCAPE Monitor B650 with CARESCAPE Patient Data Module with Nellcor OxiMax option
- CARESCAPE V100 with Nellcor OxiMax option
- Cardiolab/Mac-Lab with Tram module 451N5
- Dash 2500 with Nellcor OxiMax option
- Dash 3000/4000/5000 with Nellcor OxiMax option
- Solar 8000/8000M/8000i and 9500 with Tram module 451N5/851N5
- Solar 8000M v5/8000i v5 with CARESCAPE Patient Data Module with Nellcor OxiMax option
- E-NSATX module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors



#### OxiMax


- GE DINAMAP ProCare with Nellcor OxiMax option
- GE DINAMAP PRO Series v2 with Nellcor OxiMax option
- GE DINAMAP PRO 1000 v3 with Nellcor OxiMax option
- E-NSAT module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors



Corometrics 250 Series with Nellcor OxiMax option

### Connecting cables

#### GE



Interconnect cable with GE connector and Nellcor OxiMax connector  
**2021406-001** 3 m (10 ft)  
**2021406-002** 1.2 m (4 ft)



Interface cable Nellcor OxiMax DOC-10  
**2008773-001** 2.9 m (9.5 ft)

#### GE



Interconnect cable with GE connector and Nellcor OxiMax connector (for Corometrics only)  
**2025350-001** 3 m (10 ft)

### Sensor type



- 70124021** Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)
- 70124033** Nellcor SpO2 Multisite Sensor D-YS
- 70124034** Nellcor SpO2 Ear-Clip D-YSE, for use with 70124033 (1/box)
- 70124035** Nellcor SpO2 Sensor OXY-A/N, Neonatal/Adult (1/box)
- 414248-001** Nellcor SpO2 Sensor OXI-P/I, Pediatric/Infant (1/box)
- 70124022** Nellcor SpO2 Sensor MAX-P, Adhesive, Pediatric (24/box)
- 70124026** Nellcor SpO2 Sensor MAX-I, Adhesive, Infant (24/box)
- 70124027** Nellcor SpO2 Sensor MAX-A, Adhesive, Adult (24/box)
- 70124032** Nellcor SpO2 Sensor MAX-N, Adhesive, Neonatal (24/box)
- 407705-005** Nellcor SpO2 Sensor MAX-R, Adhesive, Nasal (24/box)
- 2028117-001** Nellcor SpO2 sensor MAX-AL, Adhesive, Adult Long (24/box)



Nellcor OxiSmart Sensors and Cables

Monitors



OxiSmart

CARESCAPE Monitor B850 with Tram module 451/451N/851/851N  
 CardioLab/Mac-Lab with Tram module 451/451N  
 Corometrics 120/2120 Series  
 Dash 3000/4000 (v5 or lower) with legacy Marquette option  
 Dash 3000/4000/5000 with GE Ohmeda option  
 Responder 3000 with Nellcor OxiSmart option  
 Solar 8000/8000M/8000i and 9500 with Tram module 451/451N/851/851N

M-NSAT module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors†  
 Configured monitor: Cardiacap/5 with N-XNSAT Nellcor option





Corometrics 118


† Refer to the CARESCAPE Monitor B850 and CARESCAPE Monitor B650 User Manuals for compatible M-NSAT module versions.

Connecting cables

GE




Sensor type



### GE connector

CARESCAPE® modular monitors and S/5, AS/3 and CS/3 modular monitors with E-(P)RESTN or E-PSM(P) modules. CARESCAPE V100 with TruSignal® option, Responder, B20 and B40 monitors.



**TruSignal TS-G3** Interconnect cable with GE connector  
3 m/10 ft

**TruSignal TS-F2-GE** Integrated Finger sensor with GE connector  
2 m/6.6 ft  
**TS-F4-GE** 4 m/13 ft


**TruSignal TS-SA4-GE** Integrated FingerTip sensor with GE connector  
4 m/13 ft

**TruSignal TS-SP3-GE** Integrated PediTip sensor with GE connector  
3 m/10 ft

**TruSignal TS-E2-GE** Integrated Ear sensors with GE connector  
2 m/6.6 ft  
**TS-E4-GE** 4 m/13 ft

### TruSat connector

GE Ohmeda TruSat® pulse oximeter.



**TruSignal TS-M3** Interconnect cable with TruSat connector  
3 m/10 ft

**TruSignal TS-F4-MC** Integrated Finger sensor with TruSat connector  
4 m/13 ft


**TruSignal TS-SA3-MC** Integrated FingerTip sensor with TruSat connector  
3 m/10 ft

**TruSignal TS-SP3-MC** Integrated PediTip sensor with TruSat connector  
3 m/10 ft

**TruSignal TS-E4-MC** Integrated Ear sensor with TruSat connector  
4 m/13 ft

### Ohmeda connector

Handheld pulse oximeters (GE Ohmeda TuffSat®, 3775), tabletop pulse oximeters (3700 series, 3800, 3900, 3900P), respiratory gas monitors (4700, 5250 RGM), CARESCAPE modular monitors and S/5, AS/3, CS/3 modular monitors with M-OSAT module, Cardiocap®/5 with N-XOSAT option.



**TruSignal TS-H3** Interconnect cable with Ohmeda connector  
3 m/10 ft


**TruSignal TS-F1-H** Integrated Finger sensor with Ohmeda connector  
1 m/3.3 ft  
**TS-F4-H** 4 m/13 ft

**TruSignal TS-E4-H** Integrated Ear sensor with Ohmeda connector  
4 m/13 ft

**TruSignal TS-W4-H** Integrated Wrap sensor with Ohmeda connector  
4 m/13 ft

### Datex connector

CARESCAPE modular monitors and S/5, AS/3 and CS/3 modular monitors with M-NE12ST(P)R, M-(N)EST(P)R or M-(P)RESTN modules; Cardiocap®/5 and S/5 Light Monitors.



**TruSignal TS-N3** Interconnect cable with Datex connector  
3 m/10 ft

**TruSignal TS-F4-N** Integrated Finger sensor with Datex connector  
4 m/13 ft

**TruSignal TS-E4-N** Integrated Ear sensor with Datex connector  
4 m/13 ft

**TruSignal TS-E-D** Ear sensor  
1 m/3.3 ft

**TruSignal TS-F-D** Finger sensor  
1 m/3.3 ft

**TruSignal TS-SA-D** FingerTip sensor  
1 m/3.3 ft

**TruSignal TS-SP-D** PediTip sensor  
1 m/3.3 ft

**TruSignal TS-W-D** Wrap sensor  
1 m/3.3 ft

**TruSignal TS-SE-3** Sensitive Skin sensor  
1 m/3.3 ft 3/box

**TruSignal TS-AF-10** AllFit sensor  
0.5 m/1.7 ft 10/box  
**TS-AF-25** 0.5 m/1.7 ft 25/box

**TruSignal TS-AAW-10** Adult Adhesive Wrap sensor  
0.3m/1ft 10/box  
**TS-AAW-25** 0.3m/1ft 25/box

**TruSignal TS-PAW-10** Pediatric Adhesive Wrap sensor  
0.3m/1ft 10/box  
**TS-PAW-25** 0.3m/1ft 25/box

**TruSignal TS-AP-10** Adult/Pediatric sensor  
0.3m/1ft 10/box  
**TS-AP-25** 0.3m/1ft 25/box



# SpO<sub>2</sub> Supplies

## GE all in one Sensor and Cable



Reference code	Product	Package size	Length	Patient weight	Patient type	Replaces
TS-F-D	TruSignal Finger sensor	1	1 m (3.3 ft)	≥20 kg (44 lb)	Adult, pediatric	OXY-F-UN* OXY-F-DB*
TS-E-D	TruSignal Ear sensor	1	1 m (3.3 ft)	≥10 kg (22 lb)	Adult, pediatric	OXY-E-UN* OXY-E-DB*
TS-W-D	TruSignal Wrap sensor	1	1 m (3.3 ft)	≥3 kg (6 lb)	Adult, pediatric	OXY-W-UN* OXY-W-DB*
TS-SE-3	TruSignal Sensitive Skin sensor	3	1 m (3.3 ft)	No limit	Adult, pediatric, infant, neonate	OXY-SE-3*
TS-SA-D	TruSignal FingerTip sensor	1	1 m (3.3 ft)	>30 kg (66 lb)	Adult	--
TS-SP-D	TruSignal PediTip sensor	1	1 m (3.3 ft)	15-30 kg (33-66 lb)	Pediatric	--

## SpO<sub>2</sub> Interconnect Cables



Reference code	Product	Package size	Length	Patient weight	Replaces
TS-G3	TruSignal interconnect cable with GE connector	1	3 m (10 ft)	No limit	OXY-ES3**
TS-M3	TruSignal interconnect cable with TruSat* connector	1	3 m (10 ft)	No limit	OXY-MC3**
TS-H3	TruSignal interconnect cable with Ohmeda connector	1	3 m (10 ft)	No limit	OXY-OL1, OXY-OL3**
TS-N3	TruSignal interconnect cable with Datex connector	1	3 m (10 ft)	No limit	OXY-SL3, OXY-C1, OXY-C3, OXY-C7 (+OXY-SLA, OXY-SLC)**

## Disposable SpO<sub>2</sub> Sensors

Reference code	Product	Package size	Length	Patient weight	Patient type	Replaces
TS-AF-10	TruSignal AllFit sensor	10	0.5 m (1.7 ft)	No limit	Adult, pediatric, infant, neonate	OXY-AF-10*
TS-AF-25	TruSignal AllFit sensor	25	0.5 m (1.7 ft)	No limit	Adult, pediatric, infant, neonate	OXY-AF-10*
TS-AP-10	TruSignal Adult/Pediatric sensor	10	0.3 m (1 ft)	≥20 kg (≥44 lb)	Adult, pediatric	OXY-AP-10*
TS-AP-25	TruSignal Adult/Pediatric sensor	25	0.3 m (1 ft)	≥20 kg (≥44 lb)	Adult, pediatric	OXY-AP-25*
TS-AAW-10	TruSignal Adult Adhesive Wrap sensor	10	0.3 m (1 ft)	≥20 kg (≥44 lb)	Adult	
TS-AAW-25	TruSignal Adult Adhesive Wrap sensor	25	0.3 m (1 ft)	≥20 kg (≥44 lb)	Adult	
TS-PAW-10	TruSignal Paediatric Adhesive Wrap sensor	10	0.3 m (1 ft)	3-20 kg (6.6-44 lb)	Paediatric, infant	
TS-PAW-25	TruSignal Paediatric Adhesive Wrap sensor	25	0.3 m (1 ft)	3-20 kg (6.6-44 lb)	Paediatric, infant	

\* Requires compatible interconnect cable \*\* Replaces if used with TruSignal sensor




### Certifications





TruSignal sensors and interconnect cables have been tested to and comply with the following standards (including but not limited to): IEC60601-1:2005 and ISO9919:2009

# SpO<sub>2</sub> Supplies

## SpO<sub>2</sub> Sensors - GE Connector

Reference code	Product	Package size	Length	Patient weight	Patient type
 TS-F2-GE	TruSignal Integrated Finger sensor	1	2 m (6.6 ft)	≥20 kg (44 lb)	Adult, pediatric
 TS-F4-GE	TruSignal Integrated Finger sensor	1	4 m (13 ft)	≥20 kg (44 lb)	Adult, pediatric
 TS-SA4-GE	TruSignal Integrated FingerTip sensor	1	4 m (13 ft)	>30 kg (66 lb)	Adult
 TS-SP3-GE	TruSignal Integrated PediTip sensor	1	3 m (10 ft)	15-30 kg (33-66 lb)	Pediatric
 TS-E2-GE	TruSignal Integrated Ear sensor	1	2 m (6.6 ft)	≥10 kg (22 lb)	Adult, pediatric
 TS-E4-GE	TruSignal Integrated Ear sensor	1	4 m (13 ft)	≥10 kg (22 lb)	Adult, pediatric

## SpO<sub>2</sub> Sensors - TruSat Connector

Reference code	Product	Package size	Length	Patient weight	Patient type
 TS-F4-MC	TruSignal Integrated Finger sensor	1	4 m (13 ft)	≥20 kg (44 lb)	Adult, pediatric
 TS-SA3-MC	TruSignal Integrated FingerTip sensor	1	3 m (10 ft)	>30 kg (66 lb)	Adult
 TS-SP3-MC	TruSignal Integrated PediTip sensor	1	3 m (10 ft)	15-30 kg (33-66 lb)	Pediatric
 TS-E4-MC	TruSignal Integrated Ear sensor	1	4 m (13 ft)	≥10 kg (22 lb)	Adult, pediatric

## SpO<sub>2</sub> Sensors - Ohmeda Connector



Reference code	Product	Package size	Length	Patient weight	Patient type
TS-F1-H	TruSignal Integrated Finger sensor	1	1 m (3.3 ft)	≥20 kg (44 lb)	Adult, pediatric
TS-F4-H	TruSignal Integrated Finger sensor	1	4 m (13 ft)	≥20 kg (44 lb)	Adult, pediatric
TS-E4-H	TruSignal Integrated Ear sensor	1	4 m (13 ft)	≥10 kg (22 lb)	Adult, pediatric
TS-W4-H	TruSignal Integrated Wrap sensor	1	4 m (13 ft)	>3kg	Adult, pediatric

## SpO<sub>2</sub> Sensors - Datex Connector



Reference code	Product	Package size	Length	Patient weight	Patient type
TS-F4-N	TruSignal Integrated Finger sensor	1	4 m (13 ft)	≥20 kg (44 lb)	Adult, pediatric
TS-E4-N	TruSignal Integrated Ear sensor	1	4 m (13 ft)	≥10 kg (22 lb)	Adult, pediatric

## SpO<sub>2</sub> Accessories



Reference code	Product	Size	Weight range	Use with
OXY-RTW	Wide adhesive replacement tape	--	No limit	Wrap, Sensitive Skin or AllFit sensor
OXY-RWL	Foam wrap replacement	Large	No limit	Wrap sensor
OXY-RWM	Foam wrap replacement	Medium	No limit	Sensitive Skin sensor
OXY-RWS	Foam wrap replacement	Small	No limit	Sensitive Skin sensor
OXY-RTB	Replacement adhesive tape (bears)	--	No limit	AllFit sensor
OXY-RT	Replacement adhesive tape (blue)	--	No limit	AllFit sensor
OXY-SND	Infant foam sandal	--	2 to 7.7 kg (4.4 to 17 lb)	Sensitive Skin sensor
OXY-HB	Replacement Headband	--	No limit	Ear sensor

### Certifications

TruSignal sensors and cables have been tested to and comply with the following standards (including but not limited to): IEC60601-1:2005 and ISO9919:2009





## Neuromuscular Transmission Sensors & Cables

NMT			
	Reference	Description	Packaging
	888414	<ul style="list-style-type: none"> <li>NMT sensor cable</li> <li>3 m (11 foot)</li> <li>For M and E series and PSM modules</li> </ul>	Each
	888415	<ul style="list-style-type: none"> <li>NMT sensor cable</li> <li>1.5 m (5 foot)</li> </ul>	Each
	888418	<ul style="list-style-type: none"> <li>E-NMT mechanosensor</li> <li>Adults</li> <li>For M and E series modules</li> </ul>	Each
	897439	<ul style="list-style-type: none"> <li>E-NMT paediatric mechanosensor</li> <li>5 to 20 kg</li> <li>For M and E series modules</li> </ul>	Each
	888416	<ul style="list-style-type: none"> <li>Electrosensor</li> <li>Research</li> <li>Other measurement sites than hand</li> </ul>	
	888417	<ul style="list-style-type: none"> <li>Regional block adapter</li> <li>0.5 m (1.5 ft)</li> <li>Stimulation needle connector ID 2 mm</li> </ul>	
	57268-HEL	<ul style="list-style-type: none"> <li>NMT electrodes</li> <li>Disposable</li> </ul>	Box of 30
	8002717	<ul style="list-style-type: none"> <li>NMT Start up kit</li> <li>Includes Mechanosensor, Sensor cable 1.5 m, NMT electrodes</li> </ul>	Each



# Anaesthesia & Respiratory Supplies

## Entropy Sensors & Cables

Entropy			
	Reference	Description	Packaging
	M1050784	<ul style="list-style-type: none"> <li>• GE Entropy cable</li> <li>• For M and E series modules</li> </ul>	Each
	M1038681	<ul style="list-style-type: none"> <li>• GE Entropy sensor</li> <li>• For M and E series modules</li> <li>• replaces 8002858 and M1038681</li> </ul>	Box of 25
	M1174413	<ul style="list-style-type: none"> <li>• GE Entropy EasyFit sensor</li> <li>• Disposable</li> <li>• Single use</li> </ul>	Box of 25
	M1174414	<ul style="list-style-type: none"> <li>• Entropy EasyFit sensor start-up kit</li> <li>• Entropy EasyFit sensor M1174413 (2 boxes), GE Entropy cable (1 Each)</li> </ul>	Each
	M1060609	<ul style="list-style-type: none"> <li>• GE Entropy start-up kit</li> <li>• GE Entropy sensor 2 Boxes and GE Entropy cable 1 Each</li> </ul>	Each





# Anaesthesia & Respiratory Supplies

## Gas Monitoring

	D-Fend Pro	D-Fend Pro+	Mini D-Fend	D-Fend	D-Fend+
					
Colour	Dark Steel Blue	Green		Black	Green
Application	Anaesthesia	Critical Care	Anaesthesia/ Critical Care		
Recommended with	E-sCAiO, E-sCAiOV	E-sCO, E-sCOV	E-MiniC	E-CAiO, E-CAiOV, E-CAiOVX	E-CO, E-COV, E-COVX
Replacing Interval	2 months	24 hours/Each new patient	Depending on use	2 months	24 hours/Each new patient
Bacterial filtration efficiency	> 99.999%	> 99.999 %			
Viral filtration efficiency	> 99.999%	> 99.999 %			
Order Code	M1182629	M1200227	8002174	876446-HEL	881319-HEL
Package Size	10/case	10/case	10/case	10/case	10/case

## Calibration Gases & Regulators

	Anaesthesia	Critical Care	Co2, E-miniC	Gas exchange in Critical Care
	<b>755534-HEL</b> Calibration Gas Regulator + <b>755583-HEL</b> Calibration Gas, blue 5% CO <sub>2</sub> , 33% N <sub>2</sub> O, 2% Desfl. 55% O <sub>2</sub> , Balance N <sub>2</sub>	<b>755534-HEL</b> Calibration Gas Regulator + <b>755581-HEL</b> Calibration Gas, green 5% CO <sub>2</sub> , 40% N <sub>2</sub> O, 55% O <sub>2</sub>	<b>755534-HEL</b> Calibration Gas Regulator + <b>755580</b> Calibration Gas, grey 5% CO <sub>2</sub> , Balance air	<b>M1006864</b> Calibration Gas Regulator + <b>755587</b> Calibration Gas 5% CO <sub>2</sub> , Balance O <sub>2</sub>
<b>E-sCO</b>		x		
<b>E-sCOV</b>		x		
<b>E-sCAiO</b>	x			
<b>E-sCAiOV</b>	x			
<b>E-CO</b>		x		
<b>E-COV</b>		x		
<b>E-COVX</b>				x
<b>E-CAiO</b>	x			
<b>E-CAiOV</b>	x			
<b>C-CAiOVX</b>	x			
<b>E-miniC</b>			x	

Use the calibration gas recommended for your module, and follow the recommended calibration intervals to ensure that the measurement accuracy remains within specifications.



# Anaesthesia & Respiratory Supplies

## Anaesthetic Breathing Circuit Spares

### Anaesthetic Breathing Circuit Spares <sup>◇</sup><sup>1</sup>

Reference	Description	Packaging
1407-7008-000	<ul style="list-style-type: none"> <li>Replacement ABS breathing circuit</li> <li>Includes all ABS parts except flow sensors, O2 cell and CO2 canisters</li> </ul>	
1407-7016-000	<ul style="list-style-type: none"> <li>Replacement ABS breathing circuit</li> <li>For Australia and New Zealand</li> </ul>	
1407-7021-000	<ul style="list-style-type: none"> <li>EZchange™ Canister manifold</li> <li>Requires rail design 1407-7004-000 reusable CO2 canister or 8003138, 8003693 disposable CO2 canisters</li> </ul>	
1407-7027-000	<ul style="list-style-type: none"> <li>EZchange Canister manifold and condenser option</li> <li>Includes 1407-7021-000 EZchange manifold, 1407-7024-000 condenser and 1407-3145-000 guard. Requires rail design 1407-7004-00 reusable CO2 canister or 8003138, 8003693 disposable CO2 canisters</li> </ul>	
1407-7004-000	<ul style="list-style-type: none"> <li>Rail design CO2 canister - reusable<sup>2</sup></li> <li>Includes qty 1, 1407-3201-000</li> </ul>	
1407-3201-000	<ul style="list-style-type: none"> <li>Rail design CO2 canister foam inserts</li> <li>For use with reusable CO2 canister</li> </ul>	
6050-0004-110	<ul style="list-style-type: none"> <li>O2 sensor</li> </ul>	
1503-3857-000	<ul style="list-style-type: none"> <li>O2 sensor plug kit</li> <li>Optional - included if no O2 monitoring present</li> </ul>	
1009-5570-000	<ul style="list-style-type: none"> <li>O2 sensor replacement cable</li> </ul>	

## Flow Sensor Spares

### Flow Sensor Spares <sup>◇</sup><sup>3</sup>

Reference	Description	Packaging
1407-7022-000	<ul style="list-style-type: none"> <li>Flow sensor module</li> <li>Includes flow sensor cover and module. Does not include flow sensors</li> </ul>	
1503-3858-000	<ul style="list-style-type: none"> <li>Offset flow sensor - disposable</li> </ul>	
1503-3856-000	<ul style="list-style-type: none"> <li>Flow sensor - disposable</li> </ul>	
1503-3244-000	<ul style="list-style-type: none"> <li>Flow sensor - autoclavable</li> </ul>	

1 One multi-absorber disposable CO2 canister is included in the supplies kit with every machine. To order additional disposable or reusable CO2 canisters, see the supplies and accessories section.

2 Effective 4/05 this part design changed to a rail design configuration

3 Two flow sensors required

\* Required

◇ Optional

## Spirometry

### D-lite flow sensors

Sensor	D-lite Reusable	Pedi-lite Reusable	D-lite+ Disposable	D-lite Disposable	Pedi-lite+ Disposable
Application	Adults	Paediatric	Adults - humid conditions	Adults	Paediatric - humid conditions
HME recommended	Yes	Yes	No	Yes	Yes
Measurement range (Tidal volume)	150-2000 ml	15-300 ml	150-2000 ml	150-2000 ml	15-300 ml
Accuracy	±6% or 30 ml	±6% or 4 ml	±6% or 30 ml	±6% or 30 ml	±6% or 4 ml
Dead space	9.5 ml	2.5 ml	9.5 ml	9.5 ml	2.5 ml
Resistance	0.5 cm H <sub>2</sub> O/30 l/min	1 cm H <sub>2</sub> O/ 15 l/min	0.5 cm H <sub>2</sub> O/30 l/min	0.5 cm H <sub>2</sub> O/30 l/min	1 cm H <sub>2</sub> O/ 15 l/min
Material	Polyphenylsulfone	Polyphenylsulfone	Polystyrene	Polystyrene	Polystyrene
Colour	Transparent yellow	Transparent yellow	Transparent	Yellow	Transparent
REF	733910-HEL	73393	896952	733950	8001948
Package Size	Each	Each	50/case	50/case	50/case



### Spirometry tubes, single use

REF	Length
890031	2 m/7 ft
884101	3 m/10 ft
891191-HEL	6 m/20 ft



Note: If you use the module for measuring metabolism, use the 2 m spirometry tube and gas sampling line.

### Spirometry kits, single use

REF	Application
894255	Critical Care
889560	Anaesthesia
8002718	Paediatric, for humid conditions (2 m/7 ft)
8004381	For humid conditions (2 m/7 ft)
8004382	For humid conditions (3 m/10 ft)
M1032634	Paediatric, for humid conditions (3 m/10 ft)

### Maintenance

	1. Cleaning		2. Disinfection		3. Sterilisation max. temp. 132°C
	Machine wash	Wipe with ordinary detergent	Machine wash	Immerse in disinfectant	Steam autoclave
D-lite flow sensor	x	o	x	o	o
Pedi-lite flow sensor	x	o	x	o	o

x recommended o if necessary

# Anaesthesia & Respiratory Supplies

## Ancillaries

Ancillaries			
	Reference	Description	Packaging
	1001-3269-000	<ul style="list-style-type: none"><li>• Caster guards</li><li>• Set of 4</li><li>• Grey</li><li>• for 5" caster</li></ul>	
	1009-3086-000	<ul style="list-style-type: none"><li>• Hook, breathing circuit</li></ul>	
	1009-3252-000	<ul style="list-style-type: none"><li>• Patient Fingers</li><li>• OEM part</li><li>• Retainer Multiple CA large</li><li>• Plastic</li></ul>	
	M1074608	<ul style="list-style-type: none"><li>• Worksurface cover, kit</li></ul>	
	1500-3378-000	<ul style="list-style-type: none"><li>• Bellows Adult</li></ul>	
	1001-3854-000	<ul style="list-style-type: none"><li>• Krytox performance lubricant</li></ul>	






# Consumables distributed by CareFusion (AU)/Obex (NZ)

---




## Diagnostic Cardiology

ECG recording paper  
ECG interface cables  
Leadwires  
Electrodes



## Temperature Management

Temperature Probes  
Connecting Cables




## Airway Management

Laryngoscope Handles  
Laryngoscope Blades



## Maternal-Infant-Care

Recording Paper  
Corometrics Transducers  
Bilipads



## Gas Monitoring

Sampling Lines  
Face masks  
Intubation  
Breathing Circuits

Filters  
CO<sub>2</sub> Absorption

### Order details for Australia

CareFusion Australia 316 Pty Ltd  
Email: [au-customerservice@carefusion.com](mailto:au-customerservice@carefusion.com)  
Phone: 1800 833 372  
Fax: 1800 833 518

### Order details for New Zealand

Obex Ltd.  
Email: [orders@obex.co.nz](mailto:orders@obex.co.nz)  
Customer Enquiries; Mary Andrew 09 631 1081  
Freephone 0800 656 239



**imagination at work**

GE Healthcare  
32 Phillip Street  
Parramatta, NSW 2150  
Australia  
1300 722 229  
[www.gehealthcare.com.au](http://www.gehealthcare.com.au)

GE Healthcare  
8 Tangihua Street  
Auckland 1010  
New Zealand  
0800 43 43 25  
[www.gehealthcare.co.nz](http://www.gehealthcare.co.nz)

GE Healthcare Limited is a subsidiary of the General Electric Company. The GE brand, logo and Lumination are trademarks of the General Electric Company.

© 2015 GE Healthcare Limited

All information in this document was correct at the time of publication and is subject to change without notice.

JB29447AU - Revision 1, April 2015