

# Mobile Detector Holders Comparison Guide



## Flexi-Holder II for GE Cassette or Flashpad™

- **Scan Flexibility** – Allows accurate radiographs to be taken from virtually any angle from floor level to 71"+ (180 cm) high
- **Versatility** – Accepts most cassettes and does not obscure cassette film area
- **Exclusive GE Design** – Designed for use with the GE Flashpad detector or GE cassette



## Mobile Adjustable Detector Holder for Flashpad

- **Lower Extremity Exams** – Maximize patient positioning during chest, spine and lower extremities exams
- **Ease of Use** – Electric height adjustment allows the detector to be positioned at any height with the touch of a button
- **Improved Positioning Accuracy** – Straps and clear acrylic panel improves the accuracy of patient positioning



## Mobile Detector Holder for Flashpad HD Detector

- **Compatibility** – Compatible with the GE FlashPad® HD detector
- **Scan Flexibility** – A stable platform enabling images to be taken at any angle
- **Easy to Use** – Counterbalanced vertical movement from floor level to 50.8" (129 cm) high



**Flashpad**



**Flashpad**



**Flashpad  
HD**



**Flexi-Holder II for GE Cassette or Flashpad™**

- Manual height and angle adjustments
- Tilt the detector/cassette 180 degrees facing straight up, down, or any angle in-between
- Telescopes up to 26" over the X-ray or operating room table
- Teflon-coated bracket will protect the detector while securing it in place
- Size: 22"W x 19"D x 75"H
- Weight: 110 lbs

**Mobile Adjustable Detector Holder for Flashpad**

- Height adjustable - detector and grid can move from 8" to 47" above the floor
- Detector and grid can be rotated 90° in either portrait or landscape
- Straps and a removable clear acrylic panel improve the accuracy of patient positioning
- Small footprint with large 4" locking wheels provides for easy movement room to room with detector and grid in place
- Size: 28.5"W x 26.25"D x 68"H
- Weight: 150 lbs

**Mobile Detector Holder for Flashpad HD Detector**

- FlashPad HD (or a cassette) can overhang the table or trolley up to 24.4" (62 cm)
- Strong and compact design
- Large, quick-locking castor wheels
- Size: 70.5" H x 26.4" D x 24" W (179 cm H x 67 cm D x 61 cm W)
- Weight: 121 lbs (55kg)

Catalog #	Description
E6314J	Flexi-Holder II for Cassette or Flashpad (not compatible with Flashpad HD)

Catalog #	Description
E6313DA	Mobile vertical detector holder for Flashpad digital detector, 120V, 60Hz, (not compatible with Flashpad HD)

Catalog #	Description
E06731BL	Mobile Detector Holder for Flashpad HD Detector

Product may not be available in all countries and regions. Full product technical specification is available upon request. Contact a GE Healthcare Representative for more information. Please visit [www.gehealthcare.com/promotional-locations](http://www.gehealthcare.com/promotional-locations)

Data subject to change.

© 2020 General Electric Company. March, 2020/JB76647US

GE and the GE Monogram are trademarks of General Electric Company.

GE Healthcare  
3000 North Grandview  
Waukesha, WI 53188  
USA  
[www.gehealthcare.com](http://www.gehealthcare.com)

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

