



SENSA-CUF

Blood Pressure Cuffs

The CRITIKON™ SENSA-CUF™ is a highly versatile solution for a wide range of extended-care settings. Made using strong, breathable, easy-to-clean material, the SENSA-CUF is designed to maximize patient comfort in low-acuity areas in a hospital where risk of contamination are diminished, and durable enough for multi-patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling a more accurate measurement

Constructed for patient safety and durability

- Soft, breathable, easy-to-clean material is not made with PVC, making them safer for the environment, or natural rubber latex, making them safer for your patients
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers, enabling a universal connection system throughout the facility

Assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



SENSA-CUF

ISO 80369-5 compliant



Size	Color	Limb Circumference	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	2470
Child	Green	12-19 cm	2472
Small Adult	Light Blue	17-25 cm	2474
Small Adult Long	Light Blue	17-25 cm	2475
Adult	Navy	23-33 cm	2476
Adult Long	Navy	23-33 cm	2477
Large Adult	Rose	31-40 cm	2478
Large Adult Long	Rose	31-40 cm	2479
Thigh	Brown	38-50 cm	2480

Note: 5 cuffs/box

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, SENA-CUF and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00474XX 8/20