



DINACLICK

GE Healthcare's universal connection system for dual-tube CRITIKON™ blood pressure cuffs

Innovative and easy-to-use design allows clinicians to quickly connect/disconnect adult and pediatric blood pressure cuffs and provides standardization of your monitors and manual NIBP devices. The audible click lets users know there is a proper connection, and the ergonomically-designed system makes the DINACLICK™ nearly effortless to use.

DINACLICK connector benefits

- **Designed for safety.** Reduces misconnections with luer-style connectors
- **Improved workflow.** Decreases manual connection time
- **Universal connectivity.** Standardization solution for single-tube and dual-tube hoses
- **Intuitive design.** Easy-to-use
- **Non-metallic components.** Material design allows for use in MR environment
- **Patient safety.** Not made with natural rubber latex or DEHP
- **ISO-80369-5 compliant.** DINACLICK connectors are compliant with ISO 80369-5 which regulates NIBP connector types to prevent misconnections and patient risks



To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

NIBP cuff part numbers

DINACLICK ISO 80369-5 compliant						
Size	Color	Limb Circumference	SOFT-CUF	CLASSIC-CUF	CLASSIC-CUF ISO	DURA-CUF
Infant	Orange	8-13 cm	SFT-P1-2A-5S	CLA-P1-2A-5S	CLA-P1-2A-X-5S	DUR-P1-2A-5S
Child	Green	12-19 cm	SFT-P2-2A-5S	CLA-P2-2A-5S	CLA-P2-2A-X-5S	DUR-P2-2A-5S
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S			
Small Adult	Light Blue	17-25 cm	SFT-A1-2A-5S	CLA-A1-2A-5S	CLA-A1-2A-X-5S	DUR-A1-2A-5S
Small Adult Long	Light Blue	17-25 cm	SFT-A1-2A-L-5S			DUR-A1-2A-L-5S
Adult	Navy	23-33 cm	SFT-A2-2A-5S	CLA-A2-2A-5S	CLA-A2-2A-X-5S	DUR-A2-2A-5S
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S	CLA-A2-2A-L-5S		DUR-A2-2A-L-5S
Large Adult	Rose	31-40 cm	SFT-A3-2A-5S	CLA-A3-2A-5S	CLA-A3-2A-X-5S	DUR-A3-2A-5S
Large Adult Long	Rose	31-40 cm	SFT-A3-2A-L-5S	CLA-A3-2A-L-5S		DUR-A3-2A-L-5S
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	CLA-T1-2A-5S	CLA-T1-2A-X-5S	DUR-T1-2A-5S

RADIAL-CUF		ISO 80369-5 compliant
Color	Forearm Circumference	2-Tube DINACLICK Connector
Purple	26-36 cm ¹	SFT-F1-2A-5S

DINACLICK connectors

	Part number	Description	Quantity/package
	2059294-001	2T DINACLICK FEMALE RECEPT	10/PACK
	2055630-001	1T DINACLICK 1/8" TUBE FEMALE RECEPT	10/PACK
	2055630-002	1T CONN 2-1 DINACLICK 5/32" TUBE FEMALE RECEPT	10/PACK

Wall mamometers and accessories

Part number	Description	Quantity/package
2059289-001	INCLUDES COILED TUBING, ADAPTER WITH DINACLICK CONNECTOR, BULB AND VALVE	5/BOX
2059289-002	INCLUDES ADAPTER WITH DINACLICK CONNECTOR AND 13" TUBING (NO COILED TUBING, NO BULB AND VALVE)	5/PACK

Adapters

	Part number	Description	Quantity/package
	2059291-001	Y-ADAPTER: 1T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-002	Y-ADAPTER: 1T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-003	Y-ADAPTER: 1T MALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-004	Y-ADAPTER: 1T UNFITTED TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-001	X-ADAPTER: 2T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-002	X-ADAPTER: 2T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-003	X-ADAPTER: 2T MATED SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-004	X-ADAPTER: 2T MALE SCREW TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059295-001	WASH PLUGS – FOR NEO-SNAP FEMALE & DINACLICK MALE	50/PACK

Adult/pediatric NIBP hoses

	Part number	Description	Length	Product compatibility
	2058203-001	CARE HOSE ADULT DINACLICK	4 FT/1.2 M	B40, B105, B125, DASH 3000, 4000, 5000, TRAM, CARESCAPE™ PATIENT DATA MODULE, COROMETRICS 120, 128, 129, 259, 259CX, E-PSM(P), E-PRESTN, E-PRETN, E-RESTN
	2058203-002	CARE HOSE ADULT DINACLICK	12 FT/3.6 M	
	2058203-003	CARE HOSE ADULT DINACLICK	24 FT/7.2 M	
	2058205-001	CRITIKON HOSE ADULT DINACLICK	12 FT/3.6 M	DINAMAP PRO SERIES, DINAMAP PROCARE, CARESCAPE V100, VC150, COMPACT, MPS SERIES MONITORS
	2058205-002	CRITIKON HOSE ADULT DINACLICK	24 FT/7.2 M	
	2087389-001	NIBP Hose CARESCAPE ONE DINACLICK	4 FT/1.2 M	CARESCAPE ONE
	2087389-002	NIBP Hose CARESCAPE ONE DINACLICK	12 FT/3.6 M	
	2087389-003	NIBP Hose CARESCAPE ONE DINACLICK	24 FT/7.2 M	



© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DINACLICK, and CARESCAPE are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB01312XX 10/20